


The Eliminator 300 SERIES: Thin Wall

APPLICATIONS: K-Middle School

EH-300

STANDARD FEATURES

- Interior Baffle Arrangement:
 - Fixed baffles (standard)
- Exhaust Collar:
 - Type 304 s/s - 10"Dia. Round (standard)
 - Other _____
- Interior Lighting Fixture:
 - Standard T-8 Fluorescent Light w/ Bulbs (standard)
- Sash Frame & Glass:
 - Painted Steel Frame w/ Laminated Safety Glass (standard)
- Sash Type:
 - Vertical (standard)



**** DECK MOUNT PLUMBING FIXTURES AVAILABLE**

- ALARM (OPTION) AS SHOWN OPPOSITE SIDE
- BLOWER SWITCH
- LIGHT SWITCH

| HOOD BASE CABINETS: | HEIGHT | | | |
|---------------------|--------|-----|-----|-----|
| | Width | 29" | 33" | 35" |
| Acid | | | | |
| Flammable | | | | |
| Vacuum Pump | | | | |
| Standard Base | | | | |
| Table Base Frame | | | | |

NOTE: 18" DEEP FLAMMABLE AND/OR ACID CABINET IS NEEDED IF DECK MOUNT FIXTURES ARE REQUIRED

- WORK SURFACE MATERIAL:**
- Epoxy (standard)
 - Type 304 Stainless Steel
 - Type 316 Stainless Steel
 - No Worksurface Required

- LINER MATERIAL:**
- POLY RESIN (standard)
 - 304 STAINLESS STEEL
 - 316 STAINLESS STEEL
 - PVC
 - SPECIAL

- BY-PASS TYPE:**
- OPEN BY-PASS (standard)
 - RESTRICTED BY PASS

- PRE-WIRING**
- Yes
 - No

LENGTH

36" 48" 60" 72" 96"

Custom sizes available; contact AMS

SPECIAL REQUIREMENTS: _____

OPTIONAL ACCESSORIES

- Lattice Assembly:
- None
 - Aluminum Lattice Assembly
 - Stainless Steel Assembly
- Interior Lighting Fixture:
- Explosion Proof Light
- Fume Hood Alarm Option:
- Alarm Tel AFA_500
 - Alarm Tel AFA_1000
- Sash Frame & Glass:
- Stainless Steel Frame w/ Laminated Safety Glass
- Ceiling Enclosure:
- Yes (Height: _____) No
- Finished Back:
- Yes No
- Sash Type:
- Combination
 - Horizontal
 - Interlocking
- Additional Options:
- | | | |
|-----------------------|------------------------------|-----------------------------|
| Blower | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Sash Stops | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Push Button Sash | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Automatic Sash Return | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Sash Locks | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Sash Alarm | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Add Air Chamber | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
- ** Specify Service Fixtures (Deck Mounted Type):**
- _____
- _____



AIR MASTER SYSTEMS CORP
 6480 NORTON CENTER DR.
 MUSKEGON, MI. 49441
 TEL. (231) 798-1111 FAX (231) 798-4000
 www.airmastersystems.com

CUSTOMER

Qty: _____

Color: _____