

# The Eliminator 600 SERIES: Radio Isotope Hood

APPLICATIONS: K-12; Univresity; Private Industry

## EH-600

### STANDARD FEATURES

- Interior Baffle Arrangement:
  - Fixed baffles (standard)
- Exhaust Collar:
  - Type 304 s/s - 10"Dia. Round (standard)
  - Other \_\_\_\_\_
- Interior Lighting Fixture:
  - Standard T-8 Fluorescent Light w/ Bulbs (standard)
- Sash Frame & Glass:
  - Stainless Steel Frame w/ Laminated Safety Glass (standard)
- Sash Type:
  - Vertical (standard)
- Additional Features:
  - Stainless Steel Airfoil
  - Type 304 Stn. Stil. one-piece integral work surface & liner
- Pre-Piping / Pre-Wiring (standard)
  - All hoods are pre-wried & pre-piped unless other wise noted

(SPECIFY FITTINGS / FIXTURES REQUIRED)

(SPECIFY FITTINGS / FIXTURES REQUIRED)

SERVICE

SERVICE

SERVICE

SERVICE

SERVICE

SERVICE

SERVICE

SERVICE

BLOWER SWITCH (STANDARD)

ALARM (OPTION)

LIGHT SWITCH (STANDARD)

SERVICE

SERVICE

120V GFI DUPLEX

120V GFI DUPLEX

#### LINER MATERIAL:


- 304 STAINLESS STEEL (standard)
- 316 STAINLESS STEEL

#### PRE-PIPING:

- PRE-PIPING UP
- PRE-PIPING DOWN

#### BY-PASS TYPE:

- OPEN BY-PASS (standard)
- RESTRICTED BY PASS



36"  48"  60"  72"  96"

SPECIAL REQUIREMENTS: \_\_\_\_\_

\_\_\_\_\_

HOOD BASE CABINETS:	HEIGHT			
	Width	29"	33"	35"
Acid				
Flammable				
Vacuum Pump				
Standard Base				
Table Base Frame				

#### WORK SURFACE MATERIAL:

- Type 304 Stainless Steel (standard)
- Type 316 Stainless Steel

### OPTIONAL ACCESSORIES

- Lattice Assembly:
  - None
  - Aluminum Lattice Assembly
  - Stainless Steel Assembly
- Interior Lighting Fixture:
  - Explosion Proof Light
- Fume Hood Alarm Option:
  - Alarm Tel AFA\_500
  - Alarm Tel AFA\_1000
  - VAV Cutout (Customer provides cutout detail)
- Sash Type:
  - Combination Sash (Horizontal/Vertical)
  - Horizontal
  - Interlocking
- Sash Frame & Glass:
  - Stainless Steel Frame w/ Laminated Safety Glass
- Airfoil:
  - Type 304 Stainless Steel
  - Type 316 Stainless Steel
- Trough:
  - None (standard)
  - Type 304 Stainless Steel
  - Type 316 Stainless Steel
- Ceiling Enclosure:
  - Yes (Height: \_\_\_\_\_)  No
- Finished Back:
  - Yes  No
- Additional Options:
 

Blower	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Push Button Sash	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Automatic Sash Return	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Sash Locks	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Sash Alarm	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Add Air Chamber	<input type="checkbox"/> Yes	<input type="checkbox"/> No



## AIR MASTER SYSTEMS CORP

6480 NORTON CENTER DR.  
MUSKEGON, MI. 49441

TEL (231) 798-1111 FAX (231) 798-4000  
www.airmastersystems.com

CUSTOMER

Qty: \_\_\_\_\_

Color: \_\_\_\_\_