


The Eliminator 900 SERIES: Demonstration Hood

APPLICATIONS: K-12



(SPECIFY FITTINGS / FIXTURES REQUIRED)

BLOWER SWITCH (OPTION)

LIGHT SWITCH (STANDARD)

120V GFI DUPLEX

Please specify Left or Right
 Left Right
 (Left-hand shown)

HOOD BASE CABINETS:	HEIGHT			
	Width	29"	33"	35"
Acid				
Flammable				
Vacuum Pump				
Standard Base				
Table Base Frame				

- LINER MATERIAL:
- POLY RESIN (standard)
 - 304 STAINLESS STEEL
 - 316 STAINLESS STEEL
 - PVC
 - SPECIAL

48" 60"
 (Other sizes available, contact AMS.)

FRONT

SPECIAL REQUIREMENTS: _____

*Note: Plumbing services must be deck mounted.

EH-900 STANDARD FEATURES

- Interior Baffle Arrangement:
 - Fixed baffles (standard)
- Exhaust Collar:
 - Type 304 s/s - 10" Dia. Round (standard)
 - Other _____
- Interior Lighting Fixture:
 - Standard T-8 Fluorescent Light w/ Bulbs (standard)
- Sash Frame & Glass:
 - Painted Steel Frame w/ Laminated Safety Glass (standard)
- Sash Type:
 - Vertical (standard)
 - Sash Interlock (standard)
- Pre-Wiring (standard)
 - All hoods are pre-wired unless other wise noted

OPTIONAL ACCESSORIES

- Lattice Assembly:
 - None
 - Plastic Lattice Assembly
 - Aluminum Lattice Assembly
 - Stainless Steel Assembly
- Interior Lighting Fixture:
 - Explosion Proof Light
- Fume Hood Alarm Option:
 - Alarm Tel AFA_500
 - Alarm Tel AFA_1000
- Sash Frame & Glass:
 - Stainless Steel Frame w/ Laminated Safety Glass
- Ceiling Enclosure:
 - Yes (Height: _____) No
- Sash Type:
 - Combination
 - Horizontal
- Additional Options:

Blower	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Sash Stops	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Sash Locks	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Sash Alarm	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Add Air Chamber	<input type="checkbox"/> Yes	<input type="checkbox"/> No
- Service Fixture (Deck Mount Type):
 - Specify _____



AIR MASTER SYSTEMS CORP
 6480 NORTON CENTER DR.
 MUSKEGON, MI. 49441
 TEL. (231) 798-1111 FAX (231) 798-4000
 www.airmastersystems.com

CUSTOMER

Qty: _____
 Color: _____