


The SOLUTION: Fume Hood (60-100 fpm)

APPLICATIONS: K-12; Higher Education; Industrial; Bio Tech / R&D



LINER MATERIAL:

- POLY RESIN (standard)
- 304 STAINLESS STEEL
- 316 STAINLESS STEEL
- PVC
- SPECIAL

PRE-PIPING:

- PRE-PIPING UP
- PRE-PIPING DOWN

Optional FH Depth:

- 39.75" deep w/ 36"d epoxy dished top (add -40"D suffix to part number)

SPECIAL REQUIREMENTS: _____

WORK SURFACE MATERIAL:

- Epoxy (standard)
- Type 304 Stainless Steel
- Type 316 Stainless Steel

EH-SOL STANDARD FEATURES

- Interior Baffle Arrangement:
 - Fixed baffles (standard)
- Exhaust Collar:
 - Type 304 s/s - 10" Dia. Round (standard)
 - Other _____
- Interior Lighting Fixture:
 - Standard T-8 Fluorescent Light w/ Bulbs (standard)
- Sash Frame & Glass:
 - Frameless w/ Laminated Safety Glass (standard)
- Airfoil:
 - Painted Steel (standard)
- Sash Type:
 - Vertical (standard)
- Pre-Piping / Pre-Wiring (standard)
 - All hoods are pre-wired & pre-piped unless other wise noted
- Sash Drive (standard)
 - Chain & Sprocket Sash Drive
- Sight Line (standard)
 - 36-inch high viewing sightline.

OPTIONAL ACCESSORIES

- Lattice Assembly:
 - None
 - Aluminum Lattice Assembly
 - Stainless Steel Assembly
- Interior Lighting Fixture:
 - Explosion Proof Light
- Fume Hood Alarm Option:
 - Alarm Tel AFA_500
 - Alarm Tel AFA_1000
 - VAV Cutout (Customer provides cutout detail)
- Sash Full-Length Finger Lift:
 - Stainless Steel
- Airfoil:
 - Type 304 Stainless Steel
 - Type 316 Stainless Steel
- Trough:
 - None (standard)
 - Painted Steel
 - Type 304 Stainless Steel
 - Type 316 Stainless Steel
- Ceiling Enclosure:
 - Yes (Height: _____) No
- Finished Back:
 - Yes No
- Sash Type:
 - Combination
 - Horizontal
 - Interlocking
- Additional Options:

Blower	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Sash Stops	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Push Button Sash	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Automatic Sash Return	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Sash Locks	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Sash Alarm	<input type="checkbox"/> Yes	<input type="checkbox"/> No
- Service Fixture Type:
 - Ball valve w/ lever handles
 - Needle valve w/ four-arm handle



AIR MASTER SYSTEMS CORP

6480 NORTON CENTER DR.
MUSKEGON, MI. 49441

TEL. (231) 798-1111 FAX (231) 798-4000
www.airmastersystems.com

CUSTOMER

Qty: _____

Color: _____