




# Electrical Accessories & Specifications

Pre-Wire service is completed in accordance with NEC and UL-1805 standards


Note: If the fume hood is not pre-wired, wiring in the field must be accomplished using UL listed electrical fixtures while observing NEC standards and local electrical codes.




**Single Pole  
15 Amp Light Switch**  
Standard Black Color  
and Cover Plate  
Part # LSW




**Red Illuminated  
Toggle Switch**  
with Black Cover Plate  
Part # BLSW




**Explosion Proof Outlet**  
Class 1, Group 1  
Part # XPO



**Explosion Proof Switch**  
Class 1, Group 1  
Part # XSW



**208 Volt, 20 Amp  
Single Pole Receptacle**  
230, 277 or other  
available upon request  
Part # 208



**120 Volt, 20 Amp Ground Fault  
Interrupter (GFI)**  
Duplex Outlet  
Standard Black color with cover  
Part # GFI



**Explosion Proof Light**  
Class 1, Group 1  
Available in  
Incandescent or fluorescent  
Part # XPL