

The Eliminator 400 SERIES: Walk-In Hood

APPLICATIONS: University; Research; Private Industry

EH-400

STANDARD FEATURES

- Interior Baffle Arrangement:
 - Fixed baffles (standard)
- Exhaust Collar:
 - Type 304 s/s - 10" Dia. Round (standard)
 - Other _____
- Interior Lighting Fixture:
 - Standard T-8 Fluorescent Light w/ Bulbs (standard)
- Sash Frame & Glass:
 - Painted Steel Frame w/ Laminated Safety Glass (standard)
- Sash Type:
 - Vertical (standard)
- Pre-Piping / Pre-Wiring (standard)
 - All hoods are pre-wired & pre-piped unless other wise noted

OPTIONAL ACCESSORIES

- Lattice Assembly:
 - None
 - Plastic Lattice Assembly
 - Aluminum Lattice Assembly
 - Stainless Steel Assembly
- Interior Lighting Fixture:
 - Explosion Proof Light
- Fume Hood Alarm Option:
 - Alarm Tel AFA_500
 - Alarm Tel AFA_1000
 - VAV Cutout (Customer provides cutout detail)
- Sash Frame & Glass:
 - Stainless Steel Frame w/ Laminated Safety Glass
- Ceiling Enclosure:
 - Yes (Height: _____) No
- Sash Type:
 - Combination
 - Horizontal Hanging Doors
 - Interlocking
- Additional Options:

| | | |
|-----------------------|------------------------------|-----------------------------|
| Blower | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Sash Stops | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Sash Locks | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Sash Alarm | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Add Air Chamber | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Panel Mounted Cupsink | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Epoxy Flooring | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Stainless Steel Ramp | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
- Service Fixture Type:
 - Ball valve w/ lever handle (ADA)
 - Needle valve w/ four-arm handle

(SPECIFY FITTINGS / FIXTURES REQUIRED)

(SPECIFY FITTINGS / FIXTURES REQUIRED)

SERVICE

SERVICE

SERVICE

ALARM (OPTION)

SERVICE

BLOWER SWITCH (OPTION)

SERVICE

120V GFI DUPLEX

SERVICE

SERVICE

SERVICE

VAV CONTROLLER CUTOUT
(CUSTOMER PROVIDE DETAIL)

SERVICE

LIGHT SWITCH (STANDARD)

SERVICE


120V GFI DUPLEX

LINER MATERIAL:

- POLY RESIN (standard)
- 304 STAINLESS STEEL
- 316 STAINLESS STEEL
- PVC
- SPECIAL

PRE-PIPING:

- PRE-PIPING UP
- PRE-PIPING DOWN



LENGTH
 48" 60" 72" 96" 120"
 (35" Standard Depth)

Special depth required Yes No
 Specify _____

SPECIAL REQUIREMENTS: _____

** For proper installation & functionality of fume hood, flooring must be leveled by owner prior to delivery of fume hoods.



AIR MASTER SYSTEMS CORP

6480 NORTON CENTER DR.
 MUSKEGON, MI. 49441

TEL. (231) 798-1111 FAX (231) 798-4000
 www.airmastersystems.com

CUSTOMER

Qty: _____

Color: _____