

# The Eliminator 500 SERIES: ADA Design

APPLICATIONS: K-12; University; Private Industry

## EH-100


### STANDARD FEATURES

- Interior Baffle Arrangement:
  - Fixed baffles (standard)
- Exhaust Collar:
  - Type 304 s/s - 10" Dia. Round (standard)
  - Other \_\_\_\_\_
- Interior Lighting Fixture:
  - Standard T-8 Fluorescent Light w/ Bulbs (standard)
- Sash Frame & Glass:
  - Painted Steel Frame w/ Laminated Safety Glass (standard)
- Airfoil:
  - Painted Steel (standard)
- Trough:
  - Painted Steel (standard)
- Sash Type:
  - Vertical Rising (standard)
- Also Included:
  - Sash Stop
  - 1/4 Turn Ball Valves (Specify fixture type)
- Pre-Piping / Pre-Wiring (standard)
  - All hoods are pre-wired & pre-piped unless other wise noted

### OPTIONAL ACCESSORIES

- Lattice Assembly:
  - None
  - Aluminum Lattice Assembly
  - Stainless Steel Assembly
- Interior Lighting Fixture:
  - Explosion Proof Light
- Fume Hood Alarm Option:
  - Alarm Tel AFA\_500
  - Alarm Tel AFA\_1000
  - VAV Cutout (Customer provides cutout detail)
- Sash Type:
  - Combination Sash (Horizontal/Vertical)
- Sash Frame & Glass:
  - Stainless Steel Frame w/ Laminated Safety Glass
- Airfoil:
  - Type 304 Stainless Steel
  - Type 316 Stainless Steel
- Trough:
  - None (standard)
  - Type 304 Stainless Steel
  - Type 316 Stainless Steel
- Ceiling Enclosure:
  - Yes (Height: \_\_\_\_\_)  No
- Finished Back:
  - Yes  No
- Sash Type:
  - Vertical
  - Horizontal
  - Interlocking
- Additional Options:
 

Blower	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Push Button Sash	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Automatic Sash Return	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Sash Locks	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Sash Alarm	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Add Air Chamber	<input type="checkbox"/> Yes	<input type="checkbox"/> No



- LINER MATERIAL:**
- POLY RESIN (standard)
  - 304 STAINLESS STEEL
  - 316 STAINLESS STEEL
  - PVC
  - SPECIAL

- PRE-PIPING:**
- PRE-PIPING UP
  - PRE-PIPING DOWN

- BY-PASS TYPE:**
- OPEN BY-PASS (standard)
  - RESTRICTED BY PASS

**SPECIAL REQUIREMENTS:** \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**HOOD BASE CABINETS:**

	Width	29"	*33"	35"
Acid				
Flammable				
Vacuum Pump				
Standard Base				
Table Base Frame				

\*PLEASE NOTE: WHERE ADA HT IS REQUIRED A LOW HT CABINET IS PROVIDED WITH A WELDED 4" APRON TO ACCOMMODATE FOR ADA HT. REQUIREMENT.

- WORK SURFACE MATERIAL:**
- Epoxy (standard)
  - Type 304 Stainless Steel
  - Type 316 Stainless Steel

Custom sizes available; contact AMS



**AIR MASTER SYSTEMS CORP**  
 6480 NORTON CENTER DR.  
 MUSKEGON, MI. 49441  
 TEL (231) 798-1111 FAX (231) 798-4000  
 www.airmastersystems.com

CUSTOMER

Qty: \_\_\_\_\_  
 Color: \_\_\_\_\_