


The Eliminator 700 SERIES: Perchloric Acid Hood

APPLICATIONS: K-12; University; Private Industry

NOTE: DRAIN IN TROUGH WILL BE SAME SIDE AS WASH DOWN VALVE

(SPECIFY FITTINGS / FIXTURES REQUIRED)

(SPECIFY FITTINGS / FIXTURES REQUIRED)



PRE-PIPING:

- PRE-PIPING UP
- PRE-PIPING DOWN

BY-PASS TYPE:

- OPEN BY-PASS (standard)
- RESTRICTED BY PASS

- 36"
- 48"
- 60"
- 72"
- 96"

SPECIAL REQUIREMENTS: _____

- Interior Baffle Arrangement:
 Fixed baffles (standard)
- Exhaust Collar:
 S/S - 10"Dia. Round (standard)
 Other _____
- Interior Lighting Fixture:
 Standard T-8 Fluorescent Light w/ Bulbs (standard)
- Sash Frame & Glass:
 Type 316 S/S Frame w/ Laminated Safety Glass (standard)
- Sash Type:
 Vertical (standard)
- Additional Features:
 Type 316 Stainless Steel Airfoil & Built-In Trough
 Type 316 Stn. Stl. one-piece integral work surface & liner
- Pre-Piping / Pre-Wiring (standard)
 All hoods are pre-wired & pre-piped unless other wise noted

OPTIONAL ACCESSORIES

- Lattice Assembly:
 None
 Aluminum Lattice Assembly
 Stainless Steel Assembly
- Interior Lighting Fixture:
 Explosion Proof Light
- Fume Hood Alarm Option:
 Alarm Tel AFA_500
 Alarm Tel AFA_1000
 VAV Cutout (Customer provides cutout detail)
- Sash Type:
 Combination Sash (Horizontal/Vertical)
- Sash Frame & Glass:
 Stainless Steel Frame w/ Laminated Safety Glass
- Ceiling Enclosure:
 Yes (Height: _____) No
- Finished Back:
 Yes No
- Sash Type:
 Combination (H/V)
 Interlocking
- Additional Options:
- | | | |
|-----------------------|------------------------------|-----------------------------|
| Blower | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Automatic Sash Return | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Sash Locks | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Sash Alarm | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
| Add Air Chamber | <input type="checkbox"/> Yes | <input type="checkbox"/> No |
- Service Fixture Type:
 Ball valve w/ lever handle (ADA)
 Needle valve w/ four-arm handle
 Front Mount Valve



AIR MASTER SYSTEMS CORP

6480 NORTON CENTER DR.
 MUSKEGON, MI. 49441
 TEL. (231) 798-1111 FAX (231) 798-4000
 www.airmastersystems.com

CUSTOMER

Qty: _____
 Color: _____