

# TITAN

---

## S E R I E S



## Heavy Duty Mobile Cart

The versatility of the Heavy Duty Mobile Cart will handle a multitude of functions in the laboratory environment. From research facilities to teaching labs, there is a configuration to meet your requirements. The cart shelving is adjustable in 1" increments of for optimum flexibility to fit your ever changing needs.

This cart is designed to handle heavy instruments not easily transported. The heavy duty casters allow for mobility needed to transport equipment from place to place.

so adaptable.

- Cart height 36" & 78"
- Shelves and work surfaces are adjustable in 1" increments
- 700 lbs locking casters with or without levelers
- Leveling casters are an optional feature
- Capabilities to attach hanging cabinets
- Widths available in : 36" 48" 60" 72"
- maximum weight capacity is 2800 lbs . load capacity breakdown is as follows. shelves 600 lbs, table frame 1000 lbs, lower shelf 600 lbs,
- Standard working height is 36" with a depth of 30"
- Counter top options available in epoxy and phenolic
- Standard shelves are 18" deep.

Air Master Systems Corp.

6480 Norton Center Drive, Muskegon, MI 49441 • 231.798.1111 • 231.798.4000 fax • [airmastersystems.com](http://airmastersystems.com)



# Ordering Guide

## Heavy Duty Mobile Cart

**Unit Options:**

- CT cart
- SF shelf
- TF table frames
- LS lower shelf

**Top:**

- E Epoxy
- P Phenolic
- N None

HDMC-CT-723078-CL

**Width:**

36"  
48"  
60"  
72"

**Depth:**

18"  
(shelf)  
30"  
(cart)

**Height:**


36"  
78"  
(cart)

**Cart:**

- C 700# caster
- CL Casters & Levelers
- X 1000# leveler caster
- CM cable management (optional)

- HDMC-CT-723078-CL
- HDMC-SF-7218-E (2)
- HDMC-TF-7230-E
- HDMC-LS-7230-E

(model as shown with all four options installed)



**Air Master Systems Corp.**

6480 Norton Center Drive, Muskegon, MI 49441 • 231.798.1111 • 231.798.4000 fax • [airmastersystems.com](http://airmastersystems.com)