

Gemini System

36.0" 42.0" and 48.0" width units

SPECIFICATION WORKSHEET - SINGLE OR DOUBLE-SIDED, TWIN-TUBE SUPPORT FRAME

RACEWAY ELECTRIC (Per Side *)

1-120VAC, 20Amp GFI Duplex

1-USB x4 charging outlet

1-20 Amp GFI / USB x2 combo

RACEWAY ELECTRIC (Per Side *)

1-120VAC, 20Amp GFI Duplex

1-USB x4 charging outlet

1-20 Amp GFI / USB x2 combo

RACEWAY ELECTRIC (Per Side *)

1-120VAC, 20Amp GFI Duplex

1-USB x4 charging outlet

1-20 Amp GFI / USB x2 combo

DATA PORTS (Left Post)
RJ45 (shielded)

x1

x2

LOWER ELECTRIC (Left Post)
pre-wired as selected below

1-120VAC, 20Amp GFI Duplex

1-USB charging outlet

1-20 Amp GFI / USB x2 combo

NONE

WIDTH OF UNIT

36.0"

42.0"

48.0"

ELECTRICAL OUTLET COLOR

WHITE

BLACK

MODESTY PANEL (non-adjustable) height at 36"

Yes

No

Single-Sided Bench

Double-Sided Bench

30.0" Table Depth (per side)

36.0" Table Depth (per side)

BENCH FRAME COLOR _____

Swap Sides, Plumbing on Left, Data on Right on GS1

POWER

PLUMBING

- ELECTRICAL WIRING OPTIONS**
- Single-circuit 20A, 1- Straight-blade plug, 1 cord
 - Double-circuit 20A, 2- Straight-blade plugs, 2 cords
 - Single-circuit 20A, 1- Twist-lock plug, L5-20P (1x120v), 1 cord
 - Double-circuit 20A, 2- Twist-lock plugs, L5-20P (2x120v), 2 cords
 - Double-circuit 20A, 1- Twist-lock plug, L14-20P (1x240v), 1 cord
 - Quad-circuit 20A, 2- Twist-lock plugs, L14-20P (2x240v), 2 cords
- Specified length (top of bench to ceiling) _____

- ELECTRICAL PIGTAIL EXIT LOCATION**
- | | |
|--|---|
| <input type="checkbox"/> Side above plumbing services | <input type="checkbox"/> GS1- Single-Sided: |
| <input type="checkbox"/> Side opposite plumbing services | <input type="checkbox"/> Exiting LEFT-SIDE |
| <input type="checkbox"/> Exit post near bottom | <input type="checkbox"/> Exiting RIGHT-SIDE |


SEE PLUMBING ILLUSTRATION ON SECOND PAGE FOR HOSE TYPES, CONNECTIONS, and FIXTURES. PLEASE SPECIFY IF BURNING GAS, OXYGEN, OR CLEAN GAS REQUIRED

-Your notes:

GS - - - P D E

example: **GS2-4830-3P2D7E** 3 Plumbed Services
 Single-Sided = 1 2 Data Ports
 Double-Sided = 2 48" wide X 30" deep (per side) 7 Electrical (inc lower)

(Sold Separately)



A
 Fixture to Ceiling
 1/2" or 5/8" PVC HOSE
 or 1/2" Braided Stainless
 w/ Male Quick-Connect


B
 Internal
 1/2" PVC HOSE
 or 1/2" Braided Stainless

Junction
 Keyed set Male/ Female
 Quick-Connects


External
 1/2" or 5/8" PVC HOSE
 w/ Male Quick-Connect

C
 Internal
 1/2" PVC HOSE
 or 1/2" Braided Stainless
 w/ Male Quick-Connect


	A	B	C
PVC HOSE	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
STAINLESS HOSE	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>




Single, Needle Valve - Black Knob
 FINISH: _____




Single, Needle Valve - ColorTech Knob
 FINISH: _____



Double, Needle Valve - Black Knobs
 FINISH: _____



Double, Needle Valve - ColorTech Knobs
 FINISH: _____



Double, Ball Valve - Lever Handles
 FINISH: _____

(There is a maximum of three fixtures per side)