

# Gemini System

36.0" 42.0" and 48.0" width units

SPECIFICATION WORKSHEET - SINGLE OR DOUBLE-SIDED, TWIN-TUBE SUPPORT FRAME

**RACEWAY ELECTRIC (Per Side \*)**

1-120VAC, 20Amp GFI Duplex

1-USB x4 charging outlet

1-20 Amp GFI / USB x2 combo

**DATA PORTS (Left Post)**  
RJ45 (shielded)

x1

x2

**LOWER ELECTRIC (Left Post)**  
pre-wired as selected below

1-120VAC, 20Amp GFI Duplex

1-USB charging outlet

1-20 Amp GFI / USB x2 combo

NONE

**WIDTH OF UNIT**

36.0"

42.0"

48.0"

**RACEWAY ELECTRIC (Per Side \*)**

1-120VAC, 20Amp GFI Duplex

1-USB x4 charging outlet

1-20 Amp GFI / USB x2 combo

---

**ELECTRICAL OUTLET COLOR**

WHITE

BLACK

\_\_\_\_\_

---

\*\*Per Side when GS2 Double-Sided Bench is selected.

---

BENCH FRAME COLOR \_\_\_\_\_

**RACEWAY ELECTRIC (Per Side \*)**

1-120VAC, 20Amp GFI Duplex

1-USB x4 charging outlet

1-20 Amp GFI / USB x2 combo

**SERVICES (Right Post)**  
pre-plumbed to top of upright  
*see below for options*

\_\_\_\_\_ indicate service required

\_\_\_\_\_ indicate service required

\_\_\_\_\_ indicate service required

**MODESTY PANEL**  
(non-adjustable) height at 36"

Yes

No

Single-Sided Bench

Double-Sided Bench

30.0" Table Depth (per side)

36.0" Table Depth (per side)

Swap Sides, Plumbing on Left, Data on Right on GS1

## POWER

## PLUMBING

- ELECTRICAL WIRING OPTIONS**
- Single-circuit 20A, 1- Straight-blade plug, 1 cord
  - Double-circuit 20A, 2- Straight-blade plugs, 2 cords
  - Single-circuit 20A, 1- Twist-lock plug, L5-20P (1x120v), 1 cord
  - Double-circuit 20A, 2- Twist-lock plugs, L5-20P (2x120v), 2 cords
  - Double-circuit 20A, 1- Twist-lock plug, L14-20P (1x240v), 1 cord
  - Quad-circuit 20A, 2- Twist-lock plugs, L14-20P (2x240v), 2 cords
- Specified length (top of bench to ceiling) \_\_\_\_\_
- ELECTRICAL PIGTAIL EXIT LOCATION**
- |  |   |
|--|---|
| <input type="checkbox"/> Side above plumbing services    | <input type="checkbox"/> GS1- Single-Sided: |
| <input type="checkbox"/> Side opposite plumbing services | <input type="checkbox"/> Exiting LEFT-SIDE  |
| <input type="checkbox"/> Exit post near bottom           | <input type="checkbox"/> Exiting RIGHT-SIDE |

SEE PLUMBING ILLUSTRATION ON SECOND PAGE FOR HOSE TYPES, CONNECTIONS, and FIXTURES. PLEASE SPECIFY IF BURNING GAS, OXYGEN, OR CLEAN GAS REQUIRED

-Your notes:

\_\_\_\_\_


\_\_\_\_\_

\_\_\_\_\_

**GS \_ - \_ - \_ P \_ D \_ E**

**example:** **GS2-4830-3P2D7E** 3 Plumbed Services  
 Single-Sided = 1 2 Data Ports  
 Double-Sided = 2 48" wide X 30" deep (per side) 7 Electrical (inc lower)

(Sold Separately)



**A**  
 Fixture to Ceiling  
 1/2" or 5/8" PVC HOSE  
 or 1/2" Braided Stainless  
 w/ Male Quick-Connect


**B**  
 Internal  
 1/2" PVC HOSE  
 or 1/2" Braided Stainless

Junction  
 Keyed set Male/ Female  
 Quick-Connects


External  
 1/2" or 5/8" PVC HOSE  
 w/ Male Quick-Connect

**C**  
 Internal  
 1/2" PVC HOSE  
 or 1/2" Braided Stainless  
 w/ Male Quick-Connect


|                |                          |                          |                          |
|----------------|--------------------------|--------------------------|--------------------------|
|                | A                        | B                        | C                        |
| PVC HOSE       | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| STAINLESS HOSE | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |




Single, Needle Valve - Black Knob  
 FINISH: \_\_\_\_\_




Single, Needle Valve - ColorTech Knob  
 FINISH: \_\_\_\_\_



Double, Needle Valve - Black Knobs  
 FINISH: \_\_\_\_\_



Double, Needle Valve - ColorTech Knobs  
 FINISH: \_\_\_\_\_



Double, Ball Valve - Lever Handles  
 FINISH: \_\_\_\_\_

(There is a maximum of three fixtures total, on one side only)