


THE ELIMINATOR 600 SERIES: RADIO ISOTOPE HOOD

APPLICATIONS: K-12; UNIVERSITY; PRIVATE INDUSTRY

**SPECIFY FITTINGS/
FIXTURES REQUIRED:**

**SPECIFY FITTINGS/
FIXTURES REQUIRED:**



STANDARD FEATURES

- **Interior Baffle Arrangement:**
Fixed Baffles (standard)
- **Exhaust Collar:**
Type 304 s/s - 10" Dia. Round (standard)
Other _____
- **Interior Lighting Fixture:**
Standard LED Fixtures (standard)
- **Sash Frame & Glass:**
304 Stainless Steel Frame w/ Laminated Safety Glass (standard)
- **Sash Type:**
Vertical (standard)
- **Additional Features:**
304 Stainless Steel Airfoil
Type 304 Stn. Stl. one-piece integral work surface & liner

OPTIONAL ACCESSORIES

- Lattice Assembly:**

- Interior Lighting Fixture:**
 Explosion Proof Light
- Sash Frame & Glass:**
 316 Stainless Steel Frame w/ Laminated Safety Glass
- Airfoil:**

- Trough:**

- Ceiling Enclosure:**
 Yes (Height: _____) No
• **Hinged Access Panel on Ceiling Enclosures?:**
 Yes No
- Finished Back:**
 Yes No
- Sash Type:**

- Additional Options:**
 Blower Sash Alarm
 Automatic Sash Return Sash Locks
 Sash Cap Sash Stops
- Service Fixture Type:**

- Pre-Wire + Electrical:**
 Yes No




AIR MASTER SYSTEMS CORP
6480 NORTON CENTER DR.
MUSKEGON, MI 49441
TEL: (231) 798-1111 FAX: (231) 798-4000


CUSTOMER: _____
QTY: _____
COLOR: _____

THE ELIMINATOR 600 SERIES: RADIO ISOTOPE HOOD

APPLICATIONS: K-12; UNIVERSITY; PRIVATE INDUSTRY



**SERVICE ROUGH-IN DETAIL
(FOR 36", 48", 60", & 72" HOODS)**



**SERVICE ROUGH-IN DETAIL
(FOR 96" HOODS)**


EPOXY WORK SURFACE (SPECIFY LOCATION OF CUPSINK)

*When utilizing 22" deep cabinets, cup sink measurement changes to 5.25" in lieu of 9"


Note: This location will have the cup sink sit slightly under the baffle.

Include Acid Cabinet Vent Hole

CUPSINK OPTIONS:



WORK SURFACE MATERIAL:



PART #	SUPERSTRUCTURE			WORKTOP			HOOD OPENING		CFM @29"	CFM @18"	SP	BLOWER MODEL #	SHIPPING WEIGHT
	W	H	D	W	H	D	W	H					
EH-613-36	36"	56"	33.25"	36"	1.25"	30"	26"	29"	567	354	.25	BL-600-36	190
EH-613-48	48"	56"	33.25"	48"	1.25"	30"	38"	29"	833	521	.25	BL-600-48	560
EH-613-60	60"	56"	33.25"	60"	1.25"	30"	50"	29"	1100	688	.34	BL-600-60	675
EH-613-72	72"	56"	33.25"	72"	1.25"	30"	62"	29"	1367	854	.34	BL-600-72	790
EH-613-96	96"	56"	33.25"	96"	1.25"	30"	86"	29"	1900	1188	.40	BL-600-96	1025