

# Gemini System

36.0" 42.0" and 48" width units

SPECIFICATION WORKSHEET - SINGLE OR DOUBLE SIDED, TWIN-TUBE SUPPORT FRAME

## RACEWAY ELECTRIC

Reference page 3 for selection sheet

## DATA PORTS (Left Post)

RJ45 (shielded)


## LOWER ELECTRIC (Left Post)

pre-wired as selected below

## WIDTH OF UNIT

\*60" & 72" can be manufactured with single full width shelf

## TABLE DEPTH



## SERVICES (Right Post)

pre-plumbed to top of upright  
see below for options

indicate service required

indicate service required

indicate service required

## MODESTY PANEL

(non-adjustable) height at 36"

Yes ☐

No ☐

## SINGLE / DOUBLE SIDED

☐ Swap Sides, Plumbing on Left, Data on Right on GS1


## POWER

ELECTRICAL WIRING OPTIONS

Specified length (top of bench to ceiling) \_\_\_\_\_

ELECTRICAL PIGTAIL EXIT LOCATION

GS1- Single-Sided:



5-20P



L5-20P



L14-20P

## DATA

DATA EXIT LOCATION

Specified length (top of bench to ceiling) \_\_\_\_\_

## PLUMBING

SEE PLUMBING ILLUSTRATION ON SECOND PAGE FOR HOSE TYPES, CONNECTIONS, and FIXTURES. PLEASE SPECIFY IF BURNING GAS, OXYGEN, OR CLEAN GAS REQUIRED

-Your notes:


GS \_ - \_ - \_ P \_ D \_ E

example: GS2-4830-3P2D7E

Single-Sided = 1 / Double-Sided = 2

7 Electrical (Including lower)

2 Data Ports



## GEMINI BENCH




example: **GSK-1535-12B**

Width      Bracket Depth

Full Shelf Depth\*      Bracket Type

\*This includes the 1" OR 3" extensions added to back edge.

## Bracket Types:



Gemini Width	Does it allow split shelving?	Part Numbers Needed
36"	Single Only	GSK-__34-__
42"	Single Only	GSK-__40-__
48"	Single Only	GSK-__46-__

### Bracket Selection:

### Bracket ONLY:

☐ ONLY Brackets\*  
No shelves

\*For purchases with alternative shelving.



Depth\*: \*Depth w/o extensions

### Shelf Depth Extension:

### Quantity:


### Width:

### Comments:

Specified length (top of bench to ceiling) \_\_\_\_\_

(Sold Separately)



PVC HOSE ☐


STAINLESS HOSE ☐

A ☐ B ☐ C ☐


A ☐ B ☐ C ☐

(There is a maximum of three fixtures per side)


\*4 Plumbing services maximum through Broen. Please contact AMS for Broen product information.




Single, Needle Valve - Black Knob  
☐ FINISH: \_\_\_\_\_




Single, Needle Valve - ColorTech Knob  
☐ FINISH: \_\_\_\_\_



Double, Needle Valve - Black Knobs  
☐ FINISH: \_\_\_\_\_



Double, Needle Valve - ColorTech Knobs  
☐ FINISH: \_\_\_\_\_



Double, Ball Valve - Lever Handles  
☐ FINISH: \_\_\_\_\_

# Gemini System

36.0" 42.0" and 48" width units

## SPECIFICATION WORKSHEET - RACEWAY ELECTRICAL



- Black
- Brown
- Gray
- Ivory
- Light Almond
- White

120VAC, 20Amp GFI Duplex



- Black
- Gray
- Ivory
- Light Almond
- White

USB Charging Outlet



- Black
- Brown
- Gray
- Ivory
- Red
- Light Almond
- White

20Amp GFI / USB x2 Combo



RJ45 Shielded Coupler

Gemini Systems 36.0", 42.0", and 48.0" Width Units allow **up to 4** GFI slots on the top raceway.

Double-sided allow **up to 8** GFI slots on the top raceway.

From Left to Right facing front of bench. ➤

SLOT 1

SLOT 2

SLOT 3

SLOT 4

**Slots 5 - 8 (Use only if Double-Sided)**

SLOT 5

SLOT 6

SLOT 7

SLOT 8