

Stainless Steel Laboratory Casework

product catalog




Contents

Cabinet Specifications	3
Ordering Guide	4
Options Menu	5
Base Cabinets — Standing Height	8
Open Cabinets	8
Door Cabinets	8
Drawer Cabinets	9
Door & Drawer Cabinets	12
Shelves & Pull Out Shelves	16
Corner Cabinets	16
Three or Five Equal Drawer Cabinet	17
Sink & Fume Hood Cabinets	18
Base Cabinets — Sitting Height	19
Open Cabinets	19
Door Cabinets	19
Drawer Cabinets	20
Door & Drawer Cabinets	23
Sink & Fume Hood Cabinets	25
Corner Cabinets	26
Base Cabinets — ADA Height	27
Open Cabinets	27
Door Cabinets	27
Drawer Cabinets	28
Door & Drawer Cabinets	29
Corner Cabinets	30
Sink & Fume Hood Cabinets	31
Base Cabinets — Corrosives & Flammables	33
Corrosives Storage Cabinets	33
Flammables Storage Cabinets	36
Vacuum Pump Cabinets	38
Mobile Cabinets	40
Suspended Cabinets	42
Overlay Cabinets	44
Tall Cabinets	46
Wall Cabinets — 24" Height	48
Wall Cabinets — 30" Height	50
Wall Cabinets — 36" Height	52
Wall Cabinets — 48" Height	54
Aprons	56
Utility Tables	57
Filler Panels	58
Knee Space Rails and Panels	60
Counter Island Reagent	62
Accessories	63
Canopy Hoods	64
Countertop Wall Brackets	64
Countertops	65
Table Tops	66
Wall Mounted Sink Enclosure	67
Scullery Sinks	68
Cylinder Racks	69
Umbilicals	70

Note:

Images in this catalog are for representation only. For complete product specifications, contact your AMS representative.

Cabinet Specifications



All surfaces shall be Stainless Steel Type 304. Exterior finish is #4 and interior finish is 2B. Specific test data is available upon request.

Gauges of base cabinet components

Aprons	18 ga.
Back Panels	18 ga.
Bottom Panels	18 ga.
Drawer/Door Outer Panel	20 ga.
Door Inner Panel	20 ga.
Drawer Bodies	20 ga.
Legs - 2"sq. tube	16 ga.
Shelves	18 ga.
Side Panels	18 ga.
Table Frames	18 ga.
Security Panels	20 ga.

For complete specifications, request our full line of product specifications. Subject to change.

IMPORTANT NOTE:

All part numbers may not reflect actual dimensions.

Ordering Guide

S		TYPE	HEIGHT	WIDTH	-	STYLE	-	OPTIONS
Cabinet Type								
B	Base (Inset)							
F	Floor (Tall)							
HB	Suspended Base							
OB	Overlay Base							
MB	Mobile Base							
W	Wall							
ASB	Acid Storage Base							
FSB	Flammable Storage							
Base Cabinet Height								
29"	Sitting Height							
32"	ADA Height							
35"	Standing Height							
Wall Cabinet Height								
24" 30" 36" 48"								
Cabinet Width								
12" 15" 18" 21" 24" 30" 36" 42" 48"								
Cabinet Style								
Reference AMS catalog for door and drawer configurations								
Option Suffix (for base cabinet)								
-18	Depth(contact AMS for approval)			-PS		1 Pull Out Shelf		
-LK	Locks			-2PS		2 Pull Out Shelves		
-LTK	Toe Kick On Left Side Of Cabinet			-3PS		3 Pull Out Shelves		
-RTK	Toe Kick On Right Side Of Cabinet			-VK		Vent Kit (only for acid cabinet)		
-SS	Stainless Steel, Type 304			-VD		Vertical Divider		
-SS316	Stainless Steel, Type 316			-FB		Finished Backs		
-BB	Bread Board(selected models only)			-EC		End Covers		
				-F		Flush Pull		
Option Suffix (for wall cabinet)								
-16	16" Depth			-TG		Tempered Glass		
-LK	Locks			-LG		Laminate (Safety) Glass		
-SS	Stainless Steel, Type 304			-THB		Titan Hanging Bracket		
-SS316	Stainless Steel, Type 316			-F		Flush Pull		
Option Suffix (for mobile cabinet)								
-18	Depth(contact AMS for approval)			-T		3/4" THK. Metal Top		
-LK	Locks			-E		3/4" THK. Epoxy Top		
-SS	Stainless Steel, Type 304			-P		3/4" THK. Phenolic Top		
-SS316	Stainless Steel, Type 316			-2"		Casters (in lieu of standard 3")		
-BB	Bread Board (selected models only)			-F		Flush Pull		
Option Suffix (for tall cabinet)								
-LK	Locks			-LG		Laminate (Safety) Glass		
-SS	Stainless Steel, Type 304			-WB		Wardrobe Bracket		
-SS316	Stainless Steel, Type 316			-PS		1 Pull Out Shelves		
-TG	Tempered Glass			-F		Flush Pull		

Options Menu

DOOR & DRAWER FRONTS



Inset Steel
(Standard Default)



Flush Overlay
Steel

DRAWER SLIDES



- ☐ 100 lb. Full Extension (Standard Default)
- ☐ 100 lb. Full Extension Soft Close
- ☐ 200 lb. Full Extension

HINGES



Inset Steel
Stainless Steel 5-Knuckle



Flush Overlay
Stainless Steel 5-Knuckle

Options Menu

MAGNETIC & ROLLER CATCH



Magnetic Catch
(Standard Default)



Roller Catch

LOCKS

Master 5-Pin
(Standard Default)



PULLS



☐ Stainless Steel Wire Pulls

Options Menu

LABEL HOLDERS



- ☐ None (Standard Default)
- ☐ One each door and drawer
- ☐ One per pair of doors and each drawer

BASE UNIT SHELF



- ☐ One piece adjustable (Standard Default)
Extra Quantity Available

GLASS FRAMED DOORS



- ☐ 3/16" Clear Float Glass (Standard Default)
- ☐ 1/4" Laminated Safety Glass
- ☐ 1/4" Tempered Safety Glass



- ☐ Without Security Panel (Standard Default)
- ☐ With Security Panel

BREADBOARD



- ☐ Sitting Height variations only

GLASS DOORS



- ☐ 1/4" Clear Float Glass (Standard Default)
- ☐ Laminated Glass
- ☐ Tempered Glass

Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Door units feature a single shelf as standard. All open units have 2 shelves. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



OPEN CABINETS



Up to 24" wide

Width	Item Number
12"	SB3512-000-SS
15"	SB3515-000-SS
18"	SB3518-000-SS
21"	SB3521-000-SS
24"	SB3524-000-SS



30" and wider

Width	Item Number
30"	SB3530-000-SS
36"	SB3536-000-SS
42"	SB3542-000-SS
48"	SB3548-000-SS

DOOR CABINETS



Example: SB3515-100R-SS

Width	Left & Right Hinged
12"	SB3512-100L-SS
15"	SB3515-100L-SS
18"	SB3518-100L-SS
21"	SB3521-100L-SS
24"	SB3524-100L-SS
12"	SB3512-100R-SS
15"	SB3515-100R-SS
18"	SB3518-100R-SS
21"	SB3521-100R-SS
24"	SB3524-100R-SS

Note: L(left) or R(right) refers to door hinge side.



30" and wider

Width	Item Number
30"	SB3530-200-SS
36"	SB3536-200-SS
42"	SB3542-200-SS
48"	SB3548-200-SS



Example: SB3530-100-LPNL-SS

Width	Left & Right Hinged
30"	SB3530-100-LPNL-SS
36"	SB3536-100-LPNL-SS
42"	SB3542-100-LPNL-SS
48"	SB3548-100-LPNL-SS
30"	SB3530-100-RPNL-SS
36"	SB3536-100-RPNL-SS
42"	SB3542-100-RPNL-SS
48"	SB3548-100-RPNL-SS

Note: L(left) or R(right) refers to door hinge side. PNL(panel) refers to panel next to door

Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB3542-034-LK**

Height: 34.75"
Doors: 2
Half Drawers: 2
Width: 42"
Drawers: 4
Lock: LK



DRAWER CABINETS



Up to 24" wide

Width	Item Number
12"	SB3512-032-SS
15"	SB3515-032-SS
18"	SB3518-032-SS
21"	SB3521-032-SS
24"	SB3524-032-SS



30" and wider

Width	Item Number
30"	SB3530-032-SS
36"	SB3536-032-SS
42"	SB3542-032-SS
48"	SB3548-032-SS



Up to 24" wide

Width	Item Number
12"	SB3512-024-SS
15"	SB3515-024-SS
18"	SB3518-024-SS
21"	SB3521-024-SS
24"	SB3524-024-SS



30" and wider

Width	Item Number
30"	SB3530-028-SS
36"	SB3536-028-SS
42"	SB3542-028-SS
48"	SB3548-028-SS



Up to 24" wide

Width	Item Number
12"	SB3512-008-SS
15"	SB3515-008-SS
18"	SB3518-008-SS
21"	SB3521-008-SS
24"	SB3524-008-SS



30" and wider

Width	Item Number
30"	SB3530-008-SS
36"	SB3536-008-SS
42"	SB3542-008-SS
48"	SB3548-008-SS



Up to 24" wide

Width	Item Number
12"	SB3512-040-SS
15"	SB3515-040-SS
18"	SB3518-040-SS
21"	SB3521-040-SS
24"	SB3524-040-SS



30" and wider

Width	Item Number
30"	SB3530-040-SS
36"	SB3536-040-SS
42"	SB3542-040-SS
48"	SB3548-040-SS

Base Cabinets — Standing Height


Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB3542-0(10)4-LK** — Lock

Height: 34.75"
Doors: 2
Half Drawers: 1
Width: 42"
Drawers: 4

Example: **SB3542-034-1F-LK** — Lock

Height: 34.75"
Doors: 2
Half Drawers: 1
Width: 48"
Drawers: 3
File Drawers: 1



DRAWER CABINETS



Width	Item Number
12"	SB3512-012-1F-SS
15"	SB3515-012-1F-SS
18"	SB3518-012-1F-SS
21"	SB3521-012-1F-SS
24"	SB3524-012-1F-SS

Up to 24" wide



Width	Item Number
12"	SB3512-000-2F-SS
15"	SB3515-000-2F-SS
18"	SB3518-000-2F-SS
21"	SB3521-000-2F-SS
24"	SB3524-000-2F-SS

Up to 24" wide



Width	Item Number
12"	SB3512-004-1F-SS
15"	SB3515-004-1F-SS
18"	SB3518-004-1F-SS
21"	SB3521-004-1F-SS
24"	SB3524-004-1F-SS

Up to 24" wide



Up to 24" wide



Width	Item Number
30"	SB3530-024-2F-SS
36"	SB3536-024-2F-SS
42"	SB3542-024-2F-SS
48"	SB3548-024-2F-SS

30" and wider



Width	Item Number
30"	SB3530-000-4F-SS
36"	SB3536-000-4F-SS
42"	SB3542-000-4F-SS
48"	SB3548-000-4F-SS

30" and wider



Width	Item Number
30"	SB3530-008-2F-SS
36"	SB3536-008-2F-SS
42"	SB3542-008-2F-SS
48"	SB3548-008-2F-SS

30" and wider



Width	Item Number
30"	SB3530-040-2F-SS
36"	SB3536-040-2F-SS
42"	SB3542-040-2F-SS
48"	SB3548-040-2F-SS


30" and wider

Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB3542-034-LK** — Lock

Height: 34.75"
Doors: 034
Half Drawers: LK
Width: 42"
Drawers: 034



DRAWER CABINETS



30" and wider

Width	Item Number
30"	SB3530-050-SS
36"	SB3536-050-SS
42"	SB3542-050-SS
48"	SB3548-050-SS



30" and wider

Width	Item Number
30"	SB3530-064-SS
36"	SB3536-064-SS
42"	SB3542-064-SS
48"	SB3548-064-SS



30" and wider

Width	Item Number
30"	SB3530-060-SS
36"	SB3536-060-SS
42"	SB3542-060-SS
48"	SB3548-060-SS



30" and wider

Width	Item Number
30"	SB3530-070-SS
36"	SB3536-070-SS
42"	SB3542-070-SS
48"	SB3548-070-SS



30" and wider

Width	Item Number
30"	SB3530-080-SS
36"	SB3536-080-SS
42"	SB3542-080-SS
48"	SB3548-080-SS



30" and wider

Width	Item Number
30"	SB3530-052-SS
36"	SB3536-052-SS
42"	SB3542-052-SS
48"	SB3548-052-SS



30" and wider

Width	Item Number
30"	SB3530-044-SS
36"	SB3536-044-SS
42"	SB3542-044-SS
48"	SB3548-044-SS




30" and wider

Width	Item Number
30"	SB3530-034-SS
36"	SB3536-034-SS
42"	SB3542-034-SS
48"	SB3548-034-SS

Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



Example: **SB3542-00(16)-LK** — Lock
 Height: 34.75"
 Doors: 00 (0 doors)
 Half Drawers: 16 (16 half drawers)
 Width: 42"
 Drawers: 00 (0 drawers)

Example: **SB3542-110R-LK** — Lock
 Height: 34.75"
 Doors: 110 (110 doors)
 Half Drawers: 00 (0 half drawers)
 Width: 42"
 Drawers: 16 (16 drawers)
 Right Hinged: R

DRAWER CABINETS



Up to 24" wide

Width	Item Number
12"	SB3512-016-SS
15"	SB3515-016-SS
18"	SB3518-016-SS
21"	SB3521-016-SS
24"	SB3524-016-SS



30" and wider

Width	Item Number
30"	SB3530-016-SS
36"	SB3536-016-SS
42"	SB3542-016-SS
48"	SB3548-016-SS



30" and wider

Width	Item Number
30"	SB3530-00(16)-SS
36"	SB3536-00(16)-SS
42"	SB3542-00(16)-SS
48"	SB3548-00(16)-SS

DOOR & DRAWER CABINETS



Right-hinged door shown

Width	Left & Right Hinged
12"	SB3512-110L-SS
15"	SB3515-110L-SS
18"	SB3518-110L-SS
21"	SB3521-110L-SS
24"	SB3524-110L-SS
12"	SB3512-110R-SS
15"	SB3515-110R-SS
18"	SB3518-110R-SS
21"	SB3521-110R-SS
24"	SB3524-110R-SS



Right-hinged door shown

Width	Left & Right Hinged
12"	SB3512-102L-SS
15"	SB3515-102L-SS
18"	SB3518-102L-SS
21"	SB3521-102L-SS
24"	SB3524-102L-SS
12"	SB3512-102R-SS
15"	SB3515-102R-SS
18"	SB3518-102R-SS
21"	SB3521-102R-SS
24"	SB3524-102R-SS

Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



Example: **SB3542-034-LK**

Height: 34.75"
Doors: 034
Half Drawers: LK
Width: 42"
Drawers: 2

DOOR & DRAWER CABINETS



30" and wider

Width	Item Number
30"	SB3530-212-SS
36"	SB3536-212-SS
42"	SB3542-212-SS
48"	SB3548-212-SS



30" and wider

Width	Item Number
30"	SB3530-220-SS
36"	SB3536-220-SS
42"	SB3542-220-SS
48"	SB3548-220-SS



30" and wider

Width	Item Number
30"	SB3530-202-A-SS
36"	SB3536-202-A-SS
42"	SB3542-202-A-SS
48"	SB3548-202-A-SS



30" and wider

Width	Left & Right Hinged
30"	SB3530-202L-SS
36"	SB3536-202L-SS
42"	SB3542-202L-SS
48"	SB3548-202L-SS
30"	SB3530-202R-SS
36"	SB3536-202R-SS
42"	SB3542-202R-SS
48"	SB3548-202R-SS



30" and wider

Width	Item Number
30"	SB3530-204-SS
36"	SB3536-204-SS
42"	SB3542-204-SS
48"	SB3548-204-SS



30" and wider

Width	Left & Right Hinged
30"	SB3530-142L-SS
36"	SB3536-142L-SS
42"	SB3542-142L-SS
48"	SB3548-142L-SS
30"	SB3530-142R-SS
36"	SB3536-142R-SS
42"	SB3542-142R-SS
48"	SB3548-142R-SS



30" and wider

Width	Left & Right Hinged
30"	SB3530-118L-SS
36"	SB3536-118L-SS
42"	SB3542-118L-SS
48"	SB3548-118L-SS
30"	SB3530-118R-SS
36"	SB3536-118R-SS
42"	SB3542-118R-SS
48"	SB3548-118R-SS



30" and wider


Width	Left & Right Hinged
30"	SB3530-118-AL-SS
36"	SB3536-118-AL-SS
42"	SB3542-118-AL-SS
48"	SB3548-118-AL-SS
30"	SB3530-118-AR-SS
36"	SB3536-118-AR-SS
42"	SB3542-118-AR-SS
48"	SB3548-118-AR-SS

Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB3542-034-LK** — Lock

Height: 34.75"
Doors: 034
Half Drawers: LK
Width: SB3542
Drawers: 034



DOOR & DRAWER CABINETS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB3530-140L-SS
36"	SB3536-140L-SS
42"	SB3542-140L-SS
48"	SB3548-140L-SS
30"	SB3530-140R-SS
36"	SB3536-140R-SS
42"	SB3542-140R-SS
48"	SB3548-140R-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB3530-140-BL-SS
36"	SB3536-140-BL-SS
42"	SB3542-140-BL-SS
48"	SB3548-140-BL-SS
30"	SB3530-140-BR-SS
36"	SB3536-140-BR-SS
42"	SB3542-140-BR-SS
48"	SB3548-140-BR-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB3530-134L-SS
36"	SB3536-134L-SS
42"	SB3542-134L-SS
48"	SB3548-134L-SS
30"	SB3530-134R-SS
36"	SB3536-134R-SS
42"	SB3542-134R-SS
48"	SB3548-134R-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB3530-134-AL-SS
36"	SB3536-134-AL-SS
42"	SB3542-134-AL-SS
48"	SB3548-134-AL-SS
30"	SB3530-134-AR-SS
36"	SB3536-134-AR-SS
42"	SB3542-134-AR-SS
48"	SB3548-134-AR-SS



30" and wider

Width	Item Number
30"	SB3530-210-SS
36"	SB3536-210-SS
42"	SB3542-210-SS
48"	SB3548-210-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB3530-210-AL-SS
36"	SB3536-210-AL-SS
42"	SB3542-210-AL-SS
48"	SB3548-210-AL-SS
30"	SB3530-210-AR-SS
36"	SB3536-210-AR-SS
42"	SB3542-210-AR-SS
48"	SB3548-210-AR-SS

Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If vertical divider is required between cup board and drawer section use the suffix "VD". If a finished back is required use the suffix "FB".



DOOR & DRAWER CABINETS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB3530-150L-SS
36"	SB3536-150L-SS
42"	SB3542-150L-SS
48"	SB3548-150L-SS
30"	SB3530-150R-SS
36"	SB3536-150R-SS
42"	SB3542-150R-SS
48"	SB3548-150R-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB3530-116L-SS
36"	SB3536-116L-SS
42"	SB3542-116L-SS
48"	SB3548-116L-SS
30"	SB3530-116R-SS
36"	SB3536-116R-SS
42"	SB3542-116R-SS
48"	SB3548-116R-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB3530-126L-SS
36"	SB3536-126L-SS
42"	SB3542-126L-SS
48"	SB3548-126L-SS
30"	SB3530-126R-SS
36"	SB3536-126R-SS
42"	SB3542-126R-SS
48"	SB3548-126R-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB3530-126-AL-SS
36"	SB3536-126-AL-SS
42"	SB3542-126-AL-SS
48"	SB3548-126-AL-SS
30"	SB3530-126-AR-SS
36"	SB3536-126-AR-SS
42"	SB3542-126-AR-SS
48"	SB3548-126-AR-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB3530-108L-SS
36"	SB3536-108L-SS
42"	SB3542-108L-SS
48"	SB3548-108L-SS
30"	SB3530-108R-SS
36"	SB3536-108R-SS
42"	SB3542-108R-SS
48"	SB3548-108R-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB3530-132L-SS
36"	SB3536-132L-SS
42"	SB3542-132L-SS
48"	SB3548-132L-SS
30"	SB3530-132R-SS
36"	SB3536-132R-SS
42"	SB3542-132R-SS
48"	SB3548-132R-SS

Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB3514-100-SVR**

Height: 34.75"
Doors: 100"
Half Drawers: SVR
Louver: Right
Width: 35"
Drawers: 100"
Sink: SVR



SHELVES & PULL OUT SHELVES



Width	Item Number	Actual Width
12"	SBS1222-SS	11 3/4"
15"	SBS1522-SS	14 3/4"
18"	SBS1822-SS	17 3/4"
21"	SBS2122-SS	20 3/4"
24"	SBS2422-SS	23 3/4"
30"	SBS3022-SS	29 3/4"
36"	SBS3622-SS	35 3/4"
42"	SBS4222-SS	41 3/4"
48"	SBS4822-SS	47 3/4"



Width	Item Number	Actual Width
15"	SBPOS-1522-SS	11 7/8"
18"	SBPOS-1822-SS	14 7/8"
21"	SBPOS-2122-SS	17 7/8"
24"	SBPOS-2422-SS	20 7/8"
30"	SBPOS-3022-SS	26 7/8"
36"	SBPOS-3622-SS	32 7/8"
42"	SBPOS-4222-SS	38 7/8"
48"	SBPOS-4822-SS	44 7/8"

CORNER CABINETS



Right-hinged door shown

Left Hinged	Right Hinged
SC3514-100L-SS	SC3514-100R-SS

Lazy Susan Add-On Suffix
"LS"




Right-hinged door shown

Left Hinged	Right Hinged
SC3514-100-SL-SS	SC3514-100-SR-SS


Left Hinged with Louvers	Right Hinged with Louvers
SC3514-100-SVL-SS	SC3514-100-SVR-SS

Lazy Susan Add-On Suffix
"LS"



Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



THREE OR FIVE EQUAL DRAWER CABINETS



Up to 24" wide

Width	Item Number
12"	SB3512-03X-SS
15"	SB3515-03X-SS
18"	SB3518-03X-SS
21"	SB3521-03X-SS
24"	SB3524-03X-SS



30" and wider

Width	Item Number
30"	SB3530-03X-SS
36"	SB3536-03X-SS
42"	SB3542-03X-SS
48"	SB3548-03X-SS



Up to 24" wide

Width	Item Number
12"	SB3512-05X-SS
15"	SB3515-05X-SS
18"	SB3518-05X-SS
21"	SB3521-05X-SS
24"	SB3524-05X-SS



30" and wider

Width	Item Number
30"	SB3530-05X-SS
36"	SB3536-05X-SS
42"	SB3542-05X-SS
48"	SB3548-05X-SS



30" and wider

Width	Item Number
30"	SB3530-0(10)X-SS
36"	SB3536-0(10)X-SS
42"	SB3542-0(10)X-SS
48"	SB3548-0(10)X-SS



Right-hinged door shown


Width	Left & Right Hinged
30"	SB3530-15X-L-SS
36"	SB3536-15X-L-SS
42"	SB3542-15X-L-SS
48"	SB3548-15X-L-SS
30"	SB3530-15X-R-SS
36"	SB3536-15X-R-SS
42"	SB3542-15X-R-SS
48"	SB3548-15X-R-SS

Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB3552-200-SI** — Island

Height: 34.75"
Doors: SB3552
Half Drawers: 200
Width: 52"
Drawers: SI
Sink



SINK & FUME HOOD CABINETS



Right-hinged door shown
Add "-SV(R/L)" for Vents

Width	Left & Right Hinged
12"	SB3512-100-SL-SS
15"	SB3515-100-SL-SS
18"	SB3518-100-SL-SS
21"	SB3521-100-SL-SS
24"	SB3524-100-SL-SS
12"	SB3512-100-SR-SS
15"	SB3515-100-SR-SS
18"	SB3518-100-SR-SS
21"	SB3521-100-SR-SS
24"	SB3524-100-SR-SS



30" and wider

Width	Item Number
30"	SB3530-200-S-SS
36"	SB3536-200-S-SS
42"	SB3542-200-S-SS
48"	SB3548-200-S-SS

Add "-SV" for Vents



Right-hinged door shown
-FHD "Full Height Doors"


Width	Left & Right Hinged
12"	SB3512-100-S-FHD-SS
15"	SB3515-100-S-FHD-SS
18"	SB3518-100-S-FHD-SS
21"	SB3521-100-S-FHD-SS
24"	SB3524-100-S-FHD-SS
12"	SB3512-100-S-FHD-SS
15"	SB3515-100-S-FHD-SS
18"	SB3518-100-S-FHD-SS
21"	SB3521-100-S-FHD-SS
24"	SB3524-100-S-FHD-SS



30" and wider

Width	Item Number
30"	SB3530-200-S-FHD-SS
36"	SB3536-200-S-FHD-SS
42"	SB3542-200-S-FHD-SS
48"	SB3548-200-S-FHD-SS

-FHD "Full Height Doors"



Sink Support Assembly

Item Number
SSS-C14
SSS-18
SSS-24
SSS-30
SSS-36
SSS-42
SSS-48
SSS-52
SSS-58



Width	Item Number
52"	SB3552-200-S-SS
58"	SB3558-200-S-SS


Available with vents:
*Change to -SV



Width	Item Number
52"	SB3552-200-SI-SS
58"	SB3558-200-SI-SS

Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 22.375" deep. Units are available in various widths as shown. Door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



OPEN CABINETS



Up to 24" wide

Width	Item Number
12"	SB2912-000-SS
15"	SB2915-000-SS
18"	SB2918-000-SS
21"	SB2921-000-SS
24"	SB2924-000-SS



30" and wider

Width	Item Number
30"	SB2930-000-SS
36"	SB2936-000-SS
42"	SB2942-000-SS
48"	SB2948-000-SS

DOOR CABINETS



Right-hinged door shown

Width	Left & Right Hinged
12"	SB2912-100L-SS
15"	SB2915-100L-SS
18"	SB2918-100L-SS
21"	SB2921-100L-SS
24"	SB2924-100L-SS
12"	SB2912-100R-SS
15"	SB2915-100R-SS
18"	SB2918-100R-SS
21"	SB2921-100R-SS
24"	SB2924-100R-SS



30" and wider

Width	Item Number
30"	SB2930-200-SS
36"	SB2936-200-SS
42"	SB2942-200-SS
48"	SB2948-200-SS



Example: SB2930-100-LPNL

Width	Left & Right Hinged
30"	SB2930-100-LPNL-SS
36"	SB2936-100-LPNL-SS
42"	SB2942-100-LPNL-SS
48"	SB2948-100-LPNL-SS
30"	SB2930-100-RPNL-SS
36"	SB2936-100-RPNL-SS
42"	SB2942-100-RPNL-SS
48"	SB2948-100-RPNL-SS


Note: L(left) or R(right) refers to door hinge side.
PNL(panel) refers to panel next to door.

Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB2930-022-LK** — Lock

Height: 28.75"
Doors: 2
Half Drawers: 0
Width: 30"
Drawers: 2



DRAWER CABINETS



Up to 24" wide

Width	Item Number
12"	SB2912-030-SS
15"	SB2915-030-SS
18"	SB2918-030-SS
21"	SB2921-030-SS
24"	SB2924-030-SS



30" and wider

Width	Item Number
30"	SB2930-030-SS
36"	SB2936-030-SS
42"	SB2942-030-SS
48"	SB2948-030-SS



Up to 24" wide

Width	Item Number
12"	SB2912-022-SS
15"	SB2915-022-SS
18"	SB2918-022-SS
21"	SB2921-022-SS
24"	SB2924-022-SS



30" and wider

Width	Item Number
30"	SB2930-022-SS
36"	SB2936-022-SS
42"	SB2942-022-SS
48"	SB2948-022-SS



Up to 24" wide

Width	Item Number
12"	SB2912-014-SS
15"	SB2915-014-SS
18"	SB2918-014-SS
21"	SB2921-014-SS
24"	SB2924-014-SS



30" and wider

Width	Item Number
30"	SB2930-014-SS
36"	SB2936-014-SS
42"	SB2942-014-SS
48"	SB2948-014-SS

Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB2930-040-2F-LK**

Height: 28.75"
Doors: 2
Half Drawers: 0
Width: 30"
Drawers: 4
File Drawers: 0
Lock: LK



DRAWER CABINETS



Up to 24" wide

Width	Item Number
12"	SB2912-006-SS
15"	SB2915-006-SS
18"	SB2918-006-SS
21"	SB2921-006-SS
24"	SB2924-006-SS



30" and wider

Width	Item Number
30"	SB2930-006-SS
36"	SB2936-006-SS
42"	SB2942-006-SS
48"	SB2948-006-SS



Up to 24" wide

Width	Item Number
12"	SB2912-002-1F-SS
15"	SB2915-002-1F-SS
18"	SB2918-002-1F-SS
21"	SB2921-002-1F-SS
24"	SB2924-002-1F-SS



30" and wider

Width	Item Number
30"	SB2930-004-2F-SS
36"	SB2936-004-2F-SS
42"	SB2942-004-2F-SS
48"	SB2948-004-2F-SS



Up to 24" wide

Width	Item Number
15"	SB2915-000-2F-SS
18"	SB2918-000-2F-SS
21"	SB2921-000-2F-SS
24"	SB2924-000-2F-SS



30" and wider

Width	Item Number
30"	SB2930-000-4F-SS
36"	SB2936-000-4F-SS
42"	SB2942-000-4F-SS
48"	SB2948-000-4F-SS



Up to 24" wide

Width	Item Number
15"	SB2915-010-1F-SS
18"	SB2918-010-1F-SS
21"	SB2921-010-1F-SS
24"	SB2924-010-1F-SS



30" and wider

Width	Item Number
30"	SB2930-020-2F-SS
36"	SB2936-020-2F-SS
42"	SB2942-020-2F-SS
48"	SB2948-020-2F-SS

Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



Height

Doors

Half Drawers

Width

Drawers

Lock

Example: **SB2930-044-LK**

DRAWER CABINETS



30" and wider

Width	Item Number
30"	SB2930-060-SS
36"	SB2936-060-SS
42"	SB2942-060-SS
48"	SB2948-060-SS



30" and wider

Width	Item Number
30"	SB2930-044-SS
36"	SB2936-044-SS
42"	SB2942-044-SS
48"	SB2948-044-SS



30" and wider

Width	Item Number
30"	SB2930-040-SS
36"	SB2936-040-SS
42"	SB2942-040-SS
48"	SB2948-040-SS



30" and wider

Width	Item Number
30"	SB2930-028-SS
36"	SB2936-028-SS
42"	SB2942-028-SS
48"	SB2948-028-SS

Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 22.375" deep. Units are available in various widths as shown. No shelves are available in sitting height base cabinets with the combination of doors and drawers. Drawer units feature 100 lb. full extension slides as standard. Specify left or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



Example: **SB2930-110R-LK** — Lock

Height: 28.75"
Doors: 110"
Half Drawers: 110"
Width: 30"
Drawers: 110"
Right Hinged: R

DOOR & DRAWER CABINETS



Right-hinged door shown

Width	Left & Right Hinged
12"	SB2912-110L-SS
15"	SB2915-110L-SS
18"	SB2918-110L-SS
21"	SB2921-110L-SS
24"	SB2924-110L-SS
12"	SB2912-110R-SS
15"	SB2915-110R-SS
18"	SB2918-110R-SS
21"	SB2921-110R-SS
24"	SB2924-110R-SS



Right-hinged door shown

Width	Left & Right Hinged
12"	SB2912-102L-SS
15"	SB2915-102L-SS
18"	SB2918-102L-SS
21"	SB2921-102L-SS
24"	SB2924-102L-SS
12"	SB2912-102R-SS
15"	SB2915-102R-SS
18"	SB2918-102R-SS
21"	SB2921-102R-SS
24"	SB2924-102R-SS



30" and wider

Width	Item Number
30"	SB2930-220-SS
36"	SB2936-220-SS
42"	SB2942-220-SS
48"	SB2948-220-SS



30" and wider

Width	Item Number
30"	SB2930-212-SS
36"	SB2936-212-SS
42"	SB2942-212-SS
48"	SB2948-212-SS



30" and wider

Width	Item Number
30"	SB2930-204-SS
36"	SB2936-204-SS
42"	SB2942-204-SS
48"	SB2948-204-SS




Right-hinged door shown

Width	Left & Right Hinged
30"	SB2930-132L-SS
36"	SB2936-132L-SS
42"	SB2942-132L-SS
48"	SB2948-132L-SS
30"	SB2930-132R-SS
36"	SB2936-132R-SS
42"	SB2942-132R-SS
48"	SB2948-132R-SS

Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 22.375" deep. Units are available in various widths as shown. No shelves are available in sitting height base cabinets with the combination of doors and drawers. Drawer units feature 100 lb. full extension slides as standard. Door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



DOOR & DRAWER CABINETS



30" and wider

Width	Item Number
30"	SB2930-210-SS
36"	SB2936-210-SS
42"	SB2942-210-SS
48"	SB2948-210-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB2930-210-AL-SS
36"	SB2936-210-AL-SS
42"	SB2942-210-AL-SS
48"	SB2948-210-AL-SS
30"	SB2930-210-AR-SS
36"	SB2936-210-AR-SS
42"	SB2942-210-AR-SS
48"	SB2948-210-AR-SS



30" and wider

Width	Item Number
30"	SB2930-202-A-SS
36"	SB2936-202-A-SS
42"	SB2942-202-A-SS
48"	SB2948-202-A-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB2930-202L-SS
36"	SB2936-202L-SS
42"	SB2942-202L-SS
48"	SB2948-202L-SS
30"	SB2930-202R-SS
36"	SB2936-202R-SS
42"	SB2942-202R-SS
48"	SB2948-202R-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB2930-130L-SS
36"	SB2936-130L-SS
42"	SB2942-130L-SS
48"	SB2948-130L-SS
30"	SB2930-130R-SS
36"	SB2936-130R-SS
42"	SB2942-130R-SS
48"	SB2948-130R-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB2930-130-AL-SS
36"	SB2936-130-AL-SS
42"	SB2942-130-AL-SS
48"	SB2948-130-AL-SS
30"	SB2930-130-AR-SS
36"	SB2936-130-AR-SS
42"	SB2942-130-AR-SS
48"	SB2948-130-AR-SS

Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 22.375" deep. Units are available in various widths as shown. No shelves are available in sitting height base cabinets with the combination of doors and drawers. Drawer units feature 100 lb. full extension slides as standard. Door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a vertical divider is required between cup board and drawer section use the suffix "VD". If a finished back is required use the suffix "FB".



DOOR & DRAWER CABINETS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB2930-116L-SS
36"	SB2936-116L-SS
42"	SB2942-116L-SS
48"	SB2948-116L-SS
30"	SB2930-116R-SS
36"	SB2936-116R-SS
42"	SB2942-116R-SS
48"	SB2948-116R-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB2930-116-AL-SS
36"	SB2936-116-AL-SS
42"	SB2942-116-AL-SS
48"	SB2948-116-AL-SS
30"	SB2930-116-AR-SS
36"	SB2936-116-AR-SS
42"	SB2942-116-AR-SS
48"	SB2948-116-AR-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB2930-124L-SS
36"	SB2936-124L-SS
42"	SB2942-124L-SS
48"	SB2948-124L-SS
30"	SB2930-124R-SS
36"	SB2936-124R-SS
42"	SB2942-124R-SS
48"	SB2948-124R-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB2930-124-AL-SS
36"	SB2936-124-AL-SS
42"	SB2942-124-AL-SS
48"	SB2948-124-AL-SS
30"	SB2930-124-AR-SS
36"	SB2936-124-AR-SS
42"	SB2942-124-AR-SS
48"	SB2948-124-AR-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB2930-114L-SS
36"	SB2936-114L-SS
42"	SB2942-114L-SS
48"	SB2948-114L-SS
30"	SB2930-114R-SS
36"	SB2936-114R-SS
42"	SB2942-114R-SS
48"	SB2948-114R-SS



Right-hinged door shown

Width	Left & Right Hinged
30"	SB2930-108L-SS
36"	SB2936-108L-SS
42"	SB2942-108L-SS
48"	SB2948-108L-SS
30"	SB2930-108R-SS
36"	SB2936-108R-SS
42"	SB2942-108R-SS
48"	SB2948-108R-SS

SINK & FUME HOOD CABINETS



Width	Left Hinged	Right Hinged
12"	SB2912-100-SL-SS	SB2912-100-SR-SS
15"	SB2915-100-SL-SS	SB2915-100-SR-SS
18"	SB2918-100-SL-SS	SB2918-100-SR-SS
21"	SB2921-100-SL-SS	SB2921-100-SR-SS
24"	SB2924-100-SL-SS	SB2924-100-SR-SS




Width	Left Hinged	Right Hinged
12"	SB2912-100-SL-FHD-SS	SB2912-100-SR-FHD-SS
15"	SB2915-100-SL-FHD-SS	SB2915-100-SR-FHD-SS
18"	SB2918-100-SL-FHD-SS	SB2918-100-SR-FHD-SS
21"	SB2921-100-SL-FHD-SS	SB2921-100-SR-FHD-SS
24"	SB2924-100-SL-FHD-SS	SB2924-100-SR-FHD-SS

-FHD "Full Height Doors"

Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Full height door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



SINK & FUME HOOD CABINETS




Width	Item Number
30"	SB2930-200-S-SS
36"	SB2936-200-S-SS
42"	SB2942-200-S-SS
48"	SB2948-200-S-SS

Add "-SV" for Vents



Width	Item Number
30"	SB2930-200-S-FHD-SS
36"	SB2936-200-S-FHD-SS
42"	SB2942-200-S-FHD-SS
48"	SB2948-200-S-FHD-SS

-FHD "Full Height Doors"



Sink Support Assembly

Item Number
SSS-C14
SSS-18
SSS-24
SSS-30
SSS-36
SSS-42
SSS-48
SSS-52
SSS-58



Width	Item Number
52"	SB2952-200-S-SS
58"	SB2958-200-S-SS



Width	Item Number
52"	SB2952-200-SI-SS
58"	SB2958-200-SI-SS

CORNER CABINETS



Right-hinged door shown

Item Number
Left & Right Hinged
SC2914-100L-SS
SC2914-100R-SS


Lazy Susan Add-On Suffix
"LS"



Right-hinged door shown

Left Hinged	Right Hinged
SC2914-100-SL-SS	SC2914-100-SR-SS

Left Hinged with Louvers	Right Hinged with Louvers
SC2914-100-SVL-SS	SC2914-100-SVR-SS



Base Cabinets — ADA Height

ADA height cabinets are 31.75" high and 22.375" deep. Units are available in various widths as shown. Door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



OPEN CABINETS



Up to 24" wide

Width	Item Number
12"	SB3212-000-SS
15"	SB3215-000-SS
18"	SB3218-000-SS
21"	SB3221-000-SS
24"	SB3224-000-SS



30" and wider

Width	Item Number
30"	SB3230-000-SS
36"	SB3236-000-SS
42"	SB3242-000-SS
48"	SB3248-000-SS

DOOR CABINETS



Right-hinged door shown

Width	Left & Right Hinged
12"	SB3212-100L-SS
15"	SB3215-100L-SS
18"	SB3218-100L-SS
21"	SB3221-100L-SS
24"	SB3224-100L-SS
12"	SB3212-100R-SS
15"	SB3215-100R-SS
18"	SB3218-100R-SS
21"	SB3221-100R-SS
24"	SB3224-100R-SS



30" and wider

Width	Item Number
30"	SB3230-200-SS
36"	SB3236-200-SS
42"	SB3242-200-SS
48"	SB3248-200-SS

Base Cabinets — ADA Height

ADA height cabinets are 31.75" high and 22.375" deep. Units are available in various widths as shown. Full height door units feature single shelf as standard. Additional shelves are available. Drawer units feature 100 lb. full extension slides as standard. If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



DRAWER CABINETS



Width	Item Number
12"	SB3212-031-SS
15"	SB3215-031-SS
18"	SB3218-031-SS
21"	SB3221-031-SS
24"	SB3224-031-SS

12" and wider



Width	Item Number
30"	SB3230-031-SS
36"	SB3236-031-SS
42"	SB3242-031-SS
48"	SB3248-031-SS

30" and wider



Width	Item Number
12"	SB3212-020-1F-SS
15"	SB3215-020-1F-SS
18"	SB3218-020-1F-SS
21"	SB3221-020-1F-SS
24"	SB3224-020-1F-SS

Up to 24" wide



Width	Item Number
30"	SB3230-040-2F-SS
36"	SB3236-040-2F-SS
42"	SB3242-040-2F-SS
48"	SB3248-040-2F-SS

30" and wider



Width	Item Number
15"	SB3215-012-1F-SS
18"	SB3218-012-1F-SS
21"	SB3221-012-1F-SS
24"	SB3224-012-1F-SS

Up to 24" wide



Width	Item Number
30"	SB3230-024-2F-SS
36"	SB3236-024-2F-SS
42"	SB3242-024-2F-SS
48"	SB3248-024-2F-SS

30" and wider



Width	Item Number
30"	SB3230-062-SS
36"	SB3236-062-SS
42"	SB3242-062-SS
48"	SB3248-062-SS

30" and wider

Base Cabinets — ADA Height

ADA height cabinets are 31.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Full height door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



DOOR & DRAWER CABINETS



Right-hinged door shown

Width	Left & Right Hinged
12"	SB3212-110L-SS
15"	SB3215-110L-SS
18"	SB3218-110L-SS
21"	SB3221-110L-SS
24"	SB3224-110L-SS
12"	SB3212-110R-SS
15"	SB3215-110R-SS
18"	SB3218-110R-SS
21"	SB3221-110R-SS
24"	SB3224-110R-SS



30" and wider

Width	Item Number
30"	SB3230-202-A-SS
36"	SB3236-202-A-SS
42"	SB3242-202-A-SS
48"	SB3248-202-A-SS



30" and wider

Width	Item Number
30"	SB3230-210-SS
36"	SB3236-210-SS
42"	SB3242-210-SS
48"	SB3248-210-SS



30" and wider

Width	Item Number
30"	SB3230-201-SS
36"	SB3236-201-SS
42"	SB3242-201-SS
48"	SB3248-201-SS



30" and wider

Width	Item Number
30"	SB3230-220-SS
36"	SB3236-220-SS
42"	SB3242-220-SS
48"	SB3248-220-SS



30" and wider

Width	Item Number
30"	SB3230-202-SS
36"	SB3236-202-SS
42"	SB3242-202-SS
48"	SB3248-202-SS

Base Cabinets — ADA Height

ADA height cabinets are 31.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



CORNER CABINETS



Right-hinged door shown

Left Hinged	Right Hinged
SC3214-100L-LS-SS	SC3214-100R-LS-SS

Lazy Susan Add-On Suffix
"LS"


Left Hinged	Right Hinged
SC3214-100L-SS	SC3214-100R-SS



Right-hinged door shown

Left Hinged	Right Hinged
SC3214-100-SL-SS	SC3214-100-SR-SS

Left Hinged with Louvers	Right Hinged with Louvers
SC3214-100-SVL-SS	SC3214-100-SVR-SS



Base Cabinets — ADA Height

ADA height cabinets are 31.75" high and 22.375" deep. Units are available in various widths as shown. Full height door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



SINK & FUME HOOD CABINETS



Right-hinged door shown
Add "-SV(R/L)" for Vents

Width	Left & Right Hinged
12"	SB3212-100-SL-SS
15"	SB3215-100-SL-SS
18"	SB3218-100-SL-SS
21"	SB3221-100-SL-SS
24"	SB3224-100-SL-SS
12"	SB3224-100-SR-SS
15"	SB3224-100-SR-SS
18"	SB3224-100-SR-SS
21"	SB3221-100-SR-SS
24"	SB3224-100-SR-SS



Add "-SV" for Vents

Width	Item Number
30"	SB3230-200-S-SS
36"	SB3236-200-S-SS
42"	SB3242-200-S-SS
48"	SB3248-200-S-SS

Add "-SV" for Vents



-FHD "Full Height Doors"

Width	Left & Right Hinged
12"	SB3212-100-SL-FHD-SS
15"	SB3215-100-SL-FHD-SS
18"	SB3218-100-SL-FHD-SS
21"	SB3221-100-SL-FHD-SS
24"	SB3224-100-SL-FHD-SS
12"	SB3224-100-SR-FHD-SS
15"	SB3224-100-SR-FHD-SS
18"	SB3224-100-SR-FHD-SS
21"	SB3221-100-SR-FHD-SS
24"	SB3224-100-SR-FHD-SS



-FHD "Full Height Doors"

Width	Item Number
30"	SB3230-200-S-FHD-SS
36"	SB3236-200-S-FHD-SS
42"	SB3242-200-S-FHD-SS
48"	SB3248-200-S-FHD-SS

-FHD "Full Height Doors"



Width	Item Number
52"	SB3252-200-S-SS
58"	SB3258-200-S-SS



Width	Item Number
52"	SB3252-200-SI-SS
58"	SB3258-200-SI-SS

Base Cabinets — ADA Height



ADA height cabinets are 31.75" high and 22.375" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".




SINK & FUME HOOD CABINETS



Width	Item Number
36"	SADAKSC36-SS
42"	SADAKSC42-SS
48"	SADAKSC48-SS



Width	Item Number
30"	SADAKSP30-SS
36"	SADAKSP36-SS
42"	SADAKSP42-SS
48"	SADAKSP48-SS



Sink Support Assembly

Item Number
SSS-C14
SSS-18
SSS-24
SSS-30
SSS-36
SSS-42
SSS-48
SSS-52
SSS-58


Base Cabinets — Corrosives, Bases, Oxidizers

Corrosives storage cabinets are reinforced steel with a hi-tech performance polymer alloy coating, which provides 100% coverage against corrosion. Cabinets are available in 18" and 22" depths.

Example: **SASB3518-22R-VK** — Vent Kit

Height
Depth
Width
Right Hinged

Removable Panel



CORROSIVES STORAGE CABINETS




SASB3524-18-SS




SASB3530-18-SS

Width	Sitting Height	ADA Height	Standing Height
18"	SASB2918-18-SS	SASB3218-18-SS	SASB3518-18-SS
24"	SASB2924-18-SS	SASB3224-18-SS	SASB3524-18-SS
30"	SASB2930-18-SS	SASB3230-18-SS	SASB3530-18-SS
36"	SASB2936-18-SS	SASB3236-18-SS	SASB3536-18-SS
42"	SASB2942-18-SS	SASB3242-18-SS	SASB3542-18-SS
48"	SASB2948-18-SS	SASB3248-18-SS	SASB3548-18-SS

Width	Sitting Height	ADA Height	Standing Height
18"	SASB2918-22-SS	SASB3218-22-SS	SASB3518-22-SS
24"	SASB2924-22-SS	SASB3224-22-SS	SASB3524-22-SS
30"	SASB2930-22-SS	SASB3230-22-SS	SASB3530-22-SS
36"	SASB2936-22-SS	SASB3236-22-SS	SASB3536-22-SS
42"	SASB2942-22-SS	SASB3242-22-SS	SASB3542-22-SS
48"	SASB2948-22-SS	SASB3248-22-SS	SASB3548-22-SS



18-3/8" ACTUAL DEPTH



22-3/8" ACTUAL DEPTH


Base Cabinets — Corrosives, Bases, Oxidizers

Corrosives storage cabinets are reinforced steel with a hi-tech performance polymer alloy coating, which provides 100% coverage against corrosion. Cabinets are available in 18" and 22" depths.

Example: **SAB3518-22R-VK** — Vent Kit


Height Depth
Width Right Hinged

Removable Panel




CORROSIVES STORAGE CABINETS

24" Wide (1 DOOR)



30" Wide (2 DOOR)



SASBM2930-22-SS

Mobile cabinets are available. Please contact your AMS representative.

OXIDIZERS

For more information regarding our Oxidizer cabinets, please contact your AMS representative.

Width	Sitting Height	ADA Height	Standing Height	Width	Sitting Height	ADA Height	Standing Height
18"	SASBM2918-18-SS	SASBM3218-18-SS	SASBM3518-18-SS	18"	SASBM2918-22-SS	SASBM3218-22-SS	SASBM3518-22-SS
24"	SASBM2924-18-SS	SASBM3224-18-SS	SASBM3524-18-SS	24"	SASBM2924-22-SS	SASBM3224-22-SS	SASBM3524-22-SS
30"	SASBM2930-18-SS	SASBM3230-18-SS	SASBM3530-18-SS	30"	SASBM2930-22-SS	SASBM3230-22-SS	SASBM3530-22-SS
36"	SASBM2936-18-SS	SASBM3236-18-SS	SASBM3536-18-SS	36"	SASBM2936-22-SS	SASBM3236-22-SS	SASBM3536-22-SS
42"	SASBM2942-18-SS	SASBM3242-18-SS	SASBM3542-18-SS	42"	SASBM2942-22-SS	SASBM3242-22-SS	SASBM3542-22-SS
48"	SASBM2948-18-SS	SASBM3248-18-SS	SASBM3548-18-SS	48"	SASBM2948-22-SS	SASBM3248-22-SS	SASBM3548-22-SS

When ordering 18" or 24" wide cabinets please place suffix "L" left or "R" right, after the part number for the hinge placement.


If a Vent Kit is required, use the suffix "VK".

Note: AMS recommends vent kits for all corrosives cabinets. Acids and bases should not be stored in the same cabinet.


Corrosive - Vent Kit

Vent systems shall be manufactured from polyolefin resin and shall conform to ASTM #1412-92 for corrosion resistance and material content. The vent system will contain no PVC or solvents used for assembly. Vent may be constructed using rigid or flexible pipe as required. Vents will be manufactured by Scientific Plastics, Inc. of Kansas City, KS, USA.

403100-VENT KIT



403100-VENT KIT-L




Base Cabinets — Flammables

Flammable storage cabinets are reinforced steel with a hi-tech performance polymer alloy coating, which provides 100% coverage against corrosion. Cabinets are available in 18" and 22" depths. All cabinets have 3 point locking systems.

Example: **SFSB3518-22R**

Height Depth
Width Right Hinged

Removable Panel



FLAMMABLES STORAGE CABINETS




SFSB2924-18-SS




SFSB2930-18-SS

Width	Sitting Height	ADA Height	Standing Height
24"	SFSB2924-18-SS	SFSB3224-18-SS	SFSB3524-18-SS
30"	SFSB2930-18-SS	SFSB3230-18-SS	SFSB3530-18-SS
36"	SFSB2936-18-SS	SFSB3236-18-SS	SFSB3536-18-SS
42"	SFSB2942-18-SS	SFSB3242-18-SS	SFSB3542-18-SS
48"	SFSB2948-18-SS	SFSB3248-18-SS	SFSB3548-18-SS

Width	Sitting Height	ADA Height	Standing Height
24"	SFSB2924-22-SS	SFSB3224-22-SS	SFSB3524-22-SS
30"	SFSB2930-22-SS	SFSB3230-22-SS	SFSB3530-22-SS
36"	SFSB2936-22-SS	SFSB3236-22-SS	SFSB3536-22-SS
42"	SFSB2942-22-SS	SFSB3242-22-SS	SFSB3542-22-SS
48"	SFSB2948-22-SS	SFSB3248-22-SS	SFSB3548-22-SS



18-3/8" ACTUAL DEPTH




22 3/8" ACTUAL DEPTH

Base Cabinets — Flammables

Flammable storage cabinets are reinforced steel with a hi-tech performance polymer alloy coating, which provides 100% coverage against corrosion. Cabinets are available in 18" and 22" depths. All cabinets have 3 point locking systems.


Example: **SFSB3518-22R**

Height Depth
Width Right Hinged




FLAMMABLES STORAGE CABINETS

24" Wide (1 DOOR)



30" Wide (2 DOOR)



SFSBM3530-18-SS

Mobile cabinets are available. Please contact your AMS representative.

Width	Sitting Height	ADA Height	Standing Height
24"	SFSBM2924-18-SS	SFSBM3224-18-SS	SFSBM3524-18-SS
30"	SFSBM2930-18-SS	SFSBM3230-18-SS	SFSBM3530-18-SS
36"	SFSBM2936-18-SS	SFSBM3236-18-SS	SFSBM3536-18-SS
42"	SFSBM2942-18-SS	SFSBM3242-18-SS	SFSBM3542-18-SS
48"	SFSBM2948-18-SS	SFSBM3248-18-SS	SFSBM3548-18-SS

Width	Sitting Height	ADA Height	Standing Height
24"	SFSBM2924-22-SS	SFSBM3224-22-SS	SFSBM3524-22-SS
30"	SFSBM2930-22-SS	SFSBM3230-22-SS	SFSBM3530-22-SS
36"	SFSBM2936-22-SS	SFSBM3236-22-SS	SFSBM3536-22-SS
42"	SFSBM2942-22-SS	SFSBM3242-22-SS	SFSBM3542-22-SS
48"	SFSBM2948-22-SS	SFSBM3248-22-SS	SFSBM3548-22-SS

Flammables storage cabinets are certified to meet NFPA 45 and are UL listed. Cabinets are available in 18" and 22" depths. For 22" depth, add suffix "-22" to the part number. For synchronized self-closing doors add "-SCC" to the part number.

*Cabinets with self-closing doors are constructed to the UL1275 standards, but cannot have a UL label

Vacuum Pump Cabinets

Vacuum Pump Cabinets come with a Polydamp Melamine Foam that has heat resistance and sound dampening qualities. This is applied to the inside of the cabinet. (Please note actual cabinet depth with finished back will be 22.875" and 18.875").



Width	Item Number
15"	SB3515-10V-SS
18"	SB3518-10V-SS
21"	SB3521-10V-SS
24"	SB3524-10V-SS

When ordering 18" or 24" wide cabinets please place suffix "L" left or "R" right, after the part number for the hinge placement.


Attention: All of these units are shipped with a loose unassembled junction box, switch and wiring. These are to be installed by your certified electrician.



Width	Item Number
30"	SB2930-20V-18-SS
36"	SB2936-20V-18-SS
42"	SB2942-20V-18-SS
48"	SB2948-20V-18-SS

Width	Item Number
30"	SB2930-20V-SS
36"	SB2936-20V-SS
42"	SB2942-20V-SS
48"	SB2948-20V-SS

Pullout Shelf Dimensions for a SB3524-10V-22:



Width	Item Number
30"	SB3230-20V-18-SS
36"	SB3236-20V-18-SS
42"	SB3242-20V-18-SS
48"	SB3248-20V-18-SS

Width	Item Number
30"	SB3230-20V-SS
36"	SB3236-20V-SS
42"	SB3242-20V-SS
48"	SB3248-20V-SS

Full Cabinet Width	Pullout Shelf Width
15"	7.4"
18"	10.4"
21"	13.4"
24"	16.4"
30"	22.4"
36"	28.4"
42"	34.4"
48"	40.4"



Width	Item Number
30"	SB3530-20V-18-SS
36"	SB3536-20V-18-SS
42"	SB3542-20V-18-SS
48"	SB3548-20V-18-SS

Width	Item Number
30"	SB3530-20V-SS
36"	SB3536-20V-SS
42"	SB3542-20V-SS
48"	SB3548-20V-SS

Full Cabinet Depth	Pullout Shelf Depth
18"	14.5"
22"	18.5"

Vacuum Pump Cabinets

Vacuum Pump Cabinets come with a Polydamp Melamine Foam that has heat resistance and sound dampening qualities. This is applied to the inside of the cabinet. (Please note actual cabinet depth with finished back will be 22.875" and 18.875").



Attention: These units are shipped with a loose unassembled junction box, switch and wiring. These are to be installed by your certified electrician.



These units include the cart pictured as part of the unit price. Vent pipe is not included.

Width	Item Number
30"	SB2930-20V-18-NB-SS
36"	SB2936-20V-18-NB-SS
42"	SB2942-20V-18-NB-SS
48"	SB2948-20V-18-NB-SS

Width	Item Number
30"	SB2930-20V-NB-SS
36"	SB2936-20V-NB-SS
42"	SB2942-20V-NB-SS
48"	SB2948-20V-NB-SS

Width	Item Number
30"	SB3230-20V-18-NB-SS
36"	SB3236-20V-18-NB-SS
42"	SB3242-20V-18-NB-SS
48"	SB3248-20V-18-NB-SS

Width	Item Number
30"	SB3230-20V-NB-SS
36"	SB3236-20V-NB-SS
42"	SB3242-20V-NB-SS
48"	SB3248-20V-NB-SS

Width	Item Number
30"	SB3530-20V-18-NB-SS
36"	SB3536-20V-18-NB-SS
42"	SB3542-20V-18-NB-SS
48"	SB3548-20V-18-NB-SS

Width	Item Number
30"	SB3530-20V-NB-SS
36"	SB3536-20V-NB-SS
42"	SB3542-20V-NB-SS
48"	SB3548-20V-NB-SS

Width	Item Number
24"	MOBILE-CART-24-SS
36"	MOBILE-CART-36-SS

(Standard depth of 22")
*Actual Depth of cart: 17.186"

Width	Item Number
24"	MOBILE-CART-24-18-SS
36"	MOBILE-CART-36-18-SS

(Suffix 18 represents 18" depth)
*Actual Depth of cart: 13.186"

Mobile Cabinets

All AMS cabinets are available in a mobile version including low profile (27.18" high), sitting height (28.75" high), ADA height (31.75" high), and standing height (34.75" high). They come standard with a finished back and optional finished top. To specify a mobile cabinet use the prefix "M" to the base cabinet model number (example: SMB2724-110). They come standard with a 3" wheel caster. There is no height changes to the cabinets when replaced toe kicks with casters. 2" casters are available by adding the suffix "-2" to the part number. Cabinet depth with finished back is 22.812". All mobile cabinets with drawers come standard with a tip resistant counter weight.

Tops: Please add the suffix "T" for steel, "E" for epoxy, "P" for phenolic. All tops are 3/4" thick and with a 7/8" overhang on front only.
(example SMB2724-110L-T)



Mobile Cabinets

Standing height mobile cabinets are 34.75" high and 22.375" deep. Units are available in various widths. Each example below is available in the standard four heights shown on the left. Door units feature a single shelf as standard for standing height cabinets. All open units have 2 shelves. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK".



DOOR & DRAWER CABINETS



Width	Item Number
12"	SMB3512-100R-SS
15"	SMB3515-100R-SS
18"	SMB3518-100R-SS
21"	SMB3521-100R-SS
24"	SMB3524-100R-SS
12"	SMB3512-100L-SS
15"	SMB3515-100L-SS
18"	SMB3518-100L-SS
21"	SMB3521-100L-SS
24"	SMB3524-100L-SS



Width	Item Number
30"	SMB3530-200-SS
36"	SMB3536-200-SS
42"	SMB3542-200-SS
48"	SMB3548-200-SS



Width	Left & Right Hinged
12"	SMB3512-110L-SS
15"	SMB3515-110L-SS
18"	SMB3518-110L-SS
21"	SMB3521-110L-SS
24"	SMB3524-110L-SS
12"	SMB3512-110R-SS
15"	SMB3515-110R-SS
18"	SMB3518-110R-SS
21"	SMB3521-110R-SS
24"	SMB3524-110R-SS

Note: L(left) or R(right) refers to door hinge side.



Width	Item Number
30"	SMB3530-220-SS
36"	SMB3536-220-SS
42"	SMB3542-220-SS
48"	SMB3548-220-SS



Width	Item Number
30"	SMB3530-210-SS
36"	SMB3536-210-SS
42"	SMB3542-210-SS
48"	SMB3548-210-SS


Suspended Cabinets

All AMS cabinets are available in a suspended version including low profile (23.18" high), sitting height (24.875" high), ADA height (27.875" high), and standing height (30.875" high). To specify a suspended cabinet use the prefix "H" to the base cabinet model number (example: SHB2324-110). Finished backs are optional and can be specified with the suffix "-FB". Low profile cabinets should be used for adjustable height work benches or table frames on casters to allow the caster to swivel if the adjustable leg is positioned in the lowest setting.



Suspended Cabinets

Standing height suspended cabinets are 30.875" high and 22.375" deep. Units are available in various widths. Each example below is available in the standard four heights shown on the left. Door units feature a single shelf as standard for standard height cabinets. All open units have 2 shelves. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



DOOR & DRAWER CABINETS



Width	Item Number
12"	SHB3112-100R-SS
15"	SHB3115-100R-SS
18"	SHB3118-100R-SS
21"	SHB3121-100R-SS
24"	SHB3124-100R-SS
12"	SHB3112-100L-SS
15"	SHB3115-100L-SS
18"	SHB3118-100L-SS
21"	SHB3121-100L-SS
24"	SHB3124-100L-SS



Width	Item Number
30"	SHB3130-200-SS
36"	SHB3136-200-SS
42"	SHB3542-200-SS
48"	SHB3148-200-SS



Width	Left & Right Hinged
12"	SHB3112-110L-SS
15"	SHB3115-110L-SS
18"	SHB3118-110L-SS
21"	SHB3121-110L-SS
24"	SHB3124-110L-SS
12"	SHB3112-110R-SS
15"	SHB3115-110R-SS
18"	SHB3118-110R-SS
21"	SHB3121-110R-SS
24"	SHB3124-110R-SS

Note: L(left) or R(right) refers to door hinge side.



Width	Item Number
30"	SHB3130-220-SS
36"	SHB3136-220-SS
42"	SHB3142-220-SS
48"	SHB3148-220-SS




Width	Item Number
30"	SHB3130-210-SS
36"	SHB3136-210-SS
42"	SHB3142-210-SS
48"	SHB3148-210-SS

Overlay Cabinets


All AMS base cabinets are available in an overlay version including low profile (27.18" high), sitting height (28.75" high), ADA height (31.75" high), and standing height (34.75" high). To specify an overlay cabinet, use the prefix "O" to the base cabinet model number (example: SOB2724-110). Overlay cabinets are also available in mobile and suspended variations.

Tops: Please add the suffix "T" for steel, "E" for epoxy, "P" for phenolic. All tops are 3/4" thick and with a 7/8" overhang on front only.
(example SOB2724-110L-T)



Overlay Cabinets

Standing height overlay base cabinets are 34.75" high and 23.248" deep (Including doors and drawers depth). Units are available in various widths. Each example below is available in the standard four heights shown on the left. Door units feature a single shelf as standard. All open units have 2 shelves. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK".



DOOR & DRAWER CABINETS



Width	Item Number
12"	SOB3512-100R-SS
15"	SOB3515-100R-SS
18"	SOB3518-100R-SS
21"	SOB3521-100R-SS
24"	SOB3524-100R-SS
12"	SOB3512-100L-SS
15"	SOB3515-100L-SS
18"	SOB3518-100L-SS
21"	SOB3521-100L-SS
24"	SOB3524-100L-SS



Width	Item Number
30"	SOB3530-200-SS
36"	SOB3536-200-SS
42"	SOB3542-200-SS
48"	SOB3548-200-SS



Width	Left & Right Hinged
12"	SOB3512-110L-SS
15"	SOB3515-110L-SS
18"	SOB3518-110L-SS
21"	SOB3521-110L-SS
24"	SOB3524-110L-SS
12"	SOB3512-110R-SS
15"	SOB3515-110R-SS
18"	SOB3518-110R-SS
21"	SOB3521-110R-SS
24"	SOB3524-110R-SS

Note: L(left) or R(right) refers to door hinge side.



Width	Item Number
30"	SOB3530-220-SS
36"	SOB3536-220-SS
42"	SOB3542-220-SS
48"	SOB3548-220-SS



Width	Item Number
30"	SOB3530-210-SS
36"	SOB3536-210-SS
42"	SOB3542-210-SS
48"	SOB3548-210-SS

Tall Cabinets

Tall cabinets are 84" high and available in 16", 18", 21" and 24" deep by using the two digit suffix. Suffix is not needed when selecting 16" deep. Units are available in various widths as shown. Door units feature five shelves as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required use the suffix "LK".

Wardrobe cabinets are available by adding the suffix "W" to any of the standard part numbers. They will include one hanging bar in lieu of shelves.



Example: **SF8418-1GHR-18** — Depth

Labels for the example part number:
 Height: 84
 # of Doors: 1
 Right Hinged: H
 Width: 18
 Door Style: G

OPEN CABINETS



Up to 24" wide

Width	Item Number
18"	SF8418-000-SS
24"	SF8424-000-SS



30" and wider

Width	Item Number
30"	SF8430-000-SS
36"	SF8436-000-SS
42"	SF8442-000-SS
48"	SF8448-000-SS

DOOR CABINETS



Right-hinged door shown

Width	Left & Right Hinged
18"	SF8418-1SHL-SS
24"	SF8424-1SHL-SS
18"	SF8418-1SHR-SS
24"	SF8424-1SHR-SS



Width	Item Number
30"	SF8430-2SH-SS
36"	SF8436-2SH-SS
42"	SF8442-2SH-SS
48"	SF8448-2SH-SS



Width	Item Number
30"	SF8430-2SS-SS
36"	SF8436-2SS-SS
42"	SF8442-2SS-SS
48"	SF8448-2SS-SS

(this unit has sliding doors)



Width	Item Number
30"	SF8430-4SH-SS
36"	SF8436-4SH-SS
42"	SF8442-4SH-SS
48"	SF8448-4SH-SS

Tall Cabinets

Tall cabinets are 84" high and available in 16", 18", 21" and 24" deep by using the two digit suffix. Suffix is not needed when selecting 16" deep. Units are available in various widths as shown. Door units feature five shelves as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required use the suffix "LK".

Wardrobe cabinets are available by adding the suffix "W" to any of the standard part numbers. They will include one hanging bar in lieu of shelves. For pull out shelves add the suffix "PS".



Example: **SF8418-1GHR-18** — Depth

Height # of Doors Right Hinged

Width Door Style

DOOR CABINETS



Width	Left & Right Hinged
18"	SF8418-1GHL-SS
24"	SF8424-1GHL-SS
18"	SF8418-1GHR-SS
24"	SF8424-1GHR-SS



Width	Item Number
30"	SF8430-2GH-SS
36"	SF8436-2GH-SS
42"	SF8442-2GH-SS
48"	SF8448-2GH-SS

Right-hinged door shown



Width	Item Number
30"	SF8430-2GSF-SS
36"	SF8436-2GSF-SS
42"	SF8442-2GSF-SS
48"	SF8448-2GSF-SS

(this unit has sliding doors)



Width	Item Number
30"	SF8430-4GH-SS
36"	SF8436-4GH-SS
42"	SF8442-4GH-SS
48"	SF8448-4GH-SS



Width	Item Number
30"	SF8430-2GH2SH-SS
36"	SF8436-2GH2SH-SS
42"	SF8442-2GH2SH-SS
48"	SF8448-2GH2SH-SS



Pull out shelves are available upon request.

Wall Cabinets — 24" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature one shelf as standard for 24" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".

Example: **SW2430-2SS-16**

Height: 24" # of Doors: 2 Depth: 16"

Width: 30" Door Style: SS



OPEN CABINETS



Up to 24" wide

Width	Item Number
18"	SW2418-000-SS
24"	SW2424-000-SS



30" and wider

Width	Item Number
30"	SW2430-000-SS
36"	SW2436-000-SS
42"	SW2442-000-SS
48"	SW2448-000-SS

STEEL-HINGED DOOR CABINETS



Right-hinged door shown

Width	Left & Right Hinged
18"	SW2418-1SHL-SS
24"	SW2424-1SHL-SS
18"	SW2418-1SHR-SS
24"	SW2424-1SHR-SS



30" and wider

Width	Item Number
30"	SW2430-2SH-SS
36"	SW2436-2SH-SS
42"	SW2442-2SH-SS
48"	SW2448-2SH-SS

STEEL SLIDING DOOR CABINETS



30" and wider

Width	Item Number
30"	SW2430-2SS-SS
36"	SW2436-2SS-SS
42"	SW2442-2SS-SS
48"	SW2448-2SS-SS

CORNER CABINETS



Item Number
SWC2418-000-SS




Left Hinged	Right Hinged
SWC2418-1SHL-SS	SWC2418-1SHR-SS

Right-hinged door shown



Left Hinged	Right Hinged
SWC2418-1GHL-SS	SWC2418-1GHR-SS



*Lazy Susan Add-On
Suffix "LS"

Wall Cabinets — 24" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature one shelf as standard for 24" units. Specify left- or right-hinged doors using the suffix "L" "R" after the two digit depth number. If locks are required, use the suffix "LK".

Example: **SW2430-2SS-16**

Height: 24" # of Doors: 2 Depth: 16"

Width: 30" Door Style: SS



13" 1"

FRAMED GLASS HINGED DOOR CABINETS



Right-hinged door shown

Width	Left & Right Hinged
18"	SW2418-1GHL-SS
24"	SW2424-1GHL-SS
18"	SW2418-1GHR-SS
24"	SW2424-1GHR-SS



30" and wider

Width	Item Number
30"	SW2430-2GH-SS
36"	SW2436-2GH-SS
42"	SW2442-2GH-SS
48"	SW2448-2GH-SS

FRAMED GLASS SLIDING DOOR CABINETS



30" and wider

Width	Item Number
30"	SW2430-2GSF-SS
36"	SW2436-2GSF-SS
42"	SW2442-2GSF-SS
48"	SW2448-2GSF-SS

GLASS SLIDING DOOR CABINETS



30" and wider

Width	Item Number
30"	SW2430-2GS-SS
36"	SW2436-2GS-SS
42"	SW2442-2GS-SS
48"	SW2448-2GS-SS


Wall Cabinets — 30" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature two shelves as standard for 30" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".

Example: **SW3030-2SS-16** — Depth

Height # of Doors

Width Door Style



OPEN CABINETS



Width	Item Number
18"	SW3018-000-SS
24"	SW3024-000-SS

Up to 24" wide



30" and wider

Width	Item Number
30"	SW3030-000-SS
36"	SW3036-000-SS
42"	SW3042-000-SS
48"	SW3048-000-SS

STEEL-HINGED DOOR CABINETS



Width	Left & Right Hinged
18"	SW3018-1SHL-SS
24"	SW3024-1SHL-SS
18"	SW3018-1SHR-SS
24"	SW3024-1SHR-SS

Right-hinged door shown



30" and wider

Width	Item Number
30"	SW3030-2SH-SS
36"	SW3036-2SH-SS
42"	SW3042-2SH-SS
48"	SW3048-2SH-SS

STEEL SLIDING DOOR CABINETS



30" and wider

Width	Item Number
30"	SW3030-2SS-SS
36"	SW3036-2SS-SS
42"	SW3042-2SS-SS
48"	SW3048-2SS-SS

CORNER CABINETS



Item Number
SWC3018-000-SS




Left Hinged	Right Hinged
SWC3018-1SHL-SS	SWC3018-1SHR-SS

Right-hinged door shown

*Lazy Susan Add-On Suffix "LS"



Left Hinged	Right Hinged
SWC3018-1GHL-SS	SWC3018-1GHR-SS




Wall Cabinets — 30" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature two shelves as standard for 30" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".

Example: **SW3030-2SS-16**

Height # of Doors Depth

Width Door Style



FRAMED GLASS HINGED DOOR CABINETS



Right-hinged door shown

Width	Left & Right Hinged
18"	SW3018-1GHL-SS
24"	SW3024-1GHL-SS
18"	SW3018-1GHR-SS
24"	SW3024-1GHR-SS



30" and wider

Width	Item Number
30"	SW3030-2GH-SS
36"	SW3036-2GH-SS
42"	SW3042-2GH-SS
48"	SW3048-2GH-SS

FRAMED GLASS SLIDING DOOR CABINETS



30" and wider

Width	Item Number
30"	SW3030-2GSF-SS
36"	SW3036-2GSF-SS
42"	SW3042-2GSF-SS
48"	SW3048-2GSF-SS

GLASS SLIDING DOOR CABINETS



30" and wider


Width	Item Number
30"	SW3030-2GS-SS
36"	SW3036-2GS-SS
42"	SW3042-2GS-SS
48"	SW3048-2GS-SS

Wall Cabinets — 36" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature two shelves as standard for 36" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".

Example: **SW3630-2SS-16**

Height # of Doors Depth
Width Door Style



OPEN CABINETS



Up to 24" wide

Width	Item Number
18"	SW3618-000-SS
24"	SW3624-000-SS



30" and wider

Width	Item Number
30"	SW3630-000-SS
36"	SW3636-000-SS
42"	SW3642-000-SS
48"	SW3648-000-SS

STEEL-HINGED DOOR CABINETS



Right-hinged door shown

Width	Left & Right Hinged
18"	SW3618-1SHL-SS
24"	SW3624-1SHL-SS
18"	SW3618-1SHR-SS
24"	SW3624-1SHR-SS



30" and wider

Width	Item Number
30"	SW3630-2SH-SS
36"	SW3636-2SH-SS
42"	SW3642-2SH-SS
48"	SW3648-2SH-SS

STEEL SLIDING DOOR CABINETS



30" and wider

Width	Item Number
30"	SW3630-2SS-SS
36"	SW3636-2SS-SS
42"	SW3642-2SS-SS
48"	SW3648-2SS-SS

CORNER CABINETS




Item Number
SWC3618-000-SS



Right-hinged door shown

*Lazy Susan Add-On Suffix "LS"




Left Hinged	Right Hinged
SWC3618-1SHL-SS	SWC3618-1SHR-SS

Left Hinged	Right Hinged
SWC3618-1GHL-SS	SWC3618-1GHR-SS

Wall Cabinets — 36" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature two shelves as standard for 36" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".



Example: **SW3630-2SS-16** — Depth

Height # of Doors
Width Door Style

FRAMED GLASS HINGED DOOR CABINETS



Width	Left & Right Hinged
18"	SW3618-1GHL-SS
24"	SW3624-1GHL-SS
18"	SW3618-1GHR-SS
24"	SW3624-1GHR-SS

Right-hinged door shown



Width	Item Number
30"	SW3630-2GH-SS
36"	SW3636-2GH-SS
42"	SW3642-2GH-SS
48"	SW3648-2GH-SS

30" and wider

FRAMED GLASS SLIDING DOOR CABINETS



Width	Item Number
30"	SW3630-2GSF-SS
36"	SW3636-2GSF-SS
42"	SW3642-2GSF-SS
48"	SW3648-2GSF-SS

30" and wider

GLASS SLIDING DOOR CABINETS




Width	Item Number
30"	SW3630-2GS-SS
36"	SW3636-2GS-SS
42"	SW3642-2GS-SS
48"	SW3648-2GS-SS

30" and wider

Wall Cabinets — 48" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature three shelves as standard for 48" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".



Example: **SW4830-2SS-16**

Height: 48" # of Doors: 2 Depth: 16"

Width: 30" Door Style: SS

OPEN CABINETS



Width	Item Number
18"	SW4818-000-SS
24"	SW4824-000-SS

Up to 24" wide



Width	Item Number
30"	SW4830-000-SS
36"	SW4836-000-SS
42"	SW4842-000-SS
48"	SW4848-000-SS

30" and wider

STEEL-HINGED DOOR CABINETS



Width	Left & Right Hinged
18"	SW4818-1SHL-SS
24"	SW4824-1SHL-SS
18"	SW4818-1SHR-SS
24"	SW4824-1SHR-SS

Right-hinged door shown



Width	Item Number
30"	SW4830-2SH-SS
36"	SW4836-2SH-SS
42"	SW4842-2SH-SS
48"	SW4848-2SH-SS

30" and wider

STEEL SLIDING DOOR CABINETS



Width	Item Number
30"	SW4830-2SS-SS
36"	SW4836-2SS-SS
42"	SW4842-2SS-SS
48"	SW4848-2SS-SS

30" and wider

CORNER CABINETS



Item Number
SWC4818-000-SS




Left Hinged	Right Hinged
SWC4818-1SHL-SS	SWC4818-1SHR-SS

Right-hinged door shown

*Lazy Susan Add-On Suffix "LS"



Left Hinged	Right Hinged
SWC4818-1GHL-SS	SWC4818-1GHR-SS




Wall Cabinets — 48" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature three shelves as standard for 48" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".

Example: **SW4830-2SS-16** — Depth

Height # of Doors Width Door Style



FRAMED GLASS HINGED DOOR CABINETS



Width	Left & Right Hinged
18"	SW4818-1GHL-SS
24"	SW4824-1GHL-SS
18"	SW4818-1GHR-SS
24"	SW4824-1GHR-SS

Right-hinged door shown



Width	Item Number
30"	SW4830-2GH-SS
36"	SW4836-2GH-SS
42"	SW4842-2GH-SS
48"	SW4848-2GH-SS

30" and wider

FRAMED GLASS SLIDING DOOR CABINETS



Width	Item Number
30"	SW4830-2GSF-SS
36"	SW4836-2GSF-SS
42"	SW4842-2GSF-SS
48"	SW4848-2GSF-SS

30" and wider

GLASS SLIDING DOOR CABINETS



Width	Item Number
30"	SW4830-2GS-SS
36"	SW4836-2GS-SS
42"	SW4842-2GS-SS
48"	SW4848-2GS-SS

30" and wider

MECHANICAL TUBE STANDARD


Wall standards are the least expensive way to install shelving for a wall. It requires full proper blocking along the whole standard to assure they are able to hold the weight on the shelves.



Height	Item Number
24"	MTS1124-SS
30"	MTS1130-SS
36"	MTS1136-SS
42"	MTS1142-SS
48"	MTS1148-SS
72"	MTS1172-SS

Example: **MTS1136**

Mechanical Tub Standard Height



Aprons

Aprons have the options to be ordered in depths of 22" and 28". The aprons can have up to three drawers if the width allows. If end covers are required, use the suffix "EC". Actual depth of an apron is 22.627 and 28".

Example: **SA043028-001** — # of Drawers

Height Depth
 |
 +-----+
 | |
 +-----+
 Width

UTILITY TABLES



Width	Item Number
24"	SA042422-000-SS
30"	SA043022-000-SS
36"	SA043622-000-SS
42"	SA044222-000-SS
48"	SA044822-000-SS
54"	SA045422-000-SS
60"	SA046022-000-SS
66"	SA046622-000-SS
72"	SA047222-000-SS



Width	Item Number	Item Number
24"	SA042422-001-SS	SA042428-001-SS
30"	SA043022-001-SS	SA043028-001-SS
36"	SA043622-001-SS	SA043628-001-SS
42"	SA044222-001-SS	SA044228-001-SS



Width	Item Number	Item Number
42"	SA044222-002-SS	SA044228-002-SS
48"	SA044822-002-SS	SA044828-002-SS
54"	SA045422-002-SS	SA045428-002-SS
60"	SA046022-002-SS	SA046028-002-SS
66"	SA046622-002-SS	SA046628-002-SS
72"	SA047222-002-SS	SA047228-002-SS



Width	Item Number	Item Number
72"	SA047222-003-SS	SA047228-003-SS

Filler Panels

Base cabinet filler panels are available in 29", 32" and 35" high. Floor case filler panels are available in 84" high. They are all available in any width required by entering the dimension in lieu of the "XX" in the part number. Note: AMS base cabinets are 22.375" deep. (Example: SLF35-6.625)

C-FILLER



Item Number	Item Number	Item Number
SCF29-XX-SS	SCF32-XX-SS	SCF35-XX-SS

C-FRONT FILLER



Item Number	Item Number	Item Number
SFFC29-XX-SS	SFFC32-XX-SS	SFFC35-XX-SS

L-FILLER



Item Number	Item Number	Item Number
SLF29-XX-SS	SLF32-XX-SS	SLF35-XX-SS

L-FRONT FILLER




Item Number	Item Number	Item Number
SFFL29-XX-SS	SFFL32-XX-SS	SFFL35-XX-SS

Filler Panels

Base cabinet filler panels are available in 29", 32" and 35" high. Floor case filler panels are available in 84" high. They are all available in any width required by entering the dimension in lieu of the "XX" in the part number. Note: AMS base cabinets are 22.375" deep. (Example: SLF35-6.625)

CORNER FILLER



Item Number	Item Number	Item Number
SICF29-XX-SS	SICF32-XX-SS	SICF35-XX-SS

FLOOR CASE FILLER



Item Number
SFFL84-XX-SS

TOE SPACE FILLER



Item Number
STSF-XX-SS

WALL CASE FRONT FILLER



Item Number	Item Number	Item Number	Item Number
SWCFF24-XX-SS	SWCFF30-XX-SS	SWCFF36-XX-SS	SWCFF48-XX-SS

Wall cabinet front filler panels are available in 24", 30", 36" and 48" high. They are all available in any width required by entering the dimension in lieu of the "XX" in the part number. To specify the depth use the suffix "-13" or "-16".

WALL CASE CORNER FILLER



Item Number	Item Number	Item Number	Item Number
SWCCF24-SS	SWCCF30-SS	SWCCF36-SS	SWCCF48-SS

Knee Space Rails And Panels

KNEE SPACE RAILS



Item Number
SKSR0418-SS
SKSR0424-SS
SKSR0430-SS
SKSR0436-SS
SKSR0442-SS
SKSR0448-SS
SKSR0452-SS
SKSR0454-SS
SKSR0460-SS
SKSR0466-SS
SKSR0472-SS
SKSR0484-SS
SKSR0496-SS

KNEE SPACE PANELS



27" Height	29" Height	32" Height	35" Height
SKSP2724-SS	SKSP2924-SS	SKSP3224-SS	SKSP3524-SS
SKSP2730-SS	SKSP2930-SS	SKSP3230-SS	SKSP3530-SS
SKSP2736-SS	SKSP2936-SS	SKSP3236-SS	SKSP3536-SS
SKSP2742-SS	SKSP2942-SS	SKSP3242-SS	SKSP3542-SS
SKSP2748-SS	SKSP2948-SS	SKSP3248-SS	SKSP3548-SS
SKSP2754-SS	SKSP2954-SS	SKSP3254-SS	SKSP3554-SS
SKSP2760-SS	SKSP2960-SS	SKSP3260-SS	SKSP3560-SS
SKSP2766-SS	SKSP2966-SS	SKSP3266-SS	SKSP3566-SS
SKSP2772-SS	SKSP2972-SS	SKSP3272-SS	SKSP3572-SS



The suffix "-2P" represents two piece panels.

27" Height	29" Height	32" Height	35" Height
SKSP2724-2P-SS	SKSP2924-2P-SS	SKSP3224-2P-SS	SKSP3524-2P-SS
SKSP2730-2P-SS	SKSP2930-2P-SS	SKSP3230-2P-SS	SKSP3530-2P-SS
SKSP2736-2P-SS	SKSP2936-2P-SS	SKSP3236-2P-SS	SKSP3536-2P-SS
SKSP2742-2P-SS	SKSP2942-2P-SS	SKSP3242-2P-SS	SKSP3542-2P-SS
SKSP2748-2P-SS	SKSP2948-2P-SS	SKSP3248-2P-SS	SKSP3548-2P-SS
SKSP2754-2P-SS	SKSP2954-2P-SS	SKSP3254-2P-SS	SKSP3554-2P-SS
SKSP2760-2P-SS	SKSP2960-2P-SS	SKSP3260-2P-SS	SKSP3560-2P-SS
SKSP2766-2P-SS	SKSP2966-2P-SS	SKSP3266-2P-SS	SKSP3566-2P-SS
SKSP2772-2P-SS	SKSP2972-2P-SS	SKSP3272-2P-SS	SKSP3572-2P-SS



The suffix "-3P" represents three piece panels.

27" Height	29" Height	32" Height	35" Height
SKSP2724-3P-SS	SKSP2924-3P-SS	SKSP3224-3P-SS	SKSP3524-3P-SS
SKSP2730-3P-SS	SKSP2930-3P-SS	SKSP3230-3P-SS	SKSP3530-3P-SS
SKSP2736-3P-SS	SKSP2936-3P-SS	SKSP3236-3P-SS	SKSP3536-3P-SS
SKSP2742-3P-SS	SKSP2942-3P-SS	SKSP3242-3P-SS	SKSP3542-3P-SS
SKSP2748-3P-SS	SKSP2948-3P-SS	SKSP3248-3P-SS	SKSP3548-3P-SS
SKSP2754-3P-SS	SKSP2954-3P-SS	SKSP3254-3P-SS	SKSP3554-3P-SS
SKSP2760-3P-SS	SKSP2960-3P-SS	SKSP3260-3P-SS	SKSP3560-3P-SS
SKSP2766-3P-SS	SKSP2966-3P-SS	SKSP3266-3P-SS	SKSP3566-3P-SS
SKSP2772-3P-SS	SKSP2972-3P-SS	SKSP3272-3P-SS	SKSP3572-3P-SS

Knee Space Rails And Panels

KNEE SPACE PANELS - APRON SIZING

To accommodate aprons with drawer housing, knee space panels must be 4 inches smaller than the cabinet height.



24" Height	25" Height	28" Height	31" Height
SKSP2424-SS	SKSP2524-SS	SKSP2824-SS	SKSP3124-SS
SKSP2430-SS	SKSP2530-SS	SKSP2830-SS	SKSP3130-SS
SKSP2436-SS	SKSP2536-SS	SKSP2836-SS	SKSP3136-SS
SKSP2442-SS	SKSP2542-SS	SKSP2842-SS	SKSP3142-SS
SKSP2448-SS	SKSP2548-SS	SKSP2848-SS	SKSP3148-SS
SKSP2454-SS	SKSP2554-SS	SKSP2854-SS	SKSP3154-SS
SKSP2460-SS	SKSP2560-SS	SKSP2860-SS	SKSP3160-SS
SKSP2466-SS	SKSP2566-SS	SKSP2866-SS	SKSP3166-SS
SKSP2472-SS	SKSP2572-SS	SKSP2872-SS	SKSP3172-SS



The suffix "-2P" represents two piece panels.

24" Height	25" Height	28" Height	31" Height
SKSP2424-2P-SS	SKSP2524-2P-SS	SKSP2824-2P-SS	SKSP3124-2P-SS
SKSP2430-2P-SS	SKSP2530-2P-SS	SKSP2830-2P-SS	SKSP3130-2P-SS
SKSP2436-2P-SS	SKSP2536-2P-SS	SKSP2836-2P-SS	SKSP3136-2P-SS
SKSP2442-2P-SS	SKSP2542-2P-SS	SKSP2842-2P-SS	SKSP3142-2P-SS
SKSP2448-2P-SS	SKSP2548-2P-SS	SKSP2848-2P-SS	SKSP3148-2P-SS
SKSP2454-2P-SS	SKSP2554-2P-SS	SKSP2854-2P-SS	SKSP3154-2P-SS
SKSP2460-2P-SS	SKSP2560-2P-SS	SKSP2860-2P-SS	SKSP3160-2P-SS
SKSP2466-2P-SS	SKSP2566-2P-SS	SKSP2866-2P-SS	SKSP3166-2P-SS
SKSP2472-2P-SS	SKSP2572-2P-SS	SKSP2872-2P-SS	SKSP3172-2P-SS



The suffix "-3P" represents three piece panels.

24" Height	25" Height	28" Height	31" Height
SKSP2424-3P-SS	SKSP2524-3P-SS	SKSP2824-3P-SS	SKSP3124-3P-SS
SKSP2430-3P-SS	SKSP2530-3P-SS	SKSP2830-3P-SS	SKSP3130-3P-SS
SKSP2436-3P-SS	SKSP2536-3P-SS	SKSP2836-3P-SS	SKSP3136-3P-SS
SKSP2442-3P-SS	SKSP2542-3P-SS	SKSP2842-3P-SS	SKSP3142-3P-SS
SKSP2448-3P-SS	SKSP2548-3P-SS	SKSP2848-3P-SS	SKSP3148-3P-SS
SKSP2454-3P-SS	SKSP2554-3P-SS	SKSP2854-3P-SS	SKSP3154-3P-SS
SKSP2460-3P-SS	SKSP2560-3P-SS	SKSP2860-3P-SS	SKSP3160-3P-SS
SKSP2466-3P-SS	SKSP2566-3P-SS	SKSP2866-3P-SS	SKSP3166-3P-SS
SKSP2472-3P-SS	SKSP2572-3P-SS	SKSP2872-3P-SS	SKSP3172-3P-SS

CABINET FINISHED BACKS

Cabinet finished back - toe kick



Cabinet finished back - no toe kick




Example: CFB2930-18-NTK-SS

height
depth
width
No Toe Kick
(toe kicks come standard on units)

Cabinet finished backs come standard with a toe kick. Finished backs can be applied to a multitude of cabinets. The example above represents the layout for the part number. The height, width, and depth must match the cabinet that the finished back is to be applied too. If you do not want a toe kick on your unit be sure to put the suffix "NTK" at the end of the part number.



Counter Island Reagent Assembly



Example: CIRA21218-2E-1D-1P

Depth
Height
Width
Electrical Cutout
Plumbing Cutout
D=Duplex
Q=Quad

Item Number
CIRA21218-SS



Accessories

END PEDESTALS



Item Number
SEP2902-22-SS
SEP2902-28-SS
SEP3502-22-SS
SEP3502-28-SS

Actual: 22.375" & 28.000"

TABLE LEGS



Item Number
SL29-22-SS
SL29-28-SS
SL35-22-SS
SL35-28-SS

Actual: 22.625"

SINGLE LEG



Item Number
SSL-29-SS
SSL-35-SS



SINGLE LEG WITH PLATE

Item Number
SSL-29P-SS
SSL-35P-SS

CORNER SLOPED TOP



Item Number
SST-Corner13-SS
SST-Corner16-SS
SST-Corner21-SS

SLOPED TOP WITH END CAPS



Item Number	Item Number	Item Number
SST13-18-SS	SST16-18-SS	SST21-18-SS
SST13-24-SS	SST16-24-SS	SST21-24-SS
SST13-30-SS	SST16-30-SS	SST21-30-SS
SST13-36-SS	SST16-36-SS	SST21-36-SS
SST13-42-SS	SST16-42-SS	SST21-42-SS
SST13-48-SS	SST16-48-SS	SST21-48-SS

*Slope height is 6"

SHELVES




Item Number	Item Number	Item Number	Item Number
SWS1813-SS	SWS1816-SS	SFS1816-SS	SFS1821-SS
SWS2413-SS	SWS2416-SS	SFS2416-SS	SFS2421-SS
SWS3013-SS	SWS3016-SS	SFS3016-SS	SFS3021-SS
SWS3613-SS	SWS3616-SS	SFS3616-SS	SFS3621-SS
SWS4213-SS	SWS4216-SS	SFS4216-SS	SFS4221-SS
SWS4813-SS	SWS4816-SS	SFS4816-SS	SFS4821-SS

Canopy Hoods

All stainless steel 16-gauge canopy hoods wall mount or island style are supplied with hanging rods and attached brackets with duct collar. Standard and custom sizes available.

ISLAND CANOPY HOODS



Item Number	Length	Depth	Height
ICH-3630-30	36"	30"	30"
ICH-4830-30	48"	30"	30"
ICH-6030-30	60"	30"	30"
ICH-7230-30	72"	30"	30"
ICH-9630-30	96"	30"	30"

WALL CANOPY HOODS




Item Number	Length	Depth	Height
WCH-3630-30	36"	30"	30"
WCH-4830-30	48"	30"	30"
WCH-6030-30	60"	30"	30"
WCH-7230-30	72"	30"	30"
WCH-9630-30	96"	30"	30"

Countertop Wall Brackets

12-gauge stainless steel wall brackets for countertops.

WALL BRACKETS




Item Number	Depth	Height
CTB-30	29.750"	18"
CTB-25	24.75"	18"
CTB-24	23.75"	18"
CTB-18	17.750"	18"

Countertops

Stainless steel countertops and tables are constructed from type 304 stainless steel, polished to a #4 finish. Type 316 stainless steel is available upon request for an additional charge. Countertops and tables are constructed from 16-gauge thick material and sink bowls from 14-gauge material.

Table frames are made using 2" square tubular construction with adjustable bullet-type feet.

COUNTERTOPS



Marine Edge



Mechanical Joint

- ☐ 14-Gauge Sink Bowls
- ☐ Drain Outlets with Basket Strainers
- ☐ (OPTIONAL) Wall Mount Brackets
- ☐ Wall Clips provided
- ☐ Pitches Available
- ☐ Partial & Full End Splashes
- ☐ Tops crated for shipping
- ☐ Price includes a submittal dwg
- ☐ Stainless Steel hat channel substrate
- ☐ Sound deadening compound applied





Stainless Steel Custom Countertops

Table Tops

Stainless steel counter tops are 1" thick 16-gauge with stainless steel hat channel substrate and sprayed with a sound deadening compound. Table tops are available in 24", 30", and 36" depths as standard. Customizations are available.

STAINLESS STEEL TABLE TOPS



	A	B	
Item Number	Length	Depth	Width
SCT-3024-1	30"	24"	1"
SCT-3624-1	36"	24"	1"
SCT-4224-1	42"	24"	1"
SCT-4824-1	48"	24"	1"
SCT-5424-1	54"	24"	1"
SCT-6024-1	60"	24"	1"
SCT-6624-1	66"	24"	1"
SCT-7224-1	72"	24"	1"
SCT-8424-1	84"	24"	1"
SCT-9624-1	96"	24"	1"

	A	B	
Item Number	Length	Depth	Width
SCT-3030-1	30"	30"	1"
SCT-3630-1	36"	30"	1"
SCT-4230-1	42"	30"	1"
SCT-4830-1	48"	30"	1"
SCT-5430-1	54"	30"	1"
SCT-6030-1	60"	30"	1"
SCT-6630-1	66"	30"	1"
SCT-7230-1	72"	30"	1"
SCT-8430-1	84"	30"	1"
SCT-9630-1	96"	30"	1"

	A	B	
Item Number	Length	Depth	Width
SCT-3036-1	30"	36"	1"
SCT-3636-1	36"	36"	1"
SCT-4236-1	42"	36"	1"
SCT-4836-1	48"	36"	1"
SCT-5436-1	54"	36"	1"
SCT-6036-1	60"	36"	1"
SCT-6636-1	66"	36"	1"
SCT-7236-1	72"	36"	1"
SCT-8436-1	84"	36"	1"
SCT-9636-1	96"	36"	1"

Wall Mounted Sink Enclosure

Stainless steel sink enclosures are supplied with 12-gauge end support brackets, an 18-gauge apron, and 18-gauge removable front panel.

WALL MOUNTED SINK ENCLOSURE



Item Number	Width	Depth	Height
SSE-30X23X22	30"	23"	22"
SSE-30X24X22	30"	24"	22"
SSE-30X29X22	30"	29"	22"
SSE-36X23X22	36"	23"	22"
SSE-36X24X22	36"	24"	22"
SSE-36X29X22	36"	29"	22"
SSE-48X23X22	48"	23"	22"
SSE-48X24X22	48"	24"	22"
SSE-48X29X22	48"	29"	22"

Scullery Sinks

Stainless steel 14-gauge scullery sinks come with 2" square tube legs on adjustable bullet feet with front to back and center stretchers. Tops are pitched for drainage. The "RFD" style comes with a lower shelf and boxed curb. Umbilical and column notches are available. Any size or configuration available. Comes with wall clips, drain outlets, and basket strainers.

SCULLERY SINKS



Item Number	Length	Depth	# of Sinks
SCUL-72-1B	72"	30"	1
SCUL-72-2B	72"	30"	2
SCUL-96-1B	96"	30"	2
SCUL-96-2B	96"	30"	3



Item Number	Length	Depth	# of Sinks
RFDS-72-1B	72"	33"	1
RFDS-72-2B	72"	33"	2
RFDS-96-2B	96"	33"	2
RFDS-96-3B	96"	33"	3
RFDS-120-2B	120"	33"	2
RFDS-120-3B	120"	33"	3

Cylinder Racks

CYLINDER RESTRAINT RACKS



SCR02-1-SS Cylinder Rack
Shown

Width	Item Number	# of Cylinders
16"	SCR01-1-SS	1
30"	SCR02-1-SS	2
44"	SCR03-1-SS	3
58"	SCR04-1-SS	4
72"	SCR05-1-SS	5

*Custom sizes available

*Includes chains

** Standard construction is welded.
Strut/Bolt version also available for a
value engineered option.



SCR02-2-SS Cylinder Rack
Shown

Width	Item Number	# of Cylinders
16"	SCR01-2-SS	2
30"	SCR02-2-SS	4
44"	SCR03-2-SS	6
58"	SCR04-2-SS	8
72"	SCR05-2-SS	10

*Custom sizes available

*Includes chains

** Standard construction is welded.
Strut/Bolt version also available for a
value engineered option.



SCR02-3-SS Cylinder Rack
Shown

Width	Item Number	# of Cylinders
16"	SCR01-3-SS	3
30"	SCR02-3-SS	6
44"	SCR03-3-SS	9
58"	SCR04-3-SS	12
72"	SCR05-3-SS	15


*Custom sizes available

*Includes chains

** Standard construction is welded.
Strut/Bolt version also available for a
value engineered option.

Umbilicals

ISLAND UMBILICALS




UMB-IR-22X5X72-SS Shown



Height	Item Number
68"	UMB-IR-X-X-68-SS
73"	UMB-IR-X-X-73-SS
78"	UMB-IR-X-X-78-SS
90"	UMB-IR-X-X-90-SS

WALL UMBILICALS



UMB-WR-22X6X78-SS Shown



Height	Item Number
68"	UMB-WR-X-X-68-SS
73"	UMB-WR-X-X-73-SS
78"	UMB-WR-X-X-78-SS
90"	UMB-WR-X-X-90-SS



Air Master Systems Corp.

6480 Norton Center Drive, Muskegon, MI 49441 • 231.798.1111 • 231.798.4000 fax • airmastersystems.com