



THE ELIMINATOR 100 SERIES: AIRFOIL FUME HOOD

APPLICATIONS: K-12; HIGHER EDUCATION; INDUSTRIAL; BIO TECH / R&D



SERVICE ROUGH-IN DETAIL
(FOR 48", 60", & 72" HOODS)



SERVICE ROUGH-IN DETAIL
(FOR 96" HOODS)

EPOXY WORK SURFACE
(SPECIFY LOCATION OF CUPSINK)

*When utilizing 22" deep cabinets, cup sink measurement changes to 5.25" in lieu of 9"

Note: At 5.25", it will have the cup sink sit slightly under the baffle.

☐ Include Acid Cabinet Vent Hole


CUPSINK OPTIONS:

NONE

☐ ADD 0.25" RAISED RIM

WORK SURFACE MATERIAL:

BLACK EPOXY (Standard)



PART #	SUPERSTRUCTURE			WORKTOP			HOOD OPENING		CFM @29"	CFM @18"	SP	BLOWER MODEL #	SHIPPING WEIGHT
	W	H	D	W	H	D	W	H					
EH-111-36	36"	56"	33.25"	36"	1.25"	30"	26"	29"	567	354	.25	BL-100-36	320
EH-111-48	48"	56"	33.25"	48"	1.25"	30"	38"	29"	833	521	.25	BL-100-48	385
EH-111-60	60"	56"	33.25"	60"	1.25"	30"	50"	29"	1100	688	.34	BL-100-60	450
EH-111-72	72"	56"	33.25"	72"	1.25"	30"	62"	29"	1367	854	.34	BL-100-72	510
EH-111-96	96"	56"	33.25"	96"	1.25"	30"	86"	29"	1900	1188	.40	BL-100-96	630