



THE ELIMINATOR 200 SERIES: DOUBLE FACED 100 SERIES

APPLICATIONS: K-12; HIGHER EDUCATION; INDUSTRIAL; BIO TECH / R&D



SERVICE ROUGH-IN DETAIL
(FOR 48", 60", & 72" HOODS)



SERVICE ROUGH-IN DETAIL
(FOR 96" HOODS)


Note: If Acid / Flammable cabinets are required - cupsink locations will vary, please contact AMS for details.

EPOXY WORK SURFACE
(SPECIFY LOCATION OF CUPSINK)

CUPSINK OPTIONS:

☐ ADD 0.25" RAISED RIM

WORK SURFACE MATERIAL:



PART #	SUPERSTRUCTURE			WORKTOP			HOOD OPENING		CFM @29"	CFM @18"	SP	BLOWER MODEL #	SHIPPING WEIGHT
	W	H	D	W	H	D	W	H					
EH-211-36	36"	56"	35.5"	36"	1.25"	29"	26"	29"	567	354	.25	BL-100-36	400
EH-211-48	48"	56"	35.5"	48"	1.25"	29"	38"	29"	833	521	.33	BL-100-48	450
EH-211-60	60"	56"	35.5"	60"	1.25"	29"	50"	29"	1100	688	.34	BL-100-60	500
EH-211-72	72"	56"	35.5"	72"	1.25"	29"	62"	29"	1367	854	.34	BL-100-72	550
EH-211-96	96"	56"	35.5"	96"	1.25"	29"	86"	29"	1900	1188	.40	BL-100-96	650