

THE ELIMINATOR 300 SERIES: THIN WALL


APPLICATIONS: K-MIDDLE SCHOOL

LINER MATERIAL:

POLY RESIN (standard)

BYPASS TYPE:

OPEN BYPASS (standard)



☐ 36" ☐ 48" ☐ 60" ☐ 72" ☐ 96"

Custom sizes available ; Contact AMS

SPECIAL REQUIREMENTS: _____

DECK MOUNT PLUMBING FIXTURES AVAILABLE

NONE

PLACEMENT OF ALARM

- ☐ AS SHOWN
☐ OPPOSITE SIDE

- ☐ BLOWER SWITCH (OPTIONAL)
☐ LIGHT SWITCH (OPTIONAL)

BASE CABINET OPTIONS:

NONE	___"W x ___"H x ___"D
NONE	___"W x ___"H x ___"D
NONE	___"W x ___"H x ___"D

NOTE: 18" DEEP FLAMMABLE AND/OR ACID CABINET IS
NEEDED IF DECK MOUNT FIXTURES ARE REQUIRED.

CABINET PULL OPTIONS:

N/A

WIRING DIAGRAM:

☐ WIRING DIAGRAM REQUIRED

STANDARD FEATURES

- **Interior Baffle Arrangement:**
Fixed Baffles (standard)
- **Exhaust Collar:**
Type 304 s/s - 10" Dia. Round (standard)
Other _____
- **Interior Lighting Fixture:**
Standard LED Fixtures (standard)
- **Sash Frame & Glass:**
Painted Steel Frame w/ Laminated
Safety Glass (standard)
- **Sash Type:**
Vertical (standard)

OPTIONAL ACCESSORIES

Lattice Assembly:

NONE

Interior Lighting Fixture:

- ☐ Explosion Proof Light

Sash Frame & Glass:

- ☐ Stainless Steel Frame w/ Laminated Safety Glass

Ceiling Enclosure:

- ☐ Yes (Height: _____) ☐ No
• **Hinged Access Panel on Ceiling Enclosures?:**
☐ Yes ☐ No

Finished Back:

- ☐ Yes ☐ No

Sash Type:

COMBINATION (Horizontal + Vertical)

Additional Options:

- ☐ Blower ☐ Sash Locks
☐ Sash Stops ☐ Sash Alarm
☐ Automatic Sash Return ☐ Sash Cap

Pass Through Cord Ports

- ☐ 1.5" Dia. Front Post Cord Ports
☐ ___" Dia. Side Panel Cord Ports (please specify size)

Service Fixture Type:

Special Plumbing Requirements

Pre-Wire + Electrical:

- ☐ Yes ☐ No


CUSTOMER: _____

QTY: _____


COLOR: _____

THE ELIMINATOR 300 SERIES: THIN WALL

APPLICATIONS: K-MIDDLE SCHOOL



SERVICE ROUGH-IN DETAIL
(FOR 36", 48", 60", & 72" HOODS)



SERVICE ROUGH-IN DETAIL
(FOR 96" HOODS)

EPOXY WORK SURFACE
(SPECIFY LOCATION OF CUPSINK)

- ☐ LEFT SIDE
- ☐ RIGHT SIDE
- ☐ Include Acid Cabinet Vent Hole


CUPSINK OPTIONS:

NONE

☐ ADD 0.25" RAISED RIM

WORK SURFACE MATERIAL:

BLACK EPOXY (Standard)



PART #	SUPERSTRUCTURE			WORKTOP			HOOD OPENING		CFM @29"	CFM @18"	SP	BLOWER MODEL #	SHIPPING WEIGHT
	W	H	D	W	H	D	W	H					
EH-311-36	36"	56"	29.5"	36"	1.25"	30"	26"	28"	646	431	.25	BL-300-36	190
EH-311-48	48"	56"	29.5"	48"	1.25"	30"	43"	28"	896	597	.25	BL-300-48	175
EH-311-60	60"	56"	29.5"	60"	1.25"	30"	55"	28"	1146	764	.39	BL-300-60	200
EH-311-72	72"	56"	29.5"	72"	1.25"	30"	67"	28"	1417	944	.37	BL-300-72	250
EH-311-96	96"	56"	29.5"	96"	1.25"	30"	91"	28"	1896	1264	.43	BL-300-96	275

*Height of open sash; CFM Calculated at 100 FPM