



THE ELIMINATOR 500 SERIES: ADA DESIGN

APPLICATIONS: K-12; UNIVERSITY; PRIVATE INDUSTRY



SERVICE ROUGH-IN DETAIL
(FOR 36", 48", 60", & 72" HOODS)



SERVICE ROUGH-IN DETAIL
(FOR 96" HOODS)

EPOXY WORK SURFACE
(SPECIFY LOCATION OF CUPSINK)

☐ Include Acid Cabinet Vent Hole


CUPSINK OPTIONS:

NONE

☐ ADD 0.25" RAISED RIM

WORK SURFACE MATERIAL:

BLACK EPOXY (Standard)



PART #	SUPERSTRUCTURE			WORKTOP			HOOD OPENING		CFM @29"	CFM @18"	SP	BLOWER MODEL #	SHIPPING WEIGHT
	W	H	D	W	H	D	W	H					
EH-511-48	48"	56"	33.25"	48"	1.25"	30"	38"	29"	833	521	.25	BL-500-48	395
EH-511-60	60"	56"	33.25"	60"	1.25"	30"	50"	29"	1100	688	.34	BL-500-60	450
EH-511-72	72"	56"	33.25"	72"	1.25"	30"	62"	29"	1367	854	.34	BL-500-72	525
EH-511-96	96"	56"	33.25"	96"	1.25"	30"	86"	29"	1900	1188	.40	BL-500-96	630