

THE ELIMINATOR 900 SERIES: DEMONSTRATION HOOD

APPLICATIONS: K-12

SPECIFY FITTINGS/
FIXTURES REQUIRED:

NONE

BLOWER SWITCH (OPTIONAL) ☐

LIGHT SWITCH (OPTIONAL) ☐


120V GFI DUPLEX (OPTIONAL) ☐

PLEASE SPECIFY
LEFT OR RIGHT:

☐ LEFT ☐ RIGHT
(Lefthand Shown)

LINER MATERIAL:

POLY RESIN (standard)



Custom sizes available ; Contact AMS

****NOTE: SERVICE FIXTURES MUST
BE DECK MOUNTED**

SPECIAL REQUIREMENTS: _____

****DECK MOUNT PLUMBING
FIXTURES AVAILABLE****

BASE CABINET OPTIONS:

NONE _____"W x _____"H x _____"D
NONE _____"W x _____"H x _____"D
NONE _____"W x _____"H x _____"D

CABINET PULL OPTIONS:

N/A

WIRING DIAGRAM

NO

STANDARD FEATURES

- **Interior Baffle Arrangement:**
Fixed Baffles (standard)
- **Exhaust Collar:**
Type 304 s/s - 10" Dia. Round (standard)
Other _____
- **Interior Lighting Fixture:**
Standard LED Fixtures (standard)
- **Sash Frame & Glass:**
Painted Steel Frame w/ Laminated Safety Glass (standard)
- **Sash Type:**
Vertical (standard)
Sash Interlock (standard)


OPTIONAL ACCESSORIES

- Lattice Assembly:
NONE
- Interior Lighting Fixture:
☐ Explosion Proof Light
- Sash Frame & Glass:
☐ Stainless Steel Frame w/ Laminated Safety Glass
- Ceiling Enclosure:
☐ Yes (Height: _____) ☐ No
• Hinged Access Panel on Ceiling Enclosures?:
☐ Yes ☐ No
- Finished Back:
☐ Yes ☐ No
- Sash Type:
COMBINATION (Horizontal + Vertical)
- Additional Options:
☐ Blower ☐ Sash Alarm
☐ Sash Stops ☐ Sash Cap
☐ Sash Locks
- Service Fixture (Deck Mount Type):

- Pre-Wire + Electrical:
☐ Yes ☐ No

THE ELIMINATOR 900 SERIES: DEMONSTRATION HOOD

APPLICATIONS: K-12



EPOXY WORK SURFACE
(SPECIFY LOCATION OF CUPSINK)

☐ LEFT SIDE ☐ RIGHT SIDE


CUPSINK OPTIONS:

NONE

☐ ADD 0.25" RAISED RIM

WORK SURFACE MATERIAL:

BLACK EPOXY (Standard)



PART #	SUPERSTRUCTURE			WORKTOP			HOOD OPENING		CFM @29"	CFM @18"	SP	BLOWER MODEL #	SHIPPING WEIGHT
	W	H	D	W	H	D	W	H					
EH-911-60	60"	48"	30"	60"	1.25"	32"	53"	26.5"	1030	722	.75	BL-900-60	500