




Electrical Accessories & Specifications

Pre-Wire service is completed in accordance with NEC and UL-1805 standards


Note: If the fume hood is not pre-wired, wiring in the field must be accomplished using UL listed electrical fixtures while observing NEC standards and local electrical codes.




**Single Pole
15 Amp Light Switch**
Standard Black Color
and Cover Plate
Part # LTSW




**Red Illuminated
Toggle Switch**
with Black Cover Plate
Part # BLSW




Explosion Proof Outlet
Class 1, Group 1
Part # XPO



Explosion Proof Switch
Class 1, Group 1
Part # XSW



**208 Volt, 20 Amp
Single Pole Receptacle**
230, 277 or other
available upon request
Part # 208



**120 Volt, 20 Amp Ground Fault
Interrupter (GFI)**
Duplex Outlet
Standard Black color with cover
Part # GFI



Explosion Proof Light
Class 1, Group 1
Available in
Incandescent or fluorescent
Part # XPL