



The **SOLUTION:** Fume Hood (6200 fpm)


APPLICATIONS: K-12; Higher Education; Industrial; Bio Tech / R&D



SERVICE ROUGH-IN DETAIL
(Specifically for 4', 5', & 6' Hoods)



SERVICE ROUGH-IN DETAIL
(Specifically for 8' Hoods)



EPOXY WORK SURFACE
(Specify Location of Cupsink)

NOTE: IF ACID / FLAMMABLE CABINETS ARE REQUIRE - CUPSINKS
LOCATIONS WILL VARY, PLEASE CONTACT AMS FOR DETAILS.

CUPSINK OPTION:

- ☐ 3" X 6" ☐ 3" X 9"
☐ OTHER _____
☐ NONE

Part No.	Superstructure			18"SASH HT.			29"SASH HT.		STATIC PRESSURE
	W	H	D	CFM @ 100 FPM	CFM @ 60 FPM	HOOD OPENING (SQ.FT)	CFM @ 100 FPM	CFM @ 60 FPM	
EH-SOL-48	48"	57.25"	33-3/4"	475	285	7.39	739	433	.25"
EH-SOL-60	60"	57.25"	33-3/4"	625	375	9.72	972	583	.33"
EH-SOL-72	72"	57.25"	33-3/4"	775	465	12.06	1206	723	.50"
EH-SOL-96	96"	57.25"	33-3/4"	1075	645	16.72	1672	1003	.26"



AIR MASTER SYSTEMS CORP
6480 NORTON CENTER DR.
MUSKEGON, MI. 49441
TEL (231) 798-1111 FAX (231) 798-4000
www.airmastersystems.com

CUSTOMER

Qty: _____
Color: _____