

# Stainless Steel Laboratory Casework

product catalog



# Contents


|                                                    |           |
|----------------------------------------------------|-----------|
| <b>Cabinet Specifications</b>                      | <b>3</b>  |
| <b>Ordering Guide</b>                              | <b>4</b>  |
| <b>Options Menu</b>                                | <b>5</b>  |
| <b>Base Cabinets — Standing Height</b>             | <b>8</b>  |
| Open Cabinets                                      | 8         |
| Door Cabinets                                      | 8         |
| Drawer Cabinets                                    | 9         |
| Door & Drawer Cabinets                             | 12        |
| Shelves & Pull Out Shelves                         | 16        |
| Corner Cabinets                                    | 16        |
| Three or Five Equal Drawer Cabinet                 | 17        |
| Sink & Fume Hood Cabinets                          | 18        |
| <b>Base Cabinets — Sitting Height</b>              | <b>19</b> |
| Open Cabinets                                      | 19        |
| Door Cabinets                                      | 19        |
| Drawer Cabinets                                    | 20        |
| Door & Drawer Cabinets                             | 23        |
| Sink & Fume Hood Cabinets                          | 25        |
| Corner Cabinets                                    | 26        |
| <b>Base Cabinets — ADA Height</b>                  | <b>27</b> |
| Open Cabinets                                      | 27        |
| Door Cabinets                                      | 27        |
| Drawer Cabinets                                    | 28        |
| Door & Drawer Cabinets                             | 29        |
| Corner Cabinets                                    | 30        |
| Sink & Fume Hood Cabinets                          | 31        |
| <b>Base Cabinets — Corrosives &amp; Flammables</b> | <b>33</b> |
| Corrosives Storage Cabinets                        | 33        |
| Flammables Storage Cabinets                        | 36        |
| <b>Vacuum Pump Cabinets</b>                        | <b>38</b> |
| <b>Mobile Cabinets</b>                             | <b>40</b> |
| <b>Suspended Cabinets</b>                          | <b>42</b> |
| <b>Overlay Cabinets</b>                            | <b>44</b> |
| <b>Tall Cabinets</b>                               | <b>46</b> |
| <b>Wall Cabinets — 24" Height</b>                  | <b>48</b> |
| <b>Wall Cabinets — 30" Height</b>                  | <b>50</b> |
| <b>Wall Cabinets — 36" Height</b>                  | <b>52</b> |
| <b>Wall Cabinets — 48" Height</b>                  | <b>54</b> |
| <b>Aprons</b>                                      | <b>56</b> |
| <b>Utility Tables</b>                              | <b>57</b> |
| <b>Filler Panels</b>                               | <b>58</b> |
| <b>Knee Space Rails and Panels</b>                 | <b>60</b> |
| <b>Counter Island Reagent</b>                      | <b>62</b> |
| <b>Accessories</b>                                 | <b>63</b> |
| <b>Canopy Hoods</b>                                | <b>64</b> |
| <b>Countertop Wall Brackets</b>                    | <b>64</b> |
| <b>Countertops</b>                                 | <b>65</b> |
| <b>Table Tops</b>                                  | <b>66</b> |
| <b>Wall Mounted Sink Enclosure</b>                 | <b>67</b> |
| <b>Scullery Sinks</b>                              | <b>68</b> |
| <b>Cylinder Racks</b>                              | <b>69</b> |
| <b>Umbilicals</b>                                  | <b>70</b> |

**Note:**

Images in this catalog are for representation only. For complete product specifications, contact your AMS representative.

# Cabinet Specifications

---



All surfaces shall be Stainless Steel Type 304. Specific test data is available upon request.

## Gauges of base cabinet components

|                         |        |
|-------------------------|--------|
| Aprons                  | 18 ga. |
| Back Panels             | 18 ga. |
| Bottom Panels           | 18 ga. |
| Drawer/Door Outer Panel | 20 ga. |
| Door Inner Panel        | 20 ga. |
| Drawer Bodies           | 20 ga. |
| Legs - 2"sq. tube       | 16 ga. |
| Shelves                 | 18 ga. |
| Side Panels             | 18 ga. |
| Table Frames            | 18 ga. |
| Security Panels         | 20 ga. |

For complete specifications, request our full line of product specifications. Subject to change.

# Ordering Guide

| S                                                        |                                    | TYPE | HEIGHT | WIDTH | - | STYLE                            | - | OPTIONS |
|----------------------------------------------------------|------------------------------------|------|--------|-------|---|----------------------------------|---|---------|
| <b>Cabinet Type</b>                                      |                                    |      |        |       |   |                                  |   |         |
| B                                                        | Base (Inset)                       |      |        |       |   |                                  |   |         |
| F                                                        | Floor (Tall)                       |      |        |       |   |                                  |   |         |
| HB                                                       | Suspended Base                     |      |        |       |   |                                  |   |         |
| OB                                                       | Overlay Base                       |      |        |       |   |                                  |   |         |
| MB                                                       | Mobile Base                        |      |        |       |   |                                  |   |         |
| W                                                        | Wall                               |      |        |       |   |                                  |   |         |
| ASB                                                      | Acid Storage Base                  |      |        |       |   |                                  |   |         |
| FSB                                                      | Flammable Storage                  |      |        |       |   |                                  |   |         |
| <b>Base Cabinet Height</b>                               |                                    |      |        |       |   |                                  |   |         |
| 29"                                                      | Sitting Height                     |      |        |       |   |                                  |   |         |
| 32"                                                      | ADA Height                         |      |        |       |   |                                  |   |         |
| 35"                                                      | Standing Height                    |      |        |       |   |                                  |   |         |
| <b>Wall Cabinet Height</b>                               |                                    |      |        |       |   |                                  |   |         |
| 24" 30" 36" 48"                                          |                                    |      |        |       |   |                                  |   |         |
| <b>Cabinet Width</b>                                     |                                    |      |        |       |   |                                  |   |         |
| 12" 15" 18" 21" 24" 30" 36" 42" 48"                      |                                    |      |        |       |   |                                  |   |         |
| <b>Cabinet Style</b>                                     |                                    |      |        |       |   |                                  |   |         |
| Reference AMS catalog for door and drawer configurations |                                    |      |        |       |   |                                  |   |         |
| <b>Option Suffix (for base cabinet)</b>                  |                                    |      |        |       |   |                                  |   |         |
| -18                                                      | Depth(contact AMS for approval)    |      |        | -PS   |   | 1 Pull Out Shelf                 |   |         |
| -LK                                                      | Locks                              |      |        | -2PS  |   | 2 Pull Out Shelves               |   |         |
| -LTK                                                     | Toe Kick On Left Side Of Cabinet   |      |        | -3PS  |   | 3 Pull Out Shelves               |   |         |
| -RTK                                                     | Toe Kick On Right Side Of Cabinet  |      |        | -VK   |   | Vent Kit (only for acid cabinet) |   |         |
| -SS                                                      | Stainless Steel, Type 304          |      |        | -VD   |   | Vertical Divider                 |   |         |
| -SS316                                                   | Stainless Steel, Type 316          |      |        | -FB   |   | Finished Backs                   |   |         |
| -BB                                                      | Bread Board(selected models only)  |      |        | -EC   |   | End Covers                       |   |         |
|                                                          |                                    |      |        | -F    |   | Flush Pull                       |   |         |
| <b>Option Suffix (for wall cabinet)</b>                  |                                    |      |        |       |   |                                  |   |         |
| -16                                                      | 16" Depth                          |      |        | -TG   |   | Tempered Glass                   |   |         |
| -LK                                                      | Locks                              |      |        | -LG   |   | Laminate (Safety) Glass          |   |         |
| -SS                                                      | Stainless Steel, Type 304          |      |        | -THB  |   | Titan Hanging Bracket            |   |         |
| -SS316                                                   | Stainless Steel, Type 316          |      |        | -F    |   | Flush Pull                       |   |         |
| <b>Option Suffix (for mobile cabinet)</b>                |                                    |      |        |       |   |                                  |   |         |
| -18                                                      | Depth(contact AMS for approval)    |      |        | -T    |   | 3/4" THK. Metal Top              |   |         |
| -LK                                                      | Locks                              |      |        | -E    |   | 3/4" THK. Epoxy Top              |   |         |
| -SS                                                      | Stainless Steel, Type 304          |      |        | -P    |   | 3/4" THK. Phenolic Top           |   |         |
| -SS316                                                   | Stainless Steel, Type 316          |      |        | -2"   |   | Casters (in lieu of standard 3") |   |         |
| -BB                                                      | Bread Board (selected models only) |      |        | -F    |   | Flush Pull                       |   |         |
| <b>Option Suffix (for tall cabinet)</b>                  |                                    |      |        |       |   |                                  |   |         |
| -LK                                                      | Locks                              |      |        | -LG   |   | Laminate (Safety) Glass          |   |         |
| -SS                                                      | Stainless Steel, Type 304          |      |        | -WB   |   | Wardrobe Bracket                 |   |         |
| -SS316                                                   | Stainless Steel, Type 316          |      |        | -PS   |   | 1 Pull Out Shelves               |   |         |
| -TG                                                      | Tempered Glass                     |      |        | -F    |   | Flush Pull                       |   |         |



# Options Menu

---

## DOOR & DRAWER FRONTS

---



Inset Steel  
(Standard Default)



Flush Overlay  
Steel

## DRAWER SLIDES

---



☐ 100 lb. Full Extension (Standard Default)

☐ 100 lb. Full Extension Soft Close

☐ 200 lb. Full Extension

## HINGES

---



Inset Steel  
Stainless Steel 5-Knuckle



Flush Overlay  
Stainless Steel 5-Knuckle

# Options Menu

---

## MAGNETIC & ROLLER CATCH

---



Magnetic Catch  
(Standard Default)



Roller Catch

## LOCKS

---

Master 5-Pin  
(Standard Default)



## PULLS

---



☐ Stainless Steel Wire Pulls

# Options Menu

---

## LABEL HOLDERS

---



- ☐ None (Standard Default)
- ☐ One each door and drawer
- ☐ One per pair of doors and each drawer

## BASE UNIT SHELF

---



- ☐ One piece adjustable (Standard Default)  
*Extra Quantity Available*

## GLASS FRAMED DOORS

---



- ☐ 3/16" Clear Float Glass (Standard Default)
- ☐ 1/4" Laminated Safety Glass
- ☐ 1/4" Tempered Safety Glass

## SECURITY PANELS

---



- ☐ Without Security Panel (Standard Default)
- ☐ With Security Panel

## BREADBOARD

---



- ☐ Sitting Height variations only

## GLASS DOORS

---



- ☐ 1/4" Clear Float Glass (Standard Default)
- ☐ Laminated Glass
- ☐ Tempered Glass

# Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 23" deep. Units are available in various widths as shown. Door units feature a single shelf as standard. All open units have 2 shelves. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



## OPEN CABINETS



Up to 24" wide

| Width | Item Number   |
|-------|---------------|
| 12"   | SB3512-000-SS |
| 15"   | SB3515-000-SS |
| 18"   | SB3518-000-SS |
| 21"   | SB3521-000-SS |
| 24"   | SB3524-000-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-000-SS |
| 36"   | SB3536-000-SS |
| 42"   | SB3542-000-SS |
| 48"   | SB3548-000-SS |

## DOOR CABINETS



Example: SB3515-100R-SS

| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SB3512-100L-SS      |
| 15"   | SB3515-100L-SS      |
| 18"   | SB3518-100L-SS      |
| 21"   | SB3521-100L-SS      |
| 24"   | SB3524-100L-SS      |
| 12"   | SB3512-100R-SS      |
| 15"   | SB3515-100R-SS      |
| 18"   | SB3518-100R-SS      |
| 21"   | SB3521-100R-SS      |
| 24"   | SB3524-100R-SS      |

Note: L(left) or R(right) refers to door hinge side.



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-200-SS |
| 36"   | SB3536-200-SS |
| 42"   | SB3542-200-SS |
| 48"   | SB3548-200-SS |



Example: SB3530-100-LPNL-SS

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-100-LPNL-SS  |
| 36"   | SB3536-100-LPNL-SS  |
| 42"   | SB3542-100-LPNL-SS  |
| 48"   | SB3548-100-LPNL-SS  |
| 30"   | SB3530-100-RPNL-SS  |
| 36"   | SB3536-100-RPNL-SS  |
| 42"   | SB3542-100-RPNL-SS  |
| 48"   | SB3548-100-RPNL-SS  |

Note: L(left) or R(right) refers to door hinge side. PNL(panel) refers to panel next to door

# Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB3542-034-LK**

Height: 34.75"  
Doors: 2  
Half Drawers: 2  
Width: 42"  
Drawers: 4  
Lock: LK



## DRAWER CABINETS



| Width | Item Number   |
|-------|---------------|
| 12"   | SB3512-032-SS |
| 15"   | SB3515-032-SS |
| 18"   | SB3518-032-SS |
| 21"   | SB3521-032-SS |
| 24"   | SB3524-032-SS |

Up to 24" wide



| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-032-SS |
| 36"   | SB3536-032-SS |
| 42"   | SB3542-032-SS |
| 48"   | SB3548-032-SS |

30" and wider



| Width | Item Number   |
|-------|---------------|
| 12"   | SB3512-024-SS |
| 15"   | SB3515-024-SS |
| 18"   | SB3518-024-SS |
| 21"   | SB3521-024-SS |
| 24"   | SB3524-024-SS |

Up to 24" wide



| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-028-SS |
| 36"   | SB3536-028-SS |
| 42"   | SB3542-028-SS |
| 48"   | SB3548-028-SS |

30" and wider



| Width | Item Number   |
|-------|---------------|
| 12"   | SB3512-008-SS |
| 15"   | SB3515-008-SS |
| 18"   | SB3518-008-SS |
| 21"   | SB3521-008-SS |
| 24"   | SB3524-008-SS |

Up to 24" wide



| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-008-SS |
| 36"   | SB3536-008-SS |
| 42"   | SB3542-008-SS |
| 48"   | SB3548-008-SS |

30" and wider



| Width | Item Number   |
|-------|---------------|
| 12"   | SB3512-040-SS |
| 15"   | SB3515-040-SS |
| 18"   | SB3518-040-SS |
| 21"   | SB3521-040-SS |
| 24"   | SB3524-040-SS |

Up to 24" wide



| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-040-SS |
| 36"   | SB3536-040-SS |
| 42"   | SB3542-040-SS |
| 48"   | SB3548-040-SS |

30" and wider

# Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB3542-0(10)4-LK** — Lock

Height: 34.75"  
Doors: 0 (10)  
Half Drawers: 4  
Width: SB3542  
Drawers: 0

Example: **SB3542-034-1F-LK** — Lock

Height: 34.75"  
Doors: 034  
Half Drawers: 1F  
Width: SB3542  
Drawers: 0  
File Drawers: 34



## DRAWER CABINETS



| Width | Item Number      |
|-------|------------------|
| 12"   | SB3512-012-1F-SS |
| 15"   | SB3515-012-1F-SS |
| 18"   | SB3518-012-1F-SS |
| 21"   | SB3521-012-1F-SS |
| 24"   | SB3524-012-1F-SS |

Up to 24" wide



| Width | Item Number      |
|-------|------------------|
| 12"   | SB3512-000-2F-SS |
| 15"   | SB3515-000-2F-SS |
| 18"   | SB3518-000-2F-SS |
| 21"   | SB3521-000-2F-SS |
| 24"   | SB3524-000-2F-SS |

Up to 24" wide



| Width | Item Number      |
|-------|------------------|
| 12"   | SB3512-004-1F-SS |
| 15"   | SB3515-004-1F-SS |
| 18"   | SB3518-004-1F-SS |
| 21"   | SB3521-004-1F-SS |
| 24"   | SB3524-004-1F-SS |

Up to 24" wide



| Width | Item Number      |
|-------|------------------|
| 12"   | SB3512-020-1F-SS |
| 15"   | SB3515-020-1F-SS |
| 18"   | SB3518-020-1F-SS |
| 21"   | SB3521-020-1F-SS |
| 24"   | SB3524-020-1F-SS |

Up to 24" wide



30" and wider

| Width | Item Number      |
|-------|------------------|
| 30"   | SB3530-024-2F-SS |
| 36"   | SB3536-024-2F-SS |
| 42"   | SB3542-024-2F-SS |
| 48"   | SB3548-024-2F-SS |



30" and wider

| Width | Item Number      |
|-------|------------------|
| 30"   | SB3530-000-4F-SS |
| 36"   | SB3536-000-4F-SS |
| 42"   | SB3542-000-4F-SS |
| 48"   | SB3548-000-4F-SS |



30" and wider

| Width | Item Number      |
|-------|------------------|
| 30"   | SB3530-008-2F-SS |
| 36"   | SB3536-008-2F-SS |
| 42"   | SB3542-008-2F-SS |
| 48"   | SB3548-008-2F-SS |



30" and wider

| Width | Item Number      |
|-------|------------------|
| 30"   | SB3530-040-2F-SS |
| 36"   | SB3536-040-2F-SS |
| 42"   | SB3542-040-2F-SS |
| 48"   | SB3548-040-2F-SS |

# Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB3542-034-LK** — Lock

Height: 34.75"  
Doors: 2  
Half Drawers: 1  
Width: 42"  
Drawers: 3



## DRAWER CABINETS



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-050-SS |
| 36"   | SB3536-050-SS |
| 42"   | SB3542-050-SS |
| 48"   | SB3548-050-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-064-SS |
| 36"   | SB3536-064-SS |
| 42"   | SB3542-064-SS |
| 48"   | SB3548-064-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-060-SS |
| 36"   | SB3536-060-SS |
| 42"   | SB3542-060-SS |
| 48"   | SB3548-060-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-070-SS |
| 36"   | SB3536-070-SS |
| 42"   | SB3542-070-SS |
| 48"   | SB3548-070-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-080-SS |
| 36"   | SB3536-080-SS |
| 42"   | SB3542-080-SS |
| 48"   | SB3548-080-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-052-SS |
| 36"   | SB3536-052-SS |
| 42"   | SB3542-052-SS |
| 48"   | SB3548-052-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-044-SS |
| 36"   | SB3536-044-SS |
| 42"   | SB3542-044-SS |
| 48"   | SB3548-044-SS |




30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-034-SS |
| 36"   | SB3536-034-SS |
| 42"   | SB3542-034-SS |
| 48"   | SB3548-034-SS |



# Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



Example: **SB3542-00(16)-LK** — Lock  
 Height: 34.75"  
 Doors: 00 (0 doors)  
 Half Drawers: 16 (16 half drawers)  
 Width: 42"  
 Drawers: 00 (0 drawers)

Example: **SB3542-110R-LK** — Lock  
 Height: 34.75"  
 Doors: 110 (110 doors)  
 Half Drawers: 00 (0 half drawers)  
 Width: 42"  
 Drawers: 16 (16 drawers)  
 Right Hinged: R

## DRAWER CABINETS



Up to 24" wide

| Width | Item Number   |
|-------|---------------|
| 12"   | SB3512-016-SS |
| 15"   | SB3515-016-SS |
| 18"   | SB3518-016-SS |
| 21"   | SB3521-016-SS |
| 24"   | SB3524-016-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-016-SS |
| 36"   | SB3536-016-SS |
| 42"   | SB3542-016-SS |
| 48"   | SB3548-016-SS |



30" and wider

| Width | Item Number      |
|-------|------------------|
| 30"   | SB3530-00(16)-SS |
| 36"   | SB3536-00(16)-SS |
| 42"   | SB3542-00(16)-SS |
| 48"   | SB3548-00(16)-SS |

## DOOR & DRAWER CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SB3512-110L-SS      |
| 15"   | SB3515-110L-SS      |
| 18"   | SB3518-110L-SS      |
| 21"   | SB3521-110L-SS      |
| 24"   | SB3524-110L-SS      |
| 12"   | SB3512-110R-SS      |
| 15"   | SB3515-110R-SS      |
| 18"   | SB3518-110R-SS      |
| 21"   | SB3521-110R-SS      |
| 24"   | SB3524-110R-SS      |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SB3512-102L-SS      |
| 15"   | SB3515-102L-SS      |
| 18"   | SB3518-102L-SS      |
| 21"   | SB3521-102L-SS      |
| 24"   | SB3524-102L-SS      |
| 12"   | SB3512-102R-SS      |
| 15"   | SB3515-102R-SS      |
| 18"   | SB3518-102R-SS      |
| 21"   | SB3521-102R-SS      |
| 24"   | SB3524-102R-SS      |

# Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



Example: **SB3542-034-LK**

Height: 34.75"  
Doors: 034  
Half Drawers: LK  
Width: 42"  
Drawers: 034

## DOOR & DRAWER CABINETS



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-212-SS |
| 36"   | SB3536-212-SS |
| 42"   | SB3542-212-SS |
| 48"   | SB3548-212-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-220-SS |
| 36"   | SB3536-220-SS |
| 42"   | SB3542-220-SS |
| 48"   | SB3548-220-SS |



30" and wider

| Width | Item Number     |
|-------|-----------------|
| 30"   | SB3530-202-A-SS |
| 36"   | SB3536-202-A-SS |
| 42"   | SB3542-202-A-SS |
| 48"   | SB3548-202-A-SS |



30" and wider

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-202L-SS      |
| 36"   | SB3536-202L-SS      |
| 42"   | SB3542-202L-SS      |
| 48"   | SB3548-202L-SS      |
| 30"   | SB3530-202R-SS      |
| 36"   | SB3536-202R-SS      |
| 42"   | SB3542-202R-SS      |
| 48"   | SB3548-202R-SS      |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-204-SS |
| 36"   | SB3536-204-SS |
| 42"   | SB3542-204-SS |
| 48"   | SB3548-204-SS |



30" and wider

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-142L-SS      |
| 36"   | SB3536-142L-SS      |
| 42"   | SB3542-142L-SS      |
| 48"   | SB3548-142L-SS      |
| 30"   | SB3530-142R-SS      |
| 36"   | SB3536-142R-SS      |
| 42"   | SB3542-142R-SS      |
| 48"   | SB3548-142R-SS      |



30" and wider

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-118L-SS      |
| 36"   | SB3536-118L-SS      |
| 42"   | SB3542-118L-SS      |
| 48"   | SB3548-118L-SS      |
| 30"   | SB3530-118R-SS      |
| 36"   | SB3536-118R-SS      |
| 42"   | SB3542-118R-SS      |
| 48"   | SB3548-118R-SS      |



30" and wider

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-118-AL-SS    |
| 36"   | SB3536-118-AL-SS    |
| 42"   | SB3542-118-AL-SS    |
| 48"   | SB3548-118-AL-SS    |
| 30"   | SB3530-118-AR-SS    |
| 36"   | SB3536-118-AR-SS    |
| 42"   | SB3542-118-AR-SS    |
| 48"   | SB3548-118-AR-SS    |

# Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard.

Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB3542-034-LK** — Lock

Height: 34.75"  
Doors: 034  
Half Drawers: LK  
Width: 42"  
Drawers: 034



## DOOR & DRAWER CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-140L-SS      |
| 36"   | SB3536-140L-SS      |
| 42"   | SB3542-140L-SS      |
| 48"   | SB3548-140L-SS      |
| 30"   | SB3530-140R-SS      |
| 36"   | SB3536-140R-SS      |
| 42"   | SB3542-140R-SS      |
| 48"   | SB3548-140R-SS      |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-140-BL-SS    |
| 36"   | SB3536-140-BL-SS    |
| 42"   | SB3542-140-BL-SS    |
| 48"   | SB3548-140-BL-SS    |
| 30"   | SB3530-140-BR-SS    |
| 36"   | SB3536-140-BR-SS    |
| 42"   | SB3542-140-BR-SS    |
| 48"   | SB3548-140-BR-SS    |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-134L-SS      |
| 36"   | SB3536-134L-SS      |
| 42"   | SB3542-134L-SS      |
| 48"   | SB3548-134L-SS      |
| 30"   | SB3530-134R-SS      |
| 36"   | SB3536-134R-SS      |
| 42"   | SB3542-134R-SS      |
| 48"   | SB3548-134R-SS      |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-134-AL-SS    |
| 36"   | SB3536-134-AL-SS    |
| 42"   | SB3542-134-AL-SS    |
| 48"   | SB3548-134-AL-SS    |
| 30"   | SB3530-134-AR-SS    |
| 36"   | SB3536-134-AR-SS    |
| 42"   | SB3542-134-AR-SS    |
| 48"   | SB3548-134-AR-SS    |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-210-SS |
| 36"   | SB3536-210-SS |
| 42"   | SB3542-210-SS |
| 48"   | SB3548-210-SS |




Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-210-AL-SS    |
| 36"   | SB3536-210-AL-SS    |
| 42"   | SB3542-210-AL-SS    |
| 48"   | SB3548-210-AL-SS    |
| 30"   | SB3530-210-AR-SS    |
| 36"   | SB3536-210-AR-SS    |
| 42"   | SB3542-210-AR-SS    |
| 48"   | SB3548-210-AR-SS    |

# Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a vertical divider is required between cup board and drawer section use the suffix "VD". If a finished back is required use the suffix "FB".



## DOOR & DRAWER CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-150L-SS      |
| 36"   | SB3536-150L-SS      |
| 42"   | SB3542-150L-SS      |
| 48"   | SB3548-150L-SS      |
| 30"   | SB3530-150R-SS      |
| 36"   | SB3536-150R-SS      |
| 42"   | SB3542-150R-SS      |
| 48"   | SB3548-150R-SS      |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-116L-SS      |
| 36"   | SB3536-116L-SS      |
| 42"   | SB3542-116L-SS      |
| 48"   | SB3548-116L-SS      |
| 30"   | SB3530-116R-SS      |
| 36"   | SB3536-116R-SS      |
| 42"   | SB3542-116R-SS      |
| 48"   | SB3548-116R-SS      |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-126L-SS      |
| 36"   | SB3536-126L-SS      |
| 42"   | SB3542-126L-SS      |
| 48"   | SB3548-126L-SS      |
| 30"   | SB3530-126R-SS      |
| 36"   | SB3536-126R-SS      |
| 42"   | SB3542-126R-SS      |
| 48"   | SB3548-126R-SS      |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-126-AL-SS    |
| 36"   | SB3536-126-AL-SS    |
| 42"   | SB3542-126-AL-SS    |
| 48"   | SB3548-126-AL-SS    |
| 30"   | SB3530-126-AR-SS    |
| 36"   | SB3536-126-AR-SS    |
| 42"   | SB3542-126-AR-SS    |
| 48"   | SB3548-126-AR-SS    |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-108L-SS      |
| 36"   | SB3536-108L-SS      |
| 42"   | SB3542-108L-SS      |
| 48"   | SB3548-108L-SS      |
| 30"   | SB3530-108R-SS      |
| 36"   | SB3536-108R-SS      |
| 42"   | SB3542-108R-SS      |
| 48"   | SB3548-108R-SS      |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-132L-SS      |
| 36"   | SB3536-132L-SS      |
| 42"   | SB3542-132L-SS      |
| 48"   | SB3548-132L-SS      |
| 30"   | SB3530-132R-SS      |
| 36"   | SB3536-132R-SS      |
| 42"   | SB3542-132R-SS      |
| 48"   | SB3548-132R-SS      |

# Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB3514-100-SVR**

Height: 34.75"  
Doors: 100"  
Half Drawers: SVR  
Louver: Right  
Width: 14"  
Drawers: 100"  
Sink: SVR



## SHELVES & PULL OUT SHELVES



| Width | Item Number |
|-------|-------------|
| 12"   | SBS1222-SS  |
| 15"   | SBS1522-SS  |
| 18"   | SBS1822-SS  |
| 21"   | SBS2122-SS  |
| 24"   | SBS2422-SS  |
| 30"   | SBS3022-SS  |
| 36"   | SBS3622-SS  |
| 42"   | SBS4222-SS  |
| 48"   | SBS4822-SS  |



| Width | Item Number   |
|-------|---------------|
| 15"   | SBPOS-1522-SS |
| 18"   | SBPOS-1822-SS |
| 21"   | SBPOS-2122-SS |
| 24"   | SBPOS-2422-SS |
| 30"   | SBPOS-3022-SS |
| 36"   | SBPOS-3622-SS |
| 42"   | SBPOS-4222-SS |
| 48"   | SBPOS-4822-SS |

## CORNER CABINETS



Right-hinged door shown

| Left Hinged    | Right Hinged   |
|----------------|----------------|
| SC3514-100L-SS | SC3514-100R-SS |

Lazy Susan Add-On Suffix  
"LS"




Right-hinged door shown

| Left Hinged      | Right Hinged     |
|------------------|------------------|
| SC3514-100-SL-SS | SC3514-100-SR-SS |

| Left Hinged with Louvers | Right Hinged with Louvers |
|--------------------------|---------------------------|
| SC3514-100-SVL-SS        | SC3514-100-SVR-SS         |

Lazy Susan Add-On Suffix  
"LS"



# Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature a single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



## THREE OR FIVE EQUAL DRAWER CABINETS



Up to 24" wide

| Width | Item Number   |
|-------|---------------|
| 12"   | SB3512-03X-SS |
| 15"   | SB3515-03X-SS |
| 18"   | SB3518-03X-SS |
| 21"   | SB3521-03X-SS |
| 24"   | SB3524-03X-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-03X-SS |
| 36"   | SB3536-03X-SS |
| 42"   | SB3542-03X-SS |
| 48"   | SB3548-03X-SS |



Up to 24" wide

| Width | Item Number   |
|-------|---------------|
| 12"   | SB3512-05X-SS |
| 15"   | SB3515-05X-SS |
| 18"   | SB3518-05X-SS |
| 21"   | SB3521-05X-SS |
| 24"   | SB3524-05X-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-05X-SS |
| 36"   | SB3536-05X-SS |
| 42"   | SB3542-05X-SS |
| 48"   | SB3548-05X-SS |



30" and wider

| Width | Item Number      |
|-------|------------------|
| 30"   | SB3530-0(10)X-SS |
| 36"   | SB3536-0(10)X-SS |
| 42"   | SB3542-0(10)X-SS |
| 48"   | SB3548-0(10)X-SS |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB3530-15X-L-SS     |
| 36"   | SB3536-15X-L-SS     |
| 42"   | SB3542-15X-L-SS     |
| 48"   | SB3548-15X-L-SS     |
| 30"   | SB3530-15X-R-SS     |
| 36"   | SB3536-15X-R-SS     |
| 42"   | SB3542-15X-R-SS     |
| 48"   | SB3548-15X-R-SS     |



# Base Cabinets — Standing Height

Standing height base cabinets are 34.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB3552-200-SI** — Island

Height: 34.75"  
Doors: SB3552  
Half Drawers: 200  
Width: 52"  
Drawers: SI  
Sink



## SINK & FUME HOOD CABINETS



Right-hinged door shown  
Add "-SV(R/L)" for Vents

| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SB3512-100-SL-SS    |
| 15"   | SB3515-100-SL-SS    |
| 18"   | SB3518-100-SL-SS    |
| 21"   | SB3521-100-SL-SS    |
| 24"   | SB3524-100-SL-SS    |
| 12"   | SB3512-100-SR-SS    |
| 15"   | SB3515-100-SR-SS    |
| 18"   | SB3518-100-SR-SS    |
| 21"   | SB3521-100-SR-SS    |
| 24"   | SB3524-100-SR-SS    |



30" and wider

| Width | Item Number     |
|-------|-----------------|
| 30"   | SB3530-200-S-SS |
| 36"   | SB3536-200-S-SS |
| 42"   | SB3542-200-S-SS |
| 48"   | SB3548-200-S-SS |

Add "-SV" for Vents



Right-hinged door shown  
-FHD "Full Height Doors"


| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SB3512-100-S-FHD-SS |
| 15"   | SB3515-100-S-FHD-SS |
| 18"   | SB3518-100-S-FHD-SS |
| 21"   | SB3521-100-S-FHD-SS |
| 24"   | SB3524-100-S-FHD-SS |
| 12"   | SB3512-100-S-FHD-SS |
| 15"   | SB3515-100-S-FHD-SS |
| 18"   | SB3518-100-S-FHD-SS |
| 21"   | SB3521-100-S-FHD-SS |
| 24"   | SB3524-100-S-FHD-SS |



30" and wider

| Width | Item Number         |
|-------|---------------------|
| 30"   | SB3530-200-S-FHD-SS |
| 36"   | SB3536-200-S-FHD-SS |
| 42"   | SB3542-200-S-FHD-SS |
| 48"   | SB3548-200-S-FHD-SS |

-FHD "Full Height Doors"



Sink Support Assembly

| Item Number |
|-------------|
| SSS-C14-SS  |
| SSS-18-SS   |
| SSS-24-SS   |
| SSS-30-SS   |
| SSS-36-SS   |
| SSS-42-SS   |
| SSS-48-SS   |
| SSS-52-SS   |
| SSS-58-SS   |



| Width | Item Number     |
|-------|-----------------|
| 52"   | SB3552-200-S-SS |
| 58"   | SB3558-200-S-SS |


Available with vents:  
\*Change to -SV



| Width | Item Number      |
|-------|------------------|
| 52"   | SB3552-200-SI-SS |
| 58"   | SB3558-200-SI-SS |

# Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 23" deep. Units are available in various widths as shown. Door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



## OPEN CABINETS



Up to 24" wide

| Width | Item Number   |
|-------|---------------|
| 12"   | SB2912-000-SS |
| 15"   | SB2915-000-SS |
| 18"   | SB2918-000-SS |
| 21"   | SB2921-000-SS |
| 24"   | SB2924-000-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-000-SS |
| 36"   | SB2936-000-SS |
| 42"   | SB2942-000-SS |
| 48"   | SB2948-000-SS |

## DOOR CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SB2912-100L-SS      |
| 15"   | SB2915-100L-SS      |
| 18"   | SB2918-100L-SS      |
| 21"   | SB2921-100L-SS      |
| 24"   | SB2924-100L-SS      |
| 12"   | SB2912-100R-SS      |
| 15"   | SB2915-100R-SS      |
| 18"   | SB2918-100R-SS      |
| 21"   | SB2921-100R-SS      |
| 24"   | SB2924-100R-SS      |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-200-SS |
| 36"   | SB2936-200-SS |
| 42"   | SB2942-200-SS |
| 48"   | SB2948-200-SS |



Example: SB2930-100-LPNL

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB2930-100-LPNL-SS  |
| 36"   | SB2936-100-LPNL-SS  |
| 42"   | SB2942-100-LPNL-SS  |
| 48"   | SB2948-100-LPNL-SS  |
| 30"   | SB2930-100-RPNL-SS  |
| 36"   | SB2936-100-RPNL-SS  |
| 42"   | SB2942-100-RPNL-SS  |
| 48"   | SB2948-100-RPNL-SS  |

Note: L(left) or R(right) refers to door hinge side.  
PNL(panel) refers to panel next to door.




# Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB2930-022-LK**

Height: 28.75"  
Doors: 2  
Half Drawers: 0  
Width: 30"  
Drawers: 2  
Lock: LK



## DRAWER CABINETS



Up to 24" wide

| Width | Item Number   |
|-------|---------------|
| 12"   | SB2912-030-SS |
| 15"   | SB2915-030-SS |
| 18"   | SB2918-030-SS |
| 21"   | SB2921-030-SS |
| 24"   | SB2924-030-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-030-SS |
| 36"   | SB2936-030-SS |
| 42"   | SB2942-030-SS |
| 48"   | SB2948-030-SS |



Up to 24" wide

| Width | Item Number   |
|-------|---------------|
| 12"   | SB2912-022-SS |
| 15"   | SB2915-022-SS |
| 18"   | SB2918-022-SS |
| 21"   | SB2921-022-SS |
| 24"   | SB2924-022-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-022-SS |
| 36"   | SB2936-022-SS |
| 42"   | SB2942-022-SS |
| 48"   | SB2948-022-SS |



Up to 24" wide

| Width | Item Number   |
|-------|---------------|
| 12"   | SB2912-014-SS |
| 15"   | SB2915-014-SS |
| 18"   | SB2918-014-SS |
| 21"   | SB2921-014-SS |
| 24"   | SB2924-014-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-014-SS |
| 36"   | SB2936-014-SS |
| 42"   | SB2942-014-SS |
| 48"   | SB2948-014-SS |

# Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".

Example: **SB2930-040-2F-LK**

Height: 28.75"  
Doors: 2  
Half Drawers: 0  
Width: 30"  
Drawers: 4  
File Drawers: 0  
Lock: LK



## DRAWER CABINETS



Up to 24" wide

| Width | Item Number   |
|-------|---------------|
| 12"   | SB2912-006-SS |
| 15"   | SB2915-006-SS |
| 18"   | SB2918-006-SS |
| 21"   | SB2921-006-SS |
| 24"   | SB2924-006-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-006-SS |
| 36"   | SB2936-006-SS |
| 42"   | SB2942-006-SS |
| 48"   | SB2948-006-SS |



Up to 24" wide

| Width | Item Number      |
|-------|------------------|
| 12"   | SB2912-002-1F-SS |
| 15"   | SB2915-002-1F-SS |
| 18"   | SB2918-002-1F-SS |
| 21"   | SB2921-002-1F-SS |
| 24"   | SB2924-002-1F-SS |



30" and wider

| Width | Item Number      |
|-------|------------------|
| 30"   | SB2930-004-2F-SS |
| 36"   | SB2936-004-2F-SS |
| 42"   | SB2942-004-2F-SS |
| 48"   | SB2948-004-2F-SS |



Up to 24" wide

| Width | Item Number      |
|-------|------------------|
| 15"   | SB2915-000-2F-SS |
| 18"   | SB2918-000-2F-SS |
| 21"   | SB2921-000-2F-SS |
| 24"   | SB2924-000-2F-SS |



30" and wider

| Width | Item Number      |
|-------|------------------|
| 30"   | SB2930-000-4F-SS |
| 36"   | SB2936-000-4F-SS |
| 42"   | SB2942-000-4F-SS |
| 48"   | SB2948-000-4F-SS |



Up to 24" wide

| Width | Item Number      |
|-------|------------------|
| 15"   | SB2915-010-1F-SS |
| 18"   | SB2918-010-1F-SS |
| 21"   | SB2921-010-1F-SS |
| 24"   | SB2924-010-1F-SS |



30" and wider

| Width | Item Number      |
|-------|------------------|
| 30"   | SB2930-020-2F-SS |
| 36"   | SB2936-020-2F-SS |
| 42"   | SB2942-020-2F-SS |
| 48"   | SB2948-020-2F-SS |

# Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



## DRAWER CABINETS



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-060-SS |
| 36"   | SB2936-060-SS |
| 42"   | SB2942-060-SS |
| 48"   | SB2948-060-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-044-SS |
| 36"   | SB2936-044-SS |
| 42"   | SB2942-044-SS |
| 48"   | SB2948-044-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-040-SS |
| 36"   | SB2936-040-SS |
| 42"   | SB2942-040-SS |
| 48"   | SB2948-040-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-028-SS |
| 36"   | SB2936-028-SS |
| 42"   | SB2942-028-SS |
| 48"   | SB2948-028-SS |

# Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 23" deep. Units are available in various widths as shown. No shelves are available in sitting height base cabinets with the combination of doors and drawers. Drawer units feature 100 lb. full extension slides as standard. Specify left or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



Example: **SB2924-110R-LK** — Lock

Height: 28.75"  
Doors: 110"  
Half Drawers: 110"  
Width: 24"  
Drawers: 110"  
Right Hinged: R

## DOOR & DRAWER CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SB2912-110L-SS      |
| 15"   | SB2915-110L-SS      |
| 18"   | SB2918-110L-SS      |
| 21"   | SB2921-110L-SS      |
| 24"   | SB2924-110L-SS      |
| 12"   | SB2912-110R-SS      |
| 15"   | SB2915-110R-SS      |
| 18"   | SB2918-110R-SS      |
| 21"   | SB2921-110R-SS      |
| 24"   | SB2924-110R-SS      |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SB2912-102L-SS      |
| 15"   | SB2915-102L-SS      |
| 18"   | SB2918-102L-SS      |
| 21"   | SB2921-102L-SS      |
| 24"   | SB2924-102L-SS      |
| 12"   | SB2912-102R-SS      |
| 15"   | SB2915-102R-SS      |
| 18"   | SB2918-102R-SS      |
| 21"   | SB2921-102R-SS      |
| 24"   | SB2924-102R-SS      |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-220-SS |
| 36"   | SB2936-220-SS |
| 42"   | SB2942-220-SS |
| 48"   | SB2948-220-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-212-SS |
| 36"   | SB2936-212-SS |
| 42"   | SB2942-212-SS |
| 48"   | SB2948-212-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-204-SS |
| 36"   | SB2936-204-SS |
| 42"   | SB2942-204-SS |
| 48"   | SB2948-204-SS |




Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB2930-132L-SS      |
| 36"   | SB2936-132L-SS      |
| 42"   | SB2942-132L-SS      |
| 48"   | SB2948-132L-SS      |
| 30"   | SB2930-132R-SS      |
| 36"   | SB2936-132R-SS      |
| 42"   | SB2942-132R-SS      |
| 48"   | SB2948-132R-SS      |

# Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 23" deep. Units are available in various widths as shown. No shelves are available in sitting height base cabinets with the combination of doors and drawers. Drawer units feature 100 lb. full extension slides as standard. Door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



## DOOR & DRAWER CABINETS



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-210-SS |
| 36"   | SB2936-210-SS |
| 42"   | SB2942-210-SS |
| 48"   | SB2948-210-SS |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB2930-210-AL-SS    |
| 36"   | SB2936-210-AL-SS    |
| 42"   | SB2942-210-AL-SS    |
| 48"   | SB2948-210-AL-SS    |
| 30"   | SB2930-210-AR-SS    |
| 36"   | SB2936-210-AR-SS    |
| 42"   | SB2942-210-AR-SS    |
| 48"   | SB2948-210-AR-SS    |



30" and wider

| Width | Item Number     |
|-------|-----------------|
| 30"   | SB2930-202-A-SS |
| 36"   | SB2936-202-A-SS |
| 42"   | SB2942-202-A-SS |
| 48"   | SB2948-202-A-SS |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB2930-202L-SS      |
| 36"   | SB2936-202L-SS      |
| 42"   | SB2942-202L-SS      |
| 48"   | SB2948-202L-SS      |
| 30"   | SB2930-202R-SS      |
| 36"   | SB2936-202R-SS      |
| 42"   | SB2942-202R-SS      |
| 48"   | SB2948-202R-SS      |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB2930-130L-SS      |
| 36"   | SB2936-130L-SS      |
| 42"   | SB2942-130L-SS      |
| 48"   | SB2948-130L-SS      |
| 30"   | SB2930-130R-SS      |
| 36"   | SB2936-130R-SS      |
| 42"   | SB2942-130R-SS      |
| 48"   | SB2948-130R-SS      |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB2930-130-AL-SS    |
| 36"   | SB2936-130-AL-SS    |
| 42"   | SB2942-130-AL-SS    |
| 48"   | SB2948-130-AL-SS    |
| 30"   | SB2930-130-AR-SS    |
| 36"   | SB2936-130-AR-SS    |
| 42"   | SB2942-130-AR-SS    |
| 48"   | SB2948-130-AR-SS    |

# Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 23" deep. Units are available in various widths as shown. No shelves are available in sitting height base cabinets with the combination of doors and drawers. Drawer units feature 100 lb. full extension slides as standard. Door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a vertical divider is required between cup board and drawer section use the suffix "VD". If a finished back is required use the suffix "FB".



## DOOR & DRAWER CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB2930-116L-SS      |
| 36"   | SB2936-116L-SS      |
| 42"   | SB2942-116L-SS      |
| 48"   | SB2948-116L-SS      |
| 30"   | SB2930-116R-SS      |
| 36"   | SB2936-116R-SS      |
| 42"   | SB2942-116R-SS      |
| 48"   | SB2948-116R-SS      |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB2930-116-AL-SS    |
| 36"   | SB2936-116-AL-SS    |
| 42"   | SB2942-116-AL-SS    |
| 48"   | SB2948-116-AL-SS    |
| 30"   | SB2930-116-AR-SS    |
| 36"   | SB2936-116-AR-SS    |
| 42"   | SB2942-116-AR-SS    |
| 48"   | SB2948-116-AR-SS    |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB2930-124L-SS      |
| 36"   | SB2936-124L-SS      |
| 42"   | SB2942-124L-SS      |
| 48"   | SB2948-124L-SS      |
| 30"   | SB2930-124R-SS      |
| 36"   | SB2936-124R-SS      |
| 42"   | SB2942-124R-SS      |
| 48"   | SB2948-124R-SS      |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB2930-124-AL-SS    |
| 36"   | SB2936-124-AL-SS    |
| 42"   | SB2942-124-AL-SS    |
| 48"   | SB2948-124-AL-SS    |
| 30"   | SB2930-124-AR-SS    |
| 36"   | SB2936-124-AR-SS    |
| 42"   | SB2942-124-AR-SS    |
| 48"   | SB2948-124-AR-SS    |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB2930-114L-SS      |
| 36"   | SB2936-114L-SS      |
| 42"   | SB2942-114L-SS      |
| 48"   | SB2948-114L-SS      |
| 30"   | SB2930-114R-SS      |
| 36"   | SB2936-114R-SS      |
| 42"   | SB2942-114R-SS      |
| 48"   | SB2948-114R-SS      |



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 30"   | SB2930-108L-SS      |
| 36"   | SB2936-108L-SS      |
| 42"   | SB2942-108L-SS      |
| 48"   | SB2948-108L-SS      |
| 30"   | SB2930-108R-SS      |
| 36"   | SB2936-108R-SS      |
| 42"   | SB2942-108R-SS      |
| 48"   | SB2948-108R-SS      |

## SINK & FUME HOOD CABINETS



| Width | Left Hinged      | Right Hinged     |
|-------|------------------|------------------|
| 12"   | SB2912-100-SL-SS | SB2912-100-SR-SS |
| 15"   | SB2915-100-SL-SS | SB2915-100-SR-SS |
| 18"   | SB2918-100-SL-SS | SB2918-100-SR-SS |
| 21"   | SB2921-100-SL-SS | SB2921-100-SR-SS |
| 24"   | SB2924-100-SL-SS | SB2924-100-SR-SS |



| Width | Left Hinged          | Right Hinged         |
|-------|----------------------|----------------------|
| 12"   | SB2912-100-SL-FHD-SS | SB2912-100-SR-FHD-SS |
| 15"   | SB2915-100-SL-FHD-SS | SB2915-100-SR-FHD-SS |
| 18"   | SB2918-100-SL-FHD-SS | SB2918-100-SR-FHD-SS |
| 21"   | SB2921-100-SL-FHD-SS | SB2921-100-SR-FHD-SS |
| 24"   | SB2924-100-SL-FHD-SS | SB2924-100-SR-FHD-SS |

-FHD "Full Height Doors"

# Base Cabinets — Sitting Height

Sitting height cabinets are 28.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Full height door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



## SINK & FUME HOOD CABINETS




| Width | Item Number     |
|-------|-----------------|
| 30"   | SB2930-200-S-SS |
| 36"   | SB2936-200-S-SS |
| 42"   | SB2942-200-S-SS |
| 48"   | SB2948-200-S-SS |

Add "-SV" for Vents



| Width | Item Number         |
|-------|---------------------|
| 30"   | SB2930-200-S-FHD-SS |
| 36"   | SB2936-200-S-FHD-SS |
| 42"   | SB2942-200-S-FHD-SS |
| 48"   | SB2948-200-S-FHD-SS |

-FHD "Full Height Doors"



Sink Support Assembly

| Item Number |
|-------------|
| SSS-C14-SS  |
| SSS-18-SS   |
| SSS-24-SS   |
| SSS-30-SS   |
| SSS-36-SS   |
| SSS-42-SS   |
| SSS-48-SS   |
| SSS-52-SS   |
| SSS-58-SS   |



| Width | Item Number     |
|-------|-----------------|
| 52"   | SB2952-200-S-SS |
| 58"   | SB2958-200-S-SS |



| Width | Item Number      |
|-------|------------------|
| 52"   | SB2952-200-SI-SS |
| 58"   | SB2958-200-SI-SS |

## CORNER CABINETS



Right-hinged door shown

| Item Number         |
|---------------------|
| Left & Right Hinged |
| SC2914-100L-SS      |
| SC2914-100R-SS      |


Lazy Susan Add-On Suffix  
"LS"



Right-hinged door shown

| Left Hinged      | Right Hinged     |
|------------------|------------------|
| SC2914-100-SL-SS | SC2914-100-SR-SS |

| Left Hinged with Louvers | Right Hinged with Louvers |
|--------------------------|---------------------------|
| SC2914-100-SVL-SS        | SC2914-100-SVR-SS         |





# Base Cabinets — ADA Height

ADA height cabinets are 31.75" high and 23" deep. Units are available in various widths as shown. Door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



## OPEN CABINETS



Up to 24" wide

| Width | Item Number   |
|-------|---------------|
| 12"   | SB3212-000-SS |
| 15"   | SB3215-000-SS |
| 18"   | SB3218-000-SS |
| 21"   | SB3221-000-SS |
| 24"   | SB3224-000-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3230-000-SS |
| 36"   | SB3236-000-SS |
| 42"   | SB3242-000-SS |
| 48"   | SB3248-000-SS |

## DOOR CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SB3212-100L-SS      |
| 15"   | SB3215-100L-SS      |
| 18"   | SB3218-100L-SS      |
| 21"   | SB3221-100L-SS      |
| 24"   | SB3224-100L-SS      |
| 12"   | SB3212-100R-SS      |
| 15"   | SB3215-100R-SS      |
| 18"   | SB3218-100R-SS      |
| 21"   | SB3221-100R-SS      |
| 24"   | SB3224-100R-SS      |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3230-200-SS |
| 36"   | SB3236-200-SS |
| 42"   | SB3242-200-SS |
| 48"   | SB3248-200-SS |

# Base Cabinets — ADA Height

ADA height cabinets are 31.75" high and 23" deep. Units are available in various widths as shown. Full height door units feature single shelf as standard. Additional shelves are available. Drawer units feature 100 lb. full extension slides as standard. If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



## DRAWER CABINETS



| Width | Item Number   |
|-------|---------------|
| 12"   | SB3212-031-SS |
| 15"   | SB3215-031-SS |
| 18"   | SB3218-031-SS |
| 21"   | SB3221-031-SS |
| 24"   | SB3224-031-SS |

12" and wider



| Width | Item Number   |
|-------|---------------|
| 30"   | SB3230-031-SS |
| 36"   | SB3236-031-SS |
| 42"   | SB3242-031-SS |
| 48"   | SB3248-031-SS |

30" and wider



| Width | Item Number      |
|-------|------------------|
| 12"   | SB3212-020-1F-SS |
| 15"   | SB3215-020-1F-SS |
| 18"   | SB3218-020-1F-SS |
| 21"   | SB3221-020-1F-SS |
| 24"   | SB3224-020-1F-SS |

Up to 24" wide



| Width | Item Number      |
|-------|------------------|
| 30"   | SB3230-040-2F-SS |
| 36"   | SB3236-040-2F-SS |
| 42"   | SB3242-040-2F-SS |
| 48"   | SB3248-040-2F-SS |

30" and wider



| Width | Item Number      |
|-------|------------------|
| 15"   | SB3215-012-1F-SS |
| 18"   | SB3218-012-1F-SS |
| 21"   | SB3221-012-1F-SS |
| 24"   | SB3224-012-1F-SS |

Up to 24" wide



| Width | Item Number      |
|-------|------------------|
| 30"   | SB3230-024-2F-SS |
| 36"   | SB3236-024-2F-SS |
| 42"   | SB3242-024-2F-SS |
| 48"   | SB3248-024-2F-SS |

30" and wider



| Width | Item Number   |
|-------|---------------|
| 30"   | SB3230-062-SS |
| 36"   | SB3236-062-SS |
| 42"   | SB3242-062-SS |
| 48"   | SB3248-062-SS |

30" and wider

# Base Cabinets — ADA Height

ADA height cabinets are 31.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Full height door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



## DOOR & DRAWER CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SB3212-110L-SS      |
| 15"   | SB3215-110L-SS      |
| 18"   | SB3218-110L-SS      |
| 21"   | SB3221-110L-SS      |
| 24"   | SB3224-110L-SS      |
| 12"   | SB3212-110R-SS      |
| 15"   | SB3215-110R-SS      |
| 18"   | SB3218-110R-SS      |
| 21"   | SB3221-110R-SS      |
| 24"   | SB3224-110R-SS      |



30" and wider

| Width | Item Number     |
|-------|-----------------|
| 30"   | SB3230-202-A-SS |
| 36"   | SB3236-202-A-SS |
| 42"   | SB3242-202-A-SS |
| 48"   | SB3248-202-A-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3230-210-SS |
| 36"   | SB3236-210-SS |
| 42"   | SB3242-210-SS |
| 48"   | SB3248-210-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3230-201-SS |
| 36"   | SB3236-201-SS |
| 42"   | SB3242-201-SS |
| 48"   | SB3248-201-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3230-220-SS |
| 36"   | SB3236-220-SS |
| 42"   | SB3242-220-SS |
| 48"   | SB3248-220-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3230-202-SS |
| 36"   | SB3236-202-SS |
| 42"   | SB3242-202-SS |
| 48"   | SB3248-202-SS |

# Base Cabinets — ADA Height

ADA height cabinets are 31.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



## CORNER CABINETS



Right-hinged door shown

| Left Hinged       | Right Hinged      |
|-------------------|-------------------|
| SC3214-100L-LS-SS | SC3214-100R-LS-SS |

Lazy Susan Add-On Suffix  
"LS"


| Left Hinged    | Right Hinged   |
|----------------|----------------|
| SC3214-100L-SS | SC3214-100R-SS |



Right-hinged door shown

| Left Hinged      | Right Hinged     |
|------------------|------------------|
| SC3214-100-SL-SS | SC3214-100-SR-SS |

| Left Hinged<br>with Louvers | Right Hinged<br>with Louvers |
|-----------------------------|------------------------------|
| SC3214-100-SVL-SS           | SC3214-100-SVR-SS            |



# Base Cabinets — ADA Height

ADA height cabinets are 31.75" high and 23" deep. Units are available in various widths as shown. Full height door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".



## SINK & FUME HOOD CABINETS



Right-hinged door shown  
Add "-SV(R/L)" for Vents

| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SB3212-100-SL-SS    |
| 15"   | SB3215-100-SL-SS    |
| 18"   | SB3218-100-SL-SS    |
| 21"   | SB3221-100-SL-SS    |
| 24"   | SB3224-100-SL-SS    |
| 12"   | SB3224-100-SR-SS    |
| 15"   | SB3224-100-SR-SS    |
| 18"   | SB3224-100-SR-SS    |
| 21"   | SB3221-100-SR-SS    |
| 24"   | SB3224-100-SR-SS    |



Add "-SV" for Vents

| Width | Item Number     |
|-------|-----------------|
| 30"   | SB3230-200-S-SS |
| 36"   | SB3236-200-S-SS |
| 42"   | SB3242-200-S-SS |
| 48"   | SB3248-200-S-SS |

Add "-SV" for Vents



-FHD "Full Height Doors"

| Width | Left & Right Hinged  |
|-------|----------------------|
| 12"   | SB3212-100-SL-FHD-SS |
| 15"   | SB3215-100-SL-FHD-SS |
| 18"   | SB3218-100-SL-FHD-SS |
| 21"   | SB3221-100-SL-FHD-SS |
| 24"   | SB3224-100-SL-FHD-SS |
| 12"   | SB3224-100-SR-FHD-SS |
| 15"   | SB3224-100-SR-FHD-SS |
| 18"   | SB3224-100-SR-FHD-SS |
| 21"   | SB3221-100-SR-FHD-SS |
| 24"   | SB3224-100-SR-FHD-SS |



-FHD "Full Height Doors"

| Width | Item Number         |
|-------|---------------------|
| 30"   | SB3230-200-S-FHD-SS |
| 36"   | SB3236-200-S-FHD-SS |
| 42"   | SB3242-200-S-FHD-SS |
| 48"   | SB3248-200-S-FHD-SS |

-FHD "Full Height Doors"



| Width | Item Number     |
|-------|-----------------|
| 52"   | SB3252-200-S-SS |
| 58"   | SB3258-200-S-SS |



| Width | Item Number      |
|-------|------------------|
| 52"   | SB3252-200-SI-SS |
| 58"   | SB3258-200-SI-SS |

# Base Cabinets — ADA Height

ADA height cabinets are 31.75" high and 23" deep. Units are available in various widths as shown. Drawer units feature 100 lb. full extension slides as standard. Door units feature single shelf as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required, use the suffix "LK". If a finished back is required use the suffix "FB".




## SINK & FUME HOOD CABINETS



| Width | Item Number  |
|-------|--------------|
| 36"   | SADAKSC36-SS |
| 42"   | SADAKSC42-SS |
| 48"   | SADAKSC48-SS |



| Width | Item Number  |
|-------|--------------|
| 30"   | SADAKSP30-SS |
| 36"   | SADAKSP36-SS |
| 42"   | SADAKSP42-SS |
| 48"   | SADAKSP48-SS |



Sink Support Assembly

| Item Number |
|-------------|
| SSS-C14-SS  |
| SSS-18-SS   |
| SSS-24-SS   |
| SSS-30-SS   |
| SSS-36-SS   |
| SSS-42-SS   |
| SSS-48-SS   |
| SSS-52-SS   |
| SSS-58-SS   |


# Base Cabinets — Corrosives, Bases, Oxidizers

Corrosives storage cabinets are reinforced steel with a hi-tech performance polymer alloy coating, which provides 100% coverage against corrosion. Cabinets are available in 18" and 22" depths.

Example: **SASB3518-22R-VK** — Vent Kit

Height  
Depth  
Width  
Right Hinged

Removable Panel



## CORROSIVES STORAGE CABINETS




SASB3524-18-SS




SASB3530-18-SS

| Width | Sitting Height | ADA Height     | Standing Height |
|-------|----------------|----------------|-----------------|
| 18"   | SASB2918-18-SS | SASB3218-18-SS | SASB3518-18-SS  |
| 24"   | SASB2924-18-SS | SASB3224-18-SS | SASB3524-18-SS  |
| 30"   | SASB2930-18-SS | SASB3230-18-SS | SASB3530-18-SS  |
| 36"   | SASB2936-18-SS | SASB3236-18-SS | SASB3536-18-SS  |
| 42"   | SASB2942-18-SS | SASB3242-18-SS | SASB3542-18-SS  |
| 48"   | SASB2948-18-SS | SASB3248-18-SS | SASB3548-18-SS  |

| Width | Sitting Height | ADA Height     | Standing Height |
|-------|----------------|----------------|-----------------|
| 18"   | SASB2918-22-SS | SASB3218-22-SS | SASB3518-22-SS  |
| 24"   | SASB2924-22-SS | SASB3224-22-SS | SASB3524-22-SS  |
| 30"   | SASB2930-22-SS | SASB3230-22-SS | SASB3530-22-SS  |
| 36"   | SASB2936-22-SS | SASB3236-22-SS | SASB3536-22-SS  |
| 42"   | SASB2942-22-SS | SASB3242-22-SS | SASB3542-22-SS  |
| 48"   | SASB2948-22-SS | SASB3248-22-SS | SASB3548-22-SS  |



18-3/8" ACTUAL DEPTH



22-3/8" ACTUAL DEPTH




# Base Cabinets — Corrosives, Bases, Oxidizers

Corrosives storage cabinets are reinforced steel with a hi-tech performance polymer alloy coating, which provides 100% coverage against corrosion. Cabinets are available in 18" and 22" depths.

Example: **SAB3518-22R-VK** — Vent Kit


Height      Depth  
Width      Right Hinged

Removable Panel




## CORROSIVES STORAGE CABINETS

24" Wide (1 DOOR)



30" Wide (2 DOOR)



SASBM2930-22-SS

Mobile cabinets are available. Please contact your AMS representative.

| Width | Sitting Height  | ADA Height      | Standing Height | Width | Sitting Height  | ADA Height      | Standing Height |
|-------|-----------------|-----------------|-----------------|-------|-----------------|-----------------|-----------------|
| 18"   | SASBM2918-18-SS | SASBM3218-18-SS | SASBM3518-18-SS | 18"   | SASBM2918-22-SS | SASBM3218-22-SS | SASBM3518-22-SS |
| 24"   | SASBM2924-18-SS | SASBM3224-18-SS | SASBM3524-18-SS | 24"   | SASBM2924-22-SS | SASBM3224-22-SS | SASBM3524-22-SS |
| 30"   | SASBM2930-18-SS | SASBM3230-18-SS | SASBM3530-18-SS | 30"   | SASBM2930-22-SS | SASBM3230-22-SS | SASBM3530-22-SS |
| 36"   | SASBM2936-18-SS | SASBM3236-18-SS | SASBM3536-18-SS | 36"   | SASBM2936-22-SS | SASBM3236-22-SS | SASBM3536-22-SS |
| 42"   | SASBM2942-18-SS | SASBM3242-18-SS | SASBM3542-18-SS | 42"   | SASBM2942-22-SS | SASBM3242-22-SS | SASBM3542-22-SS |
| 48"   | SASBM2948-18-SS | SASBM3248-18-SS | SASBM3548-18-SS | 48"   | SASBM2948-22-SS | SASBM3248-22-SS | SASBM3548-22-SS |

When ordering 18" or 24" wide cabinets please place suffix "L" left or "R" right, after the part number for the hinge placement.


If a Vent Kit is required, use the suffix "VK".

**Note:** AMS recommends vent kits for all corrosives cabinets. Acids and bases should not be stored in the same cabinet.


# Corrosive - Vent Kit

Vent systems shall be manufactured from polyolefin resin and shall conform to ASTM #1412-92 for corrosion resistance and material content. The vent system will contain no PVC or solvents used for assembly. Vent may be constructed using rigid or flexible pipe as required. Vents will be manufactured by Scientific Plastics, Inc. of Kansas City, KS, USA.

## 403100-VENT KIT



## 403100-VENT KIT-L




# Base Cabinets — Flammables

Flammable storage cabinets are reinforced steel with a hi-tech performance polymer alloy coating, which provides 100% coverage against corrosion. Cabinets are available in 18" and 22" depths. All cabinets have 3 point locking systems.

Example: **SFSB3518-22R**

Height      Depth  
Width      Right Hinged

Removable Panel



## FLAMMABLES STORAGE CABINETS




SFSB2924-18-SS




SFSB2930-18-SS

| Width | Sitting Height | ADA Height     | Standing Height |
|-------|----------------|----------------|-----------------|
| 24"   | SFSB2924-18-SS | SFSB3224-18-SS | SFSB3524-18-SS  |
| 30"   | SFSB2930-18-SS | SFSB3230-18-SS | SFSB3530-18-SS  |
| 36"   | SFSB2936-18-SS | SFSB3236-18-SS | SFSB3536-18-SS  |
| 42"   | SFSB2942-18-SS | SFSB3242-18-SS | SFSB3542-18-SS  |
| 48"   | SFSB2948-18-SS | SFSB3248-18-SS | SFSB3548-18-SS  |

| Width | Sitting Height | ADA Height     | Standing Height |
|-------|----------------|----------------|-----------------|
| 24"   | SFSB2924-22-SS | SFSB3224-22-SS | SFSB3524-22-SS  |
| 30"   | SFSB2930-22-SS | SFSB3230-22-SS | SFSB3530-22-SS  |
| 36"   | SFSB2936-22-SS | SFSB3236-22-SS | SFSB3536-22-SS  |
| 42"   | SFSB2942-22-SS | SFSB3242-22-SS | SFSB3542-22-SS  |
| 48"   | SFSB2948-22-SS | SFSB3248-22-SS | SFSB3548-22-SS  |



**18-3/8" ACTUAL DEPTH**



**22 3/8" ACTUAL DEPTH**


# Base Cabinets — Flammables

Flammable storage cabinets are reinforced **stainless steel**. Cabinets are available in 18" and 22" depths.

All cabinets have 3 point locking systems.


Example: **SFSB3518-22R**

Height      Depth  
Width      Right Hinged




## FLAMMABLES STORAGE CABINETS

24" Wide (1 DOOR)



30" Wide (2 DOOR)



SFSBM3530-18-SS

Mobile cabinets are available. Please contact your AMS representative.

| Width | Sitting Height  | ADA Height      | Standing Height |
|-------|-----------------|-----------------|-----------------|
| 24"   | SFSBM2924-18-SS | SFSBM3224-18-SS | SFSBM3524-18-SS |
| 30"   | SFSBM2930-18-SS | SFSBM3230-18-SS | SFSBM3530-18-SS |
| 36"   | SFSBM2936-18-SS | SFSBM3236-18-SS | SFSBM3536-18-SS |
| 42"   | SFSBM2942-18-SS | SFSBM3242-18-SS | SFSBM3542-18-SS |
| 48"   | SFSBM2948-18-SS | SFSBM3248-18-SS | SFSBM3548-18-SS |

| Width | Sitting Height  | ADA Height      | Standing Height |
|-------|-----------------|-----------------|-----------------|
| 24"   | SFSBM2924-22-SS | SFSBM3224-22-SS | SFSBM3524-22-SS |
| 30"   | SFSBM2930-22-SS | SFSBM3230-22-SS | SFSBM3530-22-SS |
| 36"   | SFSBM2936-22-SS | SFSBM3236-22-SS | SFSBM3536-22-SS |
| 42"   | SFSBM2942-22-SS | SFSBM3242-22-SS | SFSBM3542-22-SS |
| 48"   | SFSBM2948-22-SS | SFSBM3248-22-SS | SFSBM3548-22-SS |

Flammables storage cabinets are certified to meet NFPA 45 and are UL listed. Cabinets are available in 18" and 22" depths. For 22" depth, add suffix "-22" to the part number. For synchronized self-closing doors add "-SCC" to the part number.

\*Cabinets with self-closing doors are constructed to the UL1275 standards, but cannot have a UL label

# Vacuum Pump Cabinets

Vacuum Pump Cabinets come with a Polydamp Melamine Foam that has heat resistance and sound dampening qualities. This is applied to the inside of the cabinet. (Please note actual cabinet depth with finished back will be 23" and 18.875").



| Width | Item Number   |
|-------|---------------|
| 15"   | SB3515-10V-SS |
| 18"   | SB3518-10V-SS |
| 21"   | SB3521-10V-SS |
| 24"   | SB3524-10V-SS |

When ordering 18" or 24" wide cabinets please place suffix "L" left or "R" right, after the part number for the hinge placement.


**Attention:** All of these units are shipped with a loose unassembled junction box, switch and wiring. These are to be installed by your certified electrician.



| Width | Item Number      |
|-------|------------------|
| 30"   | SB2930-20V-18-SS |
| 36"   | SB2936-20V-18-SS |
| 42"   | SB2942-20V-18-SS |
| 48"   | SB2948-20V-18-SS |

| Width | Item Number   |
|-------|---------------|
| 30"   | SB2930-20V-SS |
| 36"   | SB2936-20V-SS |
| 42"   | SB2942-20V-SS |
| 48"   | SB2948-20V-SS |

Pullout Shelf Dimensions for a SB3524-10V-22:



| Width | Item Number      |
|-------|------------------|
| 30"   | SB3230-20V-18-SS |
| 36"   | SB3236-20V-18-SS |
| 42"   | SB3242-20V-18-SS |
| 48"   | SB3248-20V-18-SS |

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3230-20V-SS |
| 36"   | SB3236-20V-SS |
| 42"   | SB3242-20V-SS |
| 48"   | SB3248-20V-SS |

| Full Cabinet Width | Pullout Shelf Width |
|--------------------|---------------------|
| 15"                | 7.4"                |
| 18"                | 10.4"               |
| 21"                | 13.4"               |
| 24"                | 16.4"               |
| 30"                | 22.4"               |
| 36"                | 28.4"               |
| 42"                | 34.4"               |
| 48"                | 40.4"               |



| Width | Item Number      |
|-------|------------------|
| 30"   | SB3530-20V-18-SS |
| 36"   | SB3536-20V-18-SS |
| 42"   | SB3542-20V-18-SS |
| 48"   | SB3548-20V-18-SS |

| Width | Item Number   |
|-------|---------------|
| 30"   | SB3530-20V-SS |
| 36"   | SB3536-20V-SS |
| 42"   | SB3542-20V-SS |
| 48"   | SB3548-20V-SS |

| Full Cabinet Depth | Pullout Shelf Depth |
|--------------------|---------------------|
| 18"                | 14.5"               |
| 22"                | 18.5"               |

# Vacuum Pump Cabinets

Vacuum Pump Cabinets come with a Polydamp Melamine Foam that has heat resistance and sound dampening qualities. This is applied to the inside of the cabinet. (Please note actual cabinet depth with finished back will be 23" and 18.875").



**Attention:** These units are shipped with a loose unassembled junction box, switch and wiring. These are to be installed by your certified electrician.



These units include the cart pictured as part of the unit price. Vent pipe is not included.

| Width | Item Number         |
|-------|---------------------|
| 30"   | SB2930-20V-18-NB-SS |
| 36"   | SB2936-20V-18-NB-SS |
| 42"   | SB2942-20V-18-NB-SS |
| 48"   | SB2948-20V-18-NB-SS |

| Width | Item Number      |
|-------|------------------|
| 30"   | SB2930-20V-NB-SS |
| 36"   | SB2936-20V-NB-SS |
| 42"   | SB2942-20V-NB-SS |
| 48"   | SB2948-20V-NB-SS |

| Width | Item Number         |
|-------|---------------------|
| 30"   | SB3230-20V-18-NB-SS |
| 36"   | SB3236-20V-18-NB-SS |
| 42"   | SB3242-20V-18-NB-SS |
| 48"   | SB3248-20V-18-NB-SS |

| Width | Item Number      |
|-------|------------------|
| 30"   | SB3230-20V-NB-SS |
| 36"   | SB3236-20V-NB-SS |
| 42"   | SB3242-20V-NB-SS |
| 48"   | SB3248-20V-NB-SS |

| Width | Item Number         |
|-------|---------------------|
| 30"   | SB3530-20V-18-NB-SS |
| 36"   | SB3536-20V-18-NB-SS |
| 42"   | SB3542-20V-18-NB-SS |
| 48"   | SB3548-20V-18-NB-SS |

| Width | Item Number      |
|-------|------------------|
| 30"   | SB3530-20V-NB-SS |
| 36"   | SB3536-20V-NB-SS |
| 42"   | SB3542-20V-NB-SS |
| 48"   | SB3548-20V-NB-SS |

| Width | Item Number       |
|-------|-------------------|
| 24"   | MOBILE-CART-24-SS |
| 36"   | MOBILE-CART-36-SS |

(Standard depth of 22")  
\*Actual Depth of cart: 17.186"

| Width | Item Number          |
|-------|----------------------|
| 24"   | MOBILE-CART-24-18-SS |
| 36"   | MOBILE-CART-36-18-SS |

(Suffix 18 represents 18" depth)  
\*Actual Depth of cart: 13.186"

# Mobile Cabinets

All AMS cabinets are available in a mobile version including low profile (27.18" high), sitting height (28.75" high), ADA height (31.75" high), and standing height (34.75" high). They come standard with a finished back and optional finished top. To specify a mobile cabinet use the prefix "M" to the base cabinet model number (example: SMB2724-110). They come standard with a 3" wheel caster. There is no height changes to the cabinets when replaced toe kicks with casters. 2" casters are available by adding the suffix "-2" to the part number. Cabinet depth with finished back is 22.812". All mobile cabinets with drawers come standard with a tip resistant counter weight.

**Tops:** Please add the suffix "T" for steel, "E" for epoxy, "P" for phenolic. All tops are 3/4" thick and with a 7/8" overhang on front only.  
(example SMB2724-110L-T)



# Mobile Cabinets

Standing height mobile cabinets are 34.75" high and 23" deep. Units are available in various widths. Each example below is available in the standard four heights shown on the left. Door units feature a single shelf as standard for standing height cabinets. All open units have 2 shelves. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK".



## DOOR & DRAWER CABINETS



| Width | Item Number     |
|-------|-----------------|
| 12"   | SMB3512-100R-SS |
| 15"   | SMB3515-100R-SS |
| 18"   | SMB3518-100R-SS |
| 21"   | SMB3521-100R-SS |
| 24"   | SMB3524-100R-SS |
| 12"   | SMB3512-100L-SS |
| 15"   | SMB3515-100L-SS |
| 18"   | SMB3518-100L-SS |
| 21"   | SMB3521-100L-SS |
| 24"   | SMB3524-100L-SS |



| Width | Item Number    |
|-------|----------------|
| 30"   | SMB3530-200-SS |
| 36"   | SMB3536-200-SS |
| 42"   | SMB3542-200-SS |
| 48"   | SMB3548-200-SS |



| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SMB3512-110L-SS     |
| 15"   | SMB3515-110L-SS     |
| 18"   | SMB3518-110L-SS     |
| 21"   | SMB3521-110L-SS     |
| 24"   | SMB3524-110L-SS     |
| 12"   | SMB3512-110R-SS     |
| 15"   | SMB3515-110R-SS     |
| 18"   | SMB3518-110R-SS     |
| 21"   | SMB3521-110R-SS     |
| 24"   | SMB3524-110R-SS     |

Note: L(left) or R(right) refers to door hinge side.



| Width | Item Number    |
|-------|----------------|
| 30"   | SMB3530-220-SS |
| 36"   | SMB3536-220-SS |
| 42"   | SMB3542-220-SS |
| 48"   | SMB3548-220-SS |




| Width | Item Number    |
|-------|----------------|
| 30"   | SMB3530-210-SS |
| 36"   | SMB3536-210-SS |
| 42"   | SMB3542-210-SS |
| 48"   | SMB3548-210-SS |




# Suspended Cabinets

All AMS cabinets are available in a suspended version including low profile (23.18" high), sitting height (24.875" high), ADA height (27.875" high), and standing height (30.875" high). To specify a suspended cabinet use the prefix "H" to the base cabinet model number (example: SHB2324-110). Finished backs are optional and can be specified with the suffix "-FB". Low profile cabinets should be used for adjustable height work benches or table frames on casters to allow the caster to swivel if the adjustable leg is positioned in the lowest setting.



# Suspended Cabinets

Standing height suspended cabinets are 30.875" high and 23" deep. Units are available in various widths. Each example below is available in the standard four heights shown on the left. Door units feature a single shelf as standard for standard height cabinets. All open units have 2 shelves. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK". If a finished back is required use the suffix "FB".



## DOOR & DRAWER CABINETS



| Width | Item Number     |
|-------|-----------------|
| 12"   | SHB3112-100R-SS |
| 15"   | SHB3115-100R-SS |
| 18"   | SHB3118-100R-SS |
| 21"   | SHB3121-100R-SS |
| 24"   | SHB3124-100R-SS |
| 12"   | SHB3112-100L-SS |
| 15"   | SHB3115-100L-SS |
| 18"   | SHB3118-100L-SS |
| 21"   | SHB3121-100L-SS |
| 24"   | SHB3124-100L-SS |



| Width | Item Number    |
|-------|----------------|
| 30"   | SHB3130-200-SS |
| 36"   | SHB3136-200-SS |
| 42"   | SHB3542-200-SS |
| 48"   | SHB3148-200-SS |



| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SHB3112-110L-SS     |
| 15"   | SHB3115-110L-SS     |
| 18"   | SHB3118-110L-SS     |
| 21"   | SHB3121-110L-SS     |
| 24"   | SHB3124-110L-SS     |
| 12"   | SHB3112-110R-SS     |
| 15"   | SHB3115-110R-SS     |
| 18"   | SHB3118-110R-SS     |
| 21"   | SHB3121-110R-SS     |
| 24"   | SHB3124-110R-SS     |

Note: L(left) or R(right) refers to door hinge side.



| Width | Item Number    |
|-------|----------------|
| 30"   | SHB3130-220-SS |
| 36"   | SHB3136-220-SS |
| 42"   | SHB3142-220-SS |
| 48"   | SHB3148-220-SS |




| Width | Item Number    |
|-------|----------------|
| 30"   | SHB3130-210-SS |
| 36"   | SHB3136-210-SS |
| 42"   | SHB3142-210-SS |
| 48"   | SHB3148-210-SS |

# Overlay Cabinets

All AMS base cabinets are available in an overlay version including low profile (27.18" high), sitting height (28.75" high), ADA height (31.75" high), and standing height (34.75" high). To specify an overlay cabinet, use the prefix "O" to the base cabinet model number (example: SOB2724-110). Overlay cabinets are also available in mobile and suspended variations.

**Tops:** Please add the suffix "T" for steel, "E" for epoxy, "P" for phenolic. All tops are 3/4" thick and with a 7/8" overhang on front only.  
(example SOB2724-110L-T)



# Overlay Cabinets

Standing height overlay base cabinets are 34.75" high and 23.875" deep (Including doors and drawers depth). Units are available in various widths. Each example below is available in the standard four heights shown on the left. Door units feature a single shelf as standard. All open units have 2 shelves. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R". If locks are required use the suffix "LK".



## DOOR & DRAWER CABINETS



| Width | Item Number     |
|-------|-----------------|
| 12"   | SOB3512-100R-SS |
| 15"   | SOB3515-100R-SS |
| 18"   | SOB3518-100R-SS |
| 21"   | SOB3521-100R-SS |
| 24"   | SOB3524-100R-SS |
| 12"   | SOB3512-100L-SS |
| 15"   | SOB3515-100L-SS |
| 18"   | SOB3518-100L-SS |
| 21"   | SOB3521-100L-SS |
| 24"   | SOB3524-100L-SS |



| Width | Item Number    |
|-------|----------------|
| 30"   | SOB3530-200-SS |
| 36"   | SOB3536-200-SS |
| 42"   | SOB3542-200-SS |
| 48"   | SOB3548-200-SS |



| Width | Left & Right Hinged |
|-------|---------------------|
| 12"   | SOB3512-110L-SS     |
| 15"   | SOB3515-110L-SS     |
| 18"   | SOB3518-110L-SS     |
| 21"   | SOB3521-110L-SS     |
| 24"   | SOB3524-110L-SS     |
| 12"   | SOB3512-110R-SS     |
| 15"   | SOB3515-110R-SS     |
| 18"   | SOB3518-110R-SS     |
| 21"   | SOB3521-110R-SS     |
| 24"   | SOB3524-110R-SS     |

Note: L(left) or R(right) refers to door hinge side.



| Width | Item Number    |
|-------|----------------|
| 30"   | SOB3530-220-SS |
| 36"   | SOB3536-220-SS |
| 42"   | SOB3542-220-SS |
| 48"   | SOB3548-220-SS |



| Width | Item Number    |
|-------|----------------|
| 30"   | SOB3530-210-SS |
| 36"   | SOB3536-210-SS |
| 42"   | SOB3542-210-SS |
| 48"   | SOB3548-210-SS |

# Tall Cabinets

Tall cabinets are 84" high and available in 16", 18", 21" and 24" deep by using the two digit suffix. Suffix is not needed when selecting 16" deep. Units are available in various widths as shown. Door units feature five shelves as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required use the suffix "LK".

Wardrobe cabinets are available by adding the suffix "W" to any of the standard part numbers. They will include one hanging bar in lieu of shelves.



Example: **SF8418-1GHR-18** — Depth

Labels for the example part number:  
 Height: 84  
 # of Doors: 1  
 Right Hinged: H  
 Width: 18  
 Door Style: G

## OPEN CABINETS



Up to 24" wide

| Width | Item Number   |
|-------|---------------|
| 18"   | SF8418-000-SS |
| 24"   | SF8424-000-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SF8430-000-SS |
| 36"   | SF8436-000-SS |
| 42"   | SF8442-000-SS |
| 48"   | SF8448-000-SS |

## DOOR CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 18"   | SF8418-1SHL-SS      |
| 24"   | SF8424-1SHL-SS      |
| 18"   | SF8418-1SHR-SS      |
| 24"   | SF8424-1SHR-SS      |



| Width | Item Number   |
|-------|---------------|
| 30"   | SF8430-2SH-SS |
| 36"   | SF8436-2SH-SS |
| 42"   | SF8442-2SH-SS |
| 48"   | SF8448-2SH-SS |



| Width | Item Number   |
|-------|---------------|
| 30"   | SF8430-2SS-SS |
| 36"   | SF8436-2SS-SS |
| 42"   | SF8442-2SS-SS |
| 48"   | SF8448-2SS-SS |

(this unit has sliding doors)



| Width | Item Number   |
|-------|---------------|
| 30"   | SF8430-4SH-SS |
| 36"   | SF8436-4SH-SS |
| 42"   | SF8442-4SH-SS |
| 48"   | SF8448-4SH-SS |

# Tall Cabinets

Tall cabinets are 84" high and available in 16", 18", 21" and 24" deep by using the two digit suffix. Suffix is not needed when selecting 16" deep. Units are available in various widths as shown. Door units feature five shelves as standard. Additional shelves are available. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required use the suffix "LK".

Wardrobe cabinets are available by adding the suffix "W" to any of the standard part numbers. They will include one hanging bar in lieu of shelves. For pull out shelves add the suffix "PS".



Example: **SF8418-1GHR-18** — Depth

Height   # of Doors   Right Hinged

Width   Door Style

## DOOR CABINETS



| Width | Left & Right Hinged |
|-------|---------------------|
| 18"   | SF8418-1GHL-SS      |
| 24"   | SF8424-1GHL-SS      |
| 18"   | SF8418-1GHR-SS      |
| 24"   | SF8424-1GHR-SS      |



| Width | Item Number   |
|-------|---------------|
| 30"   | SF8430-2GH-SS |
| 36"   | SF8436-2GH-SS |
| 42"   | SF8442-2GH-SS |
| 48"   | SF8448-2GH-SS |

Right-hinged door shown



| Width | Item Number    |
|-------|----------------|
| 30"   | SF8430-2GSF-SS |
| 36"   | SF8436-2GSF-SS |
| 42"   | SF8442-2GSF-SS |
| 48"   | SF8448-2GSF-SS |

(this unit has sliding doors)



| Width | Item Number   |
|-------|---------------|
| 30"   | SF8430-4GH-SS |
| 36"   | SF8436-4GH-SS |
| 42"   | SF8442-4GH-SS |
| 48"   | SF8448-4GH-SS |




| Width | Item Number      |
|-------|------------------|
| 30"   | SF8430-2GH2SH-SS |
| 36"   | SF8436-2GH2SH-SS |
| 42"   | SF8442-2GH2SH-SS |
| 48"   | SF8448-2GH2SH-SS |



Pull out shelves are available upon request.

# Wall Cabinets — 24" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature one shelf as standard for 24" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".



Example: **SW2430-2SS-16** — Depth

Height    # of Doors

Width    Door Style

## OPEN CABINETS



Up to 24" wide

| Width | Item Number   |
|-------|---------------|
| 18"   | SW2418-000-SS |
| 24"   | SW2424-000-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW2430-000-SS |
| 36"   | SW2436-000-SS |
| 42"   | SW2442-000-SS |
| 48"   | SW2448-000-SS |

## STEEL-HINGED DOOR CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 18"   | SW2418-1SHL-SS      |
| 24"   | SW2424-1SHL-SS      |
| 18"   | SW2418-1SHR-SS      |
| 24"   | SW2424-1SHR-SS      |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW2430-2SH-SS |
| 36"   | SW2436-2SH-SS |
| 42"   | SW2442-2SH-SS |
| 48"   | SW2448-2SH-SS |

## STEEL SLIDING DOOR CABINETS



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW2430-2SS-SS |
| 36"   | SW2436-2SS-SS |
| 42"   | SW2442-2SS-SS |
| 48"   | SW2448-2SS-SS |

## CORNER CABINETS



| Item Number    |
|----------------|
| SWC2418-000-SS |




Right-hinged door shown

| Left Hinged     | Right Hinged    |
|-----------------|-----------------|
| SWC2418-1SHL-SS | SWC2418-1SHR-SS |




| Left Hinged     | Right Hinged    |
|-----------------|-----------------|
| SWC2418-1GHL-SS | SWC2418-1GHR-SS |



\*Lazy Susan Add-On  
Suffix "LS"

# Wall Cabinets — 24" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature one shelf as standard for 24" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".



Example: **SW2430-2SS-16** — Depth

Height    # of Doors    Width    Door Style

## FRAMED GLASS HINGED DOOR CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 18"   | SW2418-1GHL-SS      |
| 24"   | SW2424-1GHL-SS      |
| 18"   | SW2418-1GHR-SS      |
| 24"   | SW2424-1GHR-SS      |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW2430-2GH-SS |
| 36"   | SW2436-2GH-SS |
| 42"   | SW2442-2GH-SS |
| 48"   | SW2448-2GH-SS |

## FRAMED GLASS SLIDING DOOR CABINETS



30" and wider

| Width | Item Number    |
|-------|----------------|
| 30"   | SW2430-2GSF-SS |
| 36"   | SW2436-2GSF-SS |
| 42"   | SW2442-2GSF-SS |
| 48"   | SW2448-2GSF-SS |

## GLASS SLIDING DOOR CABINETS



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW2430-2GS-SS |
| 36"   | SW2436-2GS-SS |
| 42"   | SW2442-2GS-SS |
| 48"   | SW2448-2GS-SS |



# Wall Cabinets — 30" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature two shelves as standard for 30" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".

Example: **SW3030-2SS-16** — Depth

Height    # of Doors

Width    Door Style



## OPEN CABINETS



| Width | Item Number   |
|-------|---------------|
| 18"   | SW3018-000-SS |
| 24"   | SW3024-000-SS |

Up to 24" wide



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW3030-000-SS |
| 36"   | SW3036-000-SS |
| 42"   | SW3042-000-SS |
| 48"   | SW3048-000-SS |

## STEEL-HINGED DOOR CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 18"   | SW3018-1SHL-SS      |
| 24"   | SW3024-1SHL-SS      |
| 18"   | SW3018-1SHR-SS      |
| 24"   | SW3024-1SHR-SS      |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW3030-2SH-SS |
| 36"   | SW3036-2SH-SS |
| 42"   | SW3042-2SH-SS |
| 48"   | SW3048-2SH-SS |

## STEEL SLIDING DOOR CABINETS



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW3030-2SS-SS |
| 36"   | SW3036-2SS-SS |
| 42"   | SW3042-2SS-SS |
| 48"   | SW3048-2SS-SS |

## CORNER CABINETS



| Item Number    |
|----------------|
| SWC3018-000-SS |




| Left Hinged     | Right Hinged    |
|-----------------|-----------------|
| SWC3018-1SHL-SS | SWC3018-1SHR-SS |

Right-hinged door shown

\*Lazy Susan Add-On Suffix "LS"



| Left Hinged     | Right Hinged    |
|-----------------|-----------------|
| SWC3018-1GHL-SS | SWC3018-1GHR-SS |




# Wall Cabinets — 30" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature two shelves as standard for 30" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".

Example: **SW3030-2SS-16**

Height    # of Doors    Depth  
Width    Door Style



## FRAMED GLASS HINGED DOOR CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 18"   | SW3018-1GHL-SS      |
| 24"   | SW3024-1GHL-SS      |
| 18"   | SW3018-1GHR-SS      |
| 24"   | SW3024-1GHR-SS      |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW3030-2GH-SS |
| 36"   | SW3036-2GH-SS |
| 42"   | SW3042-2GH-SS |
| 48"   | SW3048-2GH-SS |

## FRAMED GLASS SLIDING DOOR CABINETS



30" and wider

| Width | Item Number    |
|-------|----------------|
| 30"   | SW3030-2GSF-SS |
| 36"   | SW3036-2GSF-SS |
| 42"   | SW3042-2GSF-SS |
| 48"   | SW3048-2GSF-SS |

## GLASS SLIDING DOOR CABINETS



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW3030-2GS-SS |
| 36"   | SW3036-2GS-SS |
| 42"   | SW3042-2GS-SS |
| 48"   | SW3048-2GS-SS |

# Wall Cabinets — 36" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature two shelves as standard for 36" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".

Example: **SW3630-2SS-16**

Height      # of Doors      Depth  
Width      Door Style



## OPEN CABINETS



Up to 24" wide

| Width | Item Number   |
|-------|---------------|
| 18"   | SW3618-000-SS |
| 24"   | SW3624-000-SS |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW3630-000-SS |
| 36"   | SW3636-000-SS |
| 42"   | SW3642-000-SS |
| 48"   | SW3648-000-SS |

## STEEL-HINGED DOOR CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 18"   | SW3618-1SHL-SS      |
| 24"   | SW3624-1SHL-SS      |
| 18"   | SW3618-1SHR-SS      |
| 24"   | SW3624-1SHR-SS      |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW3630-2SH-SS |
| 36"   | SW3636-2SH-SS |
| 42"   | SW3642-2SH-SS |
| 48"   | SW3648-2SH-SS |

## STEEL SLIDING DOOR CABINETS



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW3630-2SS-SS |
| 36"   | SW3636-2SS-SS |
| 42"   | SW3642-2SS-SS |
| 48"   | SW3648-2SS-SS |

## CORNER CABINETS




| Item Number    |
|----------------|
| SWC3618-000-SS |



Right-hinged door shown

\*Lazy Susan Add-On Suffix "LS"



| Left Hinged     | Right Hinged    |
|-----------------|-----------------|
| SWC3618-1SHL-SS | SWC3618-1SHR-SS |

| Left Hinged     | Right Hinged    |
|-----------------|-----------------|
| SWC3618-1GHL-SS | SWC3618-1GHR-SS |

# Wall Cabinets — 36" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature two shelves as standard for 36" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".



Example: **SW3630-2SS-16** — Depth

Height      # of Doors  
Width      Door Style

## FRAMED GLASS HINGED DOOR CABINETS



Right-hinged door shown

| Width | Left & Right Hinged |
|-------|---------------------|
| 18"   | SW3618-1GHL-SS      |
| 24"   | SW3624-1GHL-SS      |
| 18"   | SW3618-1GHR-SS      |
| 24"   | SW3624-1GHR-SS      |



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW3630-2GH-SS |
| 36"   | SW3636-2GH-SS |
| 42"   | SW3642-2GH-SS |
| 48"   | SW3648-2GH-SS |

## FRAMED GLASS SLIDING DOOR CABINETS



30" and wider

| Width | Item Number    |
|-------|----------------|
| 30"   | SW3630-2GSF-SS |
| 36"   | SW3636-2GSF-SS |
| 42"   | SW3642-2GSF-SS |
| 48"   | SW3648-2GSF-SS |

## GLASS SLIDING DOOR CABINETS



30" and wider

| Width | Item Number   |
|-------|---------------|
| 30"   | SW3630-2GS-SS |
| 36"   | SW3636-2GS-SS |
| 42"   | SW3642-2GS-SS |
| 48"   | SW3648-2GS-SS |

# Wall Cabinets — 48" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature three shelves as standard for 48" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".

Example: **SW4830-2SS-16**

Height: 48"    # of Doors: 2    Width: 30"    Door Style: SS    Depth: 16"



## OPEN CABINETS



| Width | Item Number   |
|-------|---------------|
| 18"   | SW4818-000-SS |
| 24"   | SW4824-000-SS |

Up to 24" wide



| Width | Item Number   |
|-------|---------------|
| 30"   | SW4830-000-SS |
| 36"   | SW4836-000-SS |
| 42"   | SW4842-000-SS |
| 48"   | SW4848-000-SS |

30" and wider

## STEEL-HINGED DOOR CABINETS



| Width | Left & Right Hinged |
|-------|---------------------|
| 18"   | SW4818-1SHL-SS      |
| 24"   | SW4824-1SHL-SS      |
| 18"   | SW4818-1SHR-SS      |
| 24"   | SW4824-1SHR-SS      |

Right-hinged door shown



| Width | Item Number   |
|-------|---------------|
| 30"   | SW4830-2SH-SS |
| 36"   | SW4836-2SH-SS |
| 42"   | SW4842-2SH-SS |
| 48"   | SW4848-2SH-SS |

30" and wider

## STEEL SLIDING DOOR CABINETS



| Width | Item Number   |
|-------|---------------|
| 30"   | SW4830-2SS-SS |
| 36"   | SW4836-2SS-SS |
| 42"   | SW4842-2SS-SS |
| 48"   | SW4848-2SS-SS |

30" and wider

## CORNER CABINETS



| Item Number    |
|----------------|
| SWC4818-000-SS |




| Left Hinged     | Right Hinged    |
|-----------------|-----------------|
| SWC4818-1SHL-SS | SWC4818-1SHR-SS |

Right-hinged door shown

\*Lazy Susan Add-On Suffix "LS"



| Left Hinged     | Right Hinged    |
|-----------------|-----------------|
| SWC4818-1GHL-SS | SWC4818-1GHR-SS |




# Wall Cabinets — 48" Height

Wall cabinets are available in 24", 30", 36" and 48" heights and come standard with a 13" depth. For 16" deep cabinets, use the suffix "16". Units are available in various widths as shown. They feature three shelves as standard for 48" units. Specify left- or right-hinged doors using the suffix "L" or "R" after the two digit depth number. If locks are required, use the suffix "LK".

Example: **SW4830-2SS-16** — Depth

Height # of Doors Width Door Style



## FRAMED GLASS HINGED DOOR CABINETS



| Width | Left & Right Hinged |
|-------|---------------------|
| 18"   | SW4818-1GHL-SS      |
| 24"   | SW4824-1GHL-SS      |
| 18"   | SW4818-1GHR-SS      |
| 24"   | SW4824-1GHR-SS      |

Right-hinged door shown



| Width | Item Number   |
|-------|---------------|
| 30"   | SW4830-2GH-SS |
| 36"   | SW4836-2GH-SS |
| 42"   | SW4842-2GH-SS |
| 48"   | SW4848-2GH-SS |

30" and wider

## FRAMED GLASS SLIDING DOOR CABINETS



| Width | Item Number    |
|-------|----------------|
| 30"   | SW4830-2GSF-SS |
| 36"   | SW4836-2GSF-SS |
| 42"   | SW4842-2GSF-SS |
| 48"   | SW4848-2GSF-SS |

30" and wider

## GLASS SLIDING DOOR CABINETS



| Width | Item Number   |
|-------|---------------|
| 30"   | SW4830-2GS-SS |
| 36"   | SW4836-2GS-SS |
| 42"   | SW4842-2GS-SS |
| 48"   | SW4848-2GS-SS |

30" and wider

## MECHANICAL TUBE STANDARD


Wall standards are the least expensive way to install shelving for a wall. It requires full proper blocking along the whole standard to assure they are able to hold the weight on the shelves.



| Height | Item Number |
|--------|-------------|
| 24"    | MTS1124-SS  |
| 30"    | MTS1130-SS  |
| 36"    | MTS1136-SS  |
| 42"    | MTS1142-SS  |
| 48"    | MTS1148-SS  |
| 72"    | MTS1172-SS  |

Example: **MTS1136**

Mechanical Tub Standard Height




# Aprons


Aprons have the options to be ordered in depths of 22" and 28". The aprons can have up to three drawers if the width allows. If end covers are required, use the suffix "EC".



## UTILITY TABLES



| Width | Item Number     |
|-------|-----------------|
| 24"   | SA042422-000-SS |
| 30"   | SA043022-000-SS |
| 36"   | SA043622-000-SS |
| 42"   | SA044222-000-SS |
| 48"   | SA044822-000-SS |
| 54"   | SA045422-000-SS |
| 60"   | SA046022-000-SS |
| 66"   | SA046622-000-SS |
| 72"   | SA047222-000-SS |



| Width | Item Number     | Item Number     |
|-------|-----------------|-----------------|
| 24"   | SA042422-001-SS | SA042428-001-SS |
| 30"   | SA043022-001-SS | SA043028-001-SS |
| 36"   | SA043622-001-SS | SA043628-001-SS |
| 42"   | SA044222-001-SS | SA044228-001-SS |



| Width | Item Number     | Item Number     |
|-------|-----------------|-----------------|
| 42"   | SA044222-002-SS | SA044228-002-SS |
| 48"   | SA044822-002-SS | SA044828-002-SS |
| 54"   | SA045422-002-SS | SA045428-002-SS |
| 60"   | SA046022-002-SS | SA046028-002-SS |
| 66"   | SA046622-002-SS | SA046628-002-SS |
| 72"   | SA047222-002-SS | SA047228-002-SS |



| Width | Item Number     | Item Number     |
|-------|-----------------|-----------------|
| 72"   | SA047222-003-SS | SA047228-003-SS |




# Utility Tables

Utility Tables can be ordered in 29", 32", and 35" height. The frame Width (A) is 2" less than nominal size. The frame depth (B) can be 22"(22.625) deep and 28" deep. SUTXXXXXX-\_\_ The number of drawers are indicated by the last three numbers. The drawers are available as 000(No Drawers), 001(One Drawer), 002 (Two Drawers), and 003(Three Drawers). Utility Tables are shipped unassembled.


Example: **SUT293024-001** — # of Drawers

Height: 29"  
Depth: 30"  
Width: 24"


## UTILITY TABLES




| Width | Item Number      | (A) | (B)    | Item Number      | (A) | (B) |
|-------|------------------|-----|--------|------------------|-----|-----|
| 24"   | SUT292424-000-SS | 22  | 22.625 | SUT292430-000-SS | 22  | 28  |
| 30"   | SUT293024-000-SS | 28  | 22.625 | SUT293030-000-SS | 28  | 28  |
| 36"   | SUT293624-000-SS | 34  | 22.625 | SUT293630-000-SS | 34  | 28  |
| 42"   | SUT294224-000-SS | 40  | 22.625 | SUT294230-000-SS | 40  | 28  |
| 48"   | SUT294824-000-SS | 46  | 22.625 | SUT294830-000-SS | 46  | 28  |
| 54"   | SUT295424-000-SS | 52  | 22.625 | SUT295430-000-SS | 52  | 28  |
| 60"   | SUT296024-000-SS | 58  | 22.625 | SUT296030-000-SS | 58  | 28  |
| 66"   | SUT296624-000-SS | 64  | 22.625 | SUT296630-000-SS | 64  | 28  |
| 72"   | SUT297224-000-SS | 70  | 22.625 | SUT297230-000-SS | 70  | 28  |



| Width | Item Number      | (A) | (B)    | Item Number      | (A) | (B) |
|-------|------------------|-----|--------|------------------|-----|-----|
| 30"   | SUT293024-001-SS | 28  | 22.625 | SUT293030-001-SS | 28  | 28  |
| 36"   | SUT293624-001-SS | 34  | 22.625 | SUT293630-001-SS | 34  | 28  |
| 42"   | SUT294224-001-SS | 40  | 22.625 | SUT294230-001-SS | 40  | 28  |



| Width | Item Number      | (A) | (B)    | Item Number      | (A) | (B) |
|-------|------------------|-----|--------|------------------|-----|-----|
| 42"   | SUT294224-002-SS | 40  | 22.625 | SUT294230-002-SS | 40  | 28  |
| 48"   | SUT294824-002-SS | 46  | 22.625 | SUT294830-002-SS | 46  | 28  |
| 54"   | SUT295424-002-SS | 52  | 22.625 | SUT295430-002-SS | 52  | 28  |
| 60"   | SUT296024-002-SS | 58  | 22.625 | SUT296030-002-SS | 58  | 28  |
| 66"   | SUT296624-002-SS | 64  | 22.625 | SUT296630-002-SS | 64  | 28  |
| 72"   | SUT297224-002-SS | 70  | 22.625 | SUT297230-002-SS | 70  | 28  |



| Width | Item Number      | (A) | (B)    | Item Number      | (A) | (B) |
|-------|------------------|-----|--------|------------------|-----|-----|
| 72"   | SUT297224-003-SS | 70  | 22.625 | SUT297230-003-SS | 70  | 28  |



# Filler Panels

Base cabinet filler panels are available in 29", 32" and 35" high. Floor case filler panels are available in 84" high. They are all available in any width required by entering the dimension in lieu of the "XX" in the part number. Note: AMS base cabinets are 23" deep. (Example: SLF35-6.625)

## C-FILLER



| Item Number | Item Number | Item Number |
|-------------|-------------|-------------|
| SCF29-XX-SS | SCF32-XX-SS | SCF35-XX-SS |

## C-FRONT FILLER



| Item Number  | Item Number  | Item Number  |
|--------------|--------------|--------------|
| SFFC29-XX-SS | SFFC32-XX-SS | SFFC35-XX-SS |

## L-FILLER



| Item Number | Item Number | Item Number |
|-------------|-------------|-------------|
| SLF29-XX-SS | SLF32-XX-SS | SLF35-XX-SS |

## L-FRONT FILLER




| Item Number  | Item Number  | Item Number  |
|--------------|--------------|--------------|
| SFFL29-XX-SS | SFFL32-XX-SS | SFFL35-XX-SS |

# Filler Panels

Base cabinet filler panels are available in 29", 32" and 35" high. Floor case filler panels are available in 84" high. They are all available in any width required by entering the dimension in lieu of the "XX" in the part number. Note: AMS base cabinets are 23" deep. (Example: SLF35-6.625)

## CORNER FILLER



| Item Number  | Item Number  | Item Number  |
|--------------|--------------|--------------|
| SICF29-XX-SS | SICF32-XX-SS | SICF35-XX-SS |

## FLOOR CASE FILLER



| Item Number  |
|--------------|
| SFFL84-XX-SS |

## TOE SPACE FILLER



| Item Number |
|-------------|
| STSF-XX-SS  |

## WALL CASE FRONT FILLER



| Item Number   | Item Number   | Item Number   | Item Number   |
|---------------|---------------|---------------|---------------|
| SWCFF24-XX-SS | SWCFF30-XX-SS | SWCFF36-XX-SS | SWCFF48-XX-SS |

Wall cabinet front filler panels are available in 24", 30", 36" and 48" high. They are all available in any width required by entering the dimension in lieu of the "XX" in the part number. To specify the depth use the suffix "-13" or "-16".




## WALL CASE CORNER FILLER



| Item Number | Item Number | Item Number | Item Number |
|-------------|-------------|-------------|-------------|
| SWCCF24-SS  | SWCCF30-SS  | SWCCF36-SS  | SWCCF48-SS  |

# Knee Space Rails And Panels

## KNEE SPACE RAILS



| Item Number |
|-------------|
| SKSR0418-SS |
| SKSR0424-SS |
| SKSR0430-SS |
| SKSR0436-SS |
| SKSR0442-SS |
| SKSR0448-SS |
| SKSR0452-SS |
| SKSR0454-SS |
| SKSR0460-SS |
| SKSR0466-SS |
| SKSR0472-SS |
| SKSR0484-SS |
| SKSR0496-SS |

## KNEE SPACE PANELS



| 27" Height  | 29" Height  | 32" Height  | 35" Height  |
|-------------|-------------|-------------|-------------|
| SKSP2724-SS | SKSP2924-SS | SKSP3224-SS | SKSP3524-SS |
| SKSP2730-SS | SKSP2930-SS | SKSP3230-SS | SKSP3530-SS |
| SKSP2736-SS | SKSP2936-SS | SKSP3236-SS | SKSP3536-SS |
| SKSP2742-SS | SKSP2942-SS | SKSP3242-SS | SKSP3542-SS |
| SKSP2748-SS | SKSP2948-SS | SKSP3248-SS | SKSP3548-SS |
| SKSP2754-SS | SKSP2954-SS | SKSP3254-SS | SKSP3554-SS |
| SKSP2760-SS | SKSP2960-SS | SKSP3260-SS | SKSP3560-SS |
| SKSP2766-SS | SKSP2966-SS | SKSP3266-SS | SKSP3566-SS |
| SKSP2772-SS | SKSP2972-SS | SKSP3272-SS | SKSP3572-SS |



The suffix "-2P" represents two piece panels.

| 27" Height     | 29" Height     | 32" Height     | 35" Height     |
|----------------|----------------|----------------|----------------|
| SKSP2724-2P-SS | SKSP2924-2P-SS | SKSP3224-2P-SS | SKSP3524-2P-SS |
| SKSP2730-2P-SS | SKSP2930-2P-SS | SKSP3230-2P-SS | SKSP3530-2P-SS |
| SKSP2736-2P-SS | SKSP2936-2P-SS | SKSP3236-2P-SS | SKSP3536-2P-SS |
| SKSP2742-2P-SS | SKSP2942-2P-SS | SKSP3242-2P-SS | SKSP3542-2P-SS |
| SKSP2748-2P-SS | SKSP2948-2P-SS | SKSP3248-2P-SS | SKSP3548-2P-SS |
| SKSP2754-2P-SS | SKSP2954-2P-SS | SKSP3254-2P-SS | SKSP3554-2P-SS |
| SKSP2760-2P-SS | SKSP2960-2P-SS | SKSP3260-2P-SS | SKSP3560-2P-SS |
| SKSP2766-2P-SS | SKSP2966-2P-SS | SKSP3266-2P-SS | SKSP3566-2P-SS |
| SKSP2772-2P-SS | SKSP2972-2P-SS | SKSP3272-2P-SS | SKSP3572-2P-SS |



The suffix "-3P" represents three piece panels.

| 27" Height     | 29" Height     | 32" Height     | 35" Height     |
|----------------|----------------|----------------|----------------|
| SKSP2724-3P-SS | SKSP2924-3P-SS | SKSP3224-3P-SS | SKSP3524-3P-SS |
| SKSP2730-3P-SS | SKSP2930-3P-SS | SKSP3230-3P-SS | SKSP3530-3P-SS |
| SKSP2736-3P-SS | SKSP2936-3P-SS | SKSP3236-3P-SS | SKSP3536-3P-SS |
| SKSP2742-3P-SS | SKSP2942-3P-SS | SKSP3242-3P-SS | SKSP3542-3P-SS |
| SKSP2748-3P-SS | SKSP2948-3P-SS | SKSP3248-3P-SS | SKSP3548-3P-SS |
| SKSP2754-3P-SS | SKSP2954-3P-SS | SKSP3254-3P-SS | SKSP3554-3P-SS |
| SKSP2760-3P-SS | SKSP2960-3P-SS | SKSP3260-3P-SS | SKSP3560-3P-SS |
| SKSP2766-3P-SS | SKSP2966-3P-SS | SKSP3266-3P-SS | SKSP3566-3P-SS |
| SKSP2772-3P-SS | SKSP2972-3P-SS | SKSP3272-3P-SS | SKSP3572-3P-SS |

# Knee Space Rails And Panels

## KNEE SPACE PANELS - APRON SIZING

To accommodate aprons with drawer housing, knee space panels must be 4 inches smaller than the cabinet height.



| 24" Height  | 25" Height  | 28" Height  | 31" Height  |
|-------------|-------------|-------------|-------------|
| SKSP2424-SS | SKSP2524-SS | SKSP2824-SS | SKSP3124-SS |
| SKSP2430-SS | SKSP2530-SS | SKSP2830-SS | SKSP3130-SS |
| SKSP2436-SS | SKSP2536-SS | SKSP2836-SS | SKSP3136-SS |
| SKSP2442-SS | SKSP2542-SS | SKSP2842-SS | SKSP3142-SS |
| SKSP2448-SS | SKSP2548-SS | SKSP2848-SS | SKSP3148-SS |
| SKSP2454-SS | SKSP2554-SS | SKSP2854-SS | SKSP3154-SS |
| SKSP2460-SS | SKSP2560-SS | SKSP2860-SS | SKSP3160-SS |
| SKSP2466-SS | SKSP2566-SS | SKSP2866-SS | SKSP3166-SS |
| SKSP2472-SS | SKSP2572-SS | SKSP2872-SS | SKSP3172-SS |



The suffix "-2P" represents two piece panels.

| 24" Height     | 25" Height     | 28" Height     | 31" Height     |
|----------------|----------------|----------------|----------------|
| SKSP2424-2P-SS | SKSP2524-2P-SS | SKSP2824-2P-SS | SKSP3124-2P-SS |
| SKSP2430-2P-SS | SKSP2530-2P-SS | SKSP2830-2P-SS | SKSP3130-2P-SS |
| SKSP2436-2P-SS | SKSP2536-2P-SS | SKSP2836-2P-SS | SKSP3136-2P-SS |
| SKSP2442-2P-SS | SKSP2542-2P-SS | SKSP2842-2P-SS | SKSP3142-2P-SS |
| SKSP2448-2P-SS | SKSP2548-2P-SS | SKSP2848-2P-SS | SKSP3148-2P-SS |
| SKSP2454-2P-SS | SKSP2554-2P-SS | SKSP2854-2P-SS | SKSP3154-2P-SS |
| SKSP2460-2P-SS | SKSP2560-2P-SS | SKSP2860-2P-SS | SKSP3160-2P-SS |
| SKSP2466-2P-SS | SKSP2566-2P-SS | SKSP2866-2P-SS | SKSP3166-2P-SS |
| SKSP2472-2P-SS | SKSP2572-2P-SS | SKSP2872-2P-SS | SKSP3172-2P-SS |



The suffix "-3P" represents three piece panels.

| 24" Height     | 25" Height     | 28" Height     | 31" Height     |
|----------------|----------------|----------------|----------------|
| SKSP2424-3P-SS | SKSP2524-3P-SS | SKSP2824-3P-SS | SKSP3124-3P-SS |
| SKSP2430-3P-SS | SKSP2530-3P-SS | SKSP2830-3P-SS | SKSP3130-3P-SS |
| SKSP2436-3P-SS | SKSP2536-3P-SS | SKSP2836-3P-SS | SKSP3136-3P-SS |
| SKSP2442-3P-SS | SKSP2542-3P-SS | SKSP2842-3P-SS | SKSP3142-3P-SS |
| SKSP2448-3P-SS | SKSP2548-3P-SS | SKSP2848-3P-SS | SKSP3148-3P-SS |
| SKSP2454-3P-SS | SKSP2554-3P-SS | SKSP2854-3P-SS | SKSP3154-3P-SS |
| SKSP2460-3P-SS | SKSP2560-3P-SS | SKSP2860-3P-SS | SKSP3160-3P-SS |
| SKSP2466-3P-SS | SKSP2566-3P-SS | SKSP2866-3P-SS | SKSP3166-3P-SS |
| SKSP2472-3P-SS | SKSP2572-3P-SS | SKSP2872-3P-SS | SKSP3172-3P-SS |

# CABINET FINISHED BACKS

Cabinet finished back - toe kick



Cabinet finished back - no toe kick




Example: CFB2930-18-NTK-SS

height  
depth  
width  
No Toe Kick  
(toe kicks come standard on units)

Cabinet finished backs come standard with a toe kick. Finished backs can be applied to a multitude of cabinets. The example above represents the layout for the part number. The height, width, and depth must match the cabinet that the finished back is to be applied too. If you do not want a toe kick on your unit be sure to put the suffix "NTK" at the end of the part number.



# Counter Island Reagent Assembly



Example: CIRA21218-2E-1D-1P

Depth  
Height  
Width  
Electrical Cutout  
D=Duplex  
Q=Quad  
Plumbing Cutout

| Item Number  |
|--------------|
| CIRA21218-SS |



# Accessories

## END PEDESTALS



| Item Number   |
|---------------|
| SEP2902-22-SS |
| SEP2902-28-SS |
| SEP3502-22-SS |
| SEP3502-28-SS |

## TABLE LEGS



| Item Number |
|-------------|
| SL29-22-SS  |
| SL29-28-SS  |
| SL35-22-SS  |
| SL35-28-SS  |

## SINGLE LEG



| Item Number |
|-------------|
| SSL-29-SS   |
| SSL-35-SS   |



## SINGLE LEG WITH PLATE

| Item Number |
|-------------|
| SSL-29P-SS  |
| SSL-35P-SS  |

## CORNER SLOPED TOP



| Item Number     |
|-----------------|
| SST-Corner13-SS |
| SST-Corner16-SS |
| SST-Corner21-SS |

## SLOPED TOP WITH END CAPS



| Item Number | Item Number | Item Number |
|-------------|-------------|-------------|
| SST13-18-SS | SST16-18-SS | SST21-18-SS |
| SST13-24-SS | SST16-24-SS | SST21-24-SS |
| SST13-30-SS | SST16-30-SS | SST21-30-SS |
| SST13-36-SS | SST16-36-SS | SST21-36-SS |
| SST13-42-SS | SST16-42-SS | SST21-42-SS |
| SST13-48-SS | SST16-48-SS | SST21-48-SS |

## SHELVES



| Item Number | Item Number | Item Number | Item Number |
|-------------|-------------|-------------|-------------|
| SWS1813-SS  | SWS1816-SS  | SFS1816-SS  | SFS1821-SS  |
| SWS2413-SS  | SWS2416-SS  | SFS2416-SS  | SFS2421-SS  |
| SWS3013-SS  | SWS3016-SS  | SFS3016-SS  | SFS3021-SS  |
| SWS3613-SS  | SWS3616-SS  | SFS3616-SS  | SFS3621-SS  |
| SWS4213-SS  | SWS4216-SS  | SFS4216-SS  | SFS4221-SS  |
| SWS4813-SS  | SWS4816-SS  | SFS4816-SS  | SFS4821-SS  |

# Canopy Hoods

All stainless steel 16-gauge canopy hoods wall mount or island style are supplied with hanging rods and attached brackets with duct collar. Standard and custom sizes available.

## ISLAND CANOPY HOODS



| Item Number    | Length | Depth | Height |
|----------------|--------|-------|--------|
| ICH-3630-30-SS | 36"    | 30"   | 30"    |
| ICH-4830-30-SS | 48"    | 30"   | 30"    |
| ICH-6030-30-SS | 60"    | 30"   | 30"    |
| ICH-7230-30-SS | 72"    | 30"   | 30"    |
| ICH-9630-30-SS | 96"    | 30"   | 30"    |

## WALL CANOPY HOODS




| Item Number    | Length | Depth | Height |
|----------------|--------|-------|--------|
| WCH-3630-30-SS | 36"    | 30"   | 30"    |
| WCH-4830-30-SS | 48"    | 30"   | 30"    |
| WCH-6030-30-SS | 60"    | 30"   | 30"    |
| WCH-7230-30-SS | 72"    | 30"   | 30"    |
| WCH-9630-30-SS | 96"    | 30"   | 30"    |

# Countertop Wall Brackets

12-gauge stainless steel wall brackets for countertops.

## WALL BRACKETS



| Item Number | Depth   | Height |
|-------------|---------|--------|
| CTB-30-SS   | 29.750" | 18"    |
| CTB-25-SS   | 24.75"  | 18"    |
| CTB-24-SS   | 23.75"  | 18"    |
| CTB-18-SS   | 17.750" | 18"    |

# Countertops

---

Stainless steel countertops and tables are constructed from type 304 stainless steel, polished to a #4 finish. Type 316 stainless steel is available upon request for an additional charge. Countertops and tables are constructed from 16-gauge thick material and sink bowls from 14-gauge material.

Table frames are made using 2" square tubular construction with adjustable bullet-type feet.

## COUNTERTOPS

---



Marine Edge



Mechanical Joint

- ☐ 14-Gauge Sink Bowls
- ☐ Drain Outlets with Basket Strainers
- ☐ (OPTIONAL) Wall Mount Brackets
- ☐ Wall Clips provided
- ☐ Pitches Available
- ☐ Partial & Full End Splashes
- ☐ Tops crated for shipping
- ☐ Price includes a submittal dwg
- ☐ Stainless Steel hat channel substrate
- ☐ Sound deadening compound applied





Stainless Steel Custom Countertops



# Table Tops

Stainless steel counter tops are 1" thick 16-gauge with stainless steel hat channel substrate and sprayed with a sound deadening compound. Table tops are available in 24", 30", and 36" depths as standard. Customizations are available.

## STAINLESS STEEL TABLE TOPS



|               | A      | B     |       |
|---------------|--------|-------|-------|
| Item Number   | Length | Depth | Width |
| SCT-3024-1-SS | 30"    | 24"   | 1"    |
| SCT-3624-1-SS | 36"    | 24"   | 1"    |
| SCT-4224-1-SS | 42"    | 24"   | 1"    |
| SCT-4824-1-SS | 48"    | 24"   | 1"    |
| SCT-5424-1-SS | 54"    | 24"   | 1"    |
| SCT-6024-1-SS | 60"    | 24"   | 1"    |
| SCT-6624-1-SS | 66"    | 24"   | 1"    |
| SCT-7224-1-SS | 72"    | 24"   | 1"    |
| SCT-8424-1-SS | 84"    | 24"   | 1"    |
| SCT-9624-1-SS | 96"    | 24"   | 1"    |

|               | A      | B     |       |
|---------------|--------|-------|-------|
| Item Number   | Length | Depth | Width |
| SCT-3030-1-SS | 30"    | 30"   | 1"    |
| SCT-3630-1-SS | 36"    | 30"   | 1"    |
| SCT-4230-1-SS | 42"    | 30"   | 1"    |
| SCT-4830-1-SS | 48"    | 30"   | 1"    |
| SCT-5430-1-SS | 54"    | 30"   | 1"    |
| SCT-6030-1-SS | 60"    | 30"   | 1"    |
| SCT-6630-1-SS | 66"    | 30"   | 1"    |
| SCT-7230-1-SS | 72"    | 30"   | 1"    |
| SCT-8430-1-SS | 84"    | 30"   | 1"    |
| SCT-9630-1-SS | 96"    | 30"   | 1"    |

|               | A      | B     |       |
|---------------|--------|-------|-------|
| Item Number   | Length | Depth | Width |
| SCT-3036-1-SS | 30"    | 36"   | 1"    |
| SCT-3636-1-SS | 36"    | 36"   | 1"    |
| SCT-4236-1-SS | 42"    | 36"   | 1"    |
| SCT-4836-1-SS | 48"    | 36"   | 1"    |
| SCT-5436-1-SS | 54"    | 36"   | 1"    |
| SCT-6036-1-SS | 60"    | 36"   | 1"    |
| SCT-6636-1-SS | 66"    | 36"   | 1"    |
| SCT-7236-1-SS | 72"    | 36"   | 1"    |
| SCT-8436-1-SS | 84"    | 36"   | 1"    |
| SCT-9636-1-SS | 96"    | 36"   | 1"    |

# Wall Mounted Sink Enclosure

Stainless steel sink enclosures are supplied with 12-gauge end support brackets, an 18-gauge apron, and 18-gauge removable front panel.

## WALL MOUNTED SINK ENCLOSURE



| Item Number     | Width | Depth | Height |
|-----------------|-------|-------|--------|
| SSE-30X23X22-SS | 30"   | 23"   | 22"    |
| SSE-30X24X22-SS | 30"   | 24"   | 22"    |
| SSE-30X29X22-SS | 30"   | 29"   | 22"    |
| SSE-36X23X22-SS | 36"   | 23"   | 22"    |
| SSE-36X24X22-SS | 36"   | 24"   | 22"    |
| SSE-36X29X22-SS | 36"   | 29"   | 22"    |
| SSE-48X23X22-SS | 48"   | 23"   | 22"    |
| SSE-48X24X22-SS | 48"   | 24"   | 22"    |
| SSE-48X29X22-SS | 48"   | 29"   | 22"    |

# Scullery Sinks

Stainless steel 14-gauge scullery sinks come with 2" square tube legs on adjustable bullet feet with front to back and center stretchers. Tops are pitched for drainage. The "RFD" style comes with a lower shelf and boxed curb. Umbilical and column notches are available. Any size or configuration available. Comes with wall clips, drain outlets, and basket strainers.

## SCULLERY SINKS



| Item Number   | Length | Depth | # of Sinks |
|---------------|--------|-------|------------|
| SCUL-72-1B-SS | 72"    | 30"   | 1          |
| SCUL-72-2B-SS | 72"    | 30"   | 2          |
| SCUL-96-1B-SS | 96"    | 30"   | 2          |
| SCUL-96-2B-SS | 96"    | 30"   | 3          |



| Item Number    | Length | Depth | # of Sinks |
|----------------|--------|-------|------------|
| RFDS-72-1B-SS  | 72"    | 33"   | 1          |
| RFDS-72-2B-SS  | 72"    | 33"   | 2          |
| RFDS-96-2B-SS  | 96"    | 33"   | 2          |
| RFDS-96-3B-SS  | 96"    | 33"   | 3          |
| RFDS-120-2B-SS | 120"   | 33"   | 2          |
| RFDS-120-3B-SS | 120"   | 33"   | 3          |

# Cylinder Racks

## CYLINDER RESTRAINT RACKS



SCR02-1-SS Cylinder Rack  
Shown

| Width | Item Number | # of Cylinders |
|-------|-------------|----------------|
| 16"   | SCR01-1-SS  | 1              |
| 30"   | SCR02-1-SS  | 2              |
| 44"   | SCR03-1-SS  | 3              |
| 58"   | SCR04-1-SS  | 4              |
| 72"   | SCR05-1-SS  | 5              |

\*Custom sizes available

\*Includes chains

\*\* Standard construction is welded.  
Strut/Bolt version also available for a  
value engineered option.



SCR02-2-SS Cylinder Rack  
Shown

| Width | Item Number | # of Cylinders |
|-------|-------------|----------------|
| 16"   | SCR01-2-SS  | 2              |
| 30"   | SCR02-2-SS  | 4              |
| 44"   | SCR03-2-SS  | 6              |
| 58"   | SCR04-2-SS  | 8              |
| 72"   | SCR05-2-SS  | 10             |

\*Custom sizes available

\*Includes chains

\*\* Standard construction is welded.  
Strut/Bolt version also available for a  
value engineered option.



SCR02-3-SS Cylinder Rack  
Shown

| Width | Item Number | # of Cylinders |
|-------|-------------|----------------|
| 16"   | SCR01-3-SS  | 3              |
| 30"   | SCR02-3-SS  | 6              |
| 44"   | SCR03-3-SS  | 9              |
| 58"   | SCR04-3-SS  | 12             |
| 72"   | SCR05-3-SS  | 15             |


\*Custom sizes available

\*Includes chains

\*\* Standard construction is welded.  
Strut/Bolt version also available for a  
value engineered option.

# Umbilicals

## ISLAND UMBILICALS




UMB-IR-22X5X72-SS Shown



| Height | Item Number      |
|--------|------------------|
| 68"    | UMB-IR-X-X-68-SS |
| 73"    | UMB-IR-X-X-73-SS |
| 78"    | UMB-IR-X-X-78-SS |
| 90"    | UMB-IR-X-X-90-SS |

## WALL UMBILICALS



UMB-WR-22X6X78-SS Shown



| Height | Item Number      |
|--------|------------------|
| 68"    | UMB-WR-X-X-68-SS |
| 73"    | UMB-WR-X-X-73-SS |
| 78"    | UMB-WR-X-X-78-SS |
| 90"    | UMB-WR-X-X-90-SS |



**Air Master Systems Corp.**

6480 Norton Center Drive, Muskegon, MI 49441 • 231.798.1111 • 231.798.4000 fax • [airmastersystems.com](http://airmastersystems.com)