

Round Tube Legs

14 Gauge Steel Inner Tube

16 Gauge Steel Outer
Tube

Powder Coated Steel
Frame


Adjustable Height

MODULAR COMPACT TABLE

OVERVIEW

The modular compact table is the perfect addition to any classroom. Available with round tube legs with levelers. The telescoping legs, allow for easy height adjustment from 29" to 37" in 1" increments. Telescoping legs are secured to a box frame for added strength and security. The modular compact table comes with a 3/4" thick black lab grade phenolic resin surface. When selecting a butcher block top the minimum depth is 25" and the thickness is 1.5"

A variety of width and depth options to maximize your classroom and provide students with ample workspace. Add casters to the square leg tables for even more flexibility and easy reconfigurations.



MATERIALS

The telescoping legs are constructed of 16 gauge 2" diameter outer round tube with 14-gauge 1.75" diameter inner round tube steel. The round leg has a box frame secured to the legs to support the surface. The tables can be knocked down for cost effective shipping, or shipped fully assembled. If shipped knocked down minimal assembly required on site, hardware provided.

** Available colors:

Ink Black, Summer Green, Engine Red, Royal Purple, Sky Blue, Tiger Orange, Petal White, Bankers Gray, Model Gray



Modular Compact Table Foot Depth 3/4" Black Phenolic Surface

MCT-RL4830-PBLKL

Leg Shape Width

Key:

Legs:

R - Round Legs

Feet:

L - Levelers

Surfaces:

B - Butcher Block (1.5" Thick)

P - Phenolic

E - Epoxy

L - Laminate

Materials:

BBWT - Butcher Block Wood Top (1.5" Thick)

BLKL - Black (Phenolic, Epoxy, Laminate)

ASHL - Ash (Laminate)

GRPL - Storm (Laminate)

FWLL - Formal Walnut (Laminate)

NTAL - Natural Ash (Laminate)

BKOM - Bianco Kos (Laminate)

*Butcher Block surface depths are **25"**, 30" and 36"

Part Number	Leg Shape	Foot	Surface
MCT-RL4824-PBLKL	Round	Leveler	3/4" Phenolic
MCT-RL4830-PBLKL	Round	Leveler	3/4" Phenolic
MCT-RL4836-PBLKL	Round	Leveler	3/4" Phenolic
MCT-RL6024-PBLKL	Round	Leveler	3/4" Phenolic
MCT-RL6030-PBLKL	Round	Leveler	3/4" Phenolic
MCT-RL6036-PBLKL	Round	Leveler	3/4" Phenolic
MCT-RL7224-PBLKL	Round	Leveler	3/4" Phenolic
MCT-RL7230-PBLKL	Round	Leveler	3/4" Phenolic
MCT-RL7236-PBLKL	Round	Leveler	3/4" Phenolic


REV-003



Black



Ash



Storm



Formal Walnut






Natural Ash



Bianco Kos

*Additional Laminates Available Upon Request.



Scan for more info