

Square 2x2 Steel Legs

Casters

Levelers

Powder Coated Steel
Frame

Adjustable Height

MODULAR COMPACT TABLE

OVERVIEW

The modular compact table is the perfect addition to any classroom. Available with square tube legs and on casters or levelers. The telescoping legs, allow for easy height adjustment from 29" to 37" in 1" increments. Telescoping legs are secured to a box frame for added strength and security. The modular compact table comes with a 3/4" thick black lab grade phenolic resin surface. When a butcher block top is selected the minimum depth is 25" and the thickness will be 1.5"

A variety of width and depth options to maximize your classroom and provide students with ample workspace. Add casters to the square leg tables for even more flexibility and easy reconfigurations. Note: Casters only available on square leg table.

MATERIALS

The telescoping legs are constructed of 16 gauge, 2" x 2" square tube outer legs with 14 gauge 1.75" x 1.75" square tube inner legs. The square leg tables have a box frame secured to the legs for stability and full support of the surface. The modular compact table can be knocked down for cost effective shipping, or shipped fully assembled. Minimal assembly required on site with knockdown option, hardware provided.



****MCT-SC4830-PBLKL**

** Available colors:

Ink Black, Summer Green, Engine Red, Royal Purple, Sky Blue, Tiger Orange, Petal White, Bankers Gray, Model Gray



Modular Compact Table Foot Depth 3/4" Black Phenolic Surface

MCT-SC4830-PBLKL

Leg Shape Width

Key:

Legs:

S - Square Legs

Feet:

C - Casters

L - Levelers

Surfaces:

B - Butcher Block (1.5" Thick)

P - Phenolic

E - Epoxy

L - Laminate

Materials:

BBWT - Butcher Block Wood Top (1.5" Thick)

BLKL - Black (Phenolic, Epoxy, Laminate)

ASHL - Ash (Laminate)

GRPL - Storm (Laminate)

FWLL - Formal Walnut (Laminate)


NTAL - Natural Ash (Laminate)

BKOM - Bianco Kos (Laminate)

Part Number	Leg Shape	Foot	Surface
MCT-SL4824-PBLKL	Square	Leveler	3/4" Phenolic
MCT-SL4830-PBLKL	Square	Leveler	3/4" Phenolic
MCT-SL4836-PBLKL	Square	Leveler	3/4" Phenolic
MCT-SL6024-PBLKL	Square	Leveler	3/4" Phenolic
MCT-SL6030-PBLKL	Square	Leveler	3/4" Phenolic
MCT-SL6036-PBLKL	Square	Leveler	3/4" Phenolic
MCT-SL7224-PBLKL	Square	Leveler	3/4" Phenolic
MCT-SL7230-PBLKL	Square	Leveler	3/4" Phenolic
MCT-SL7236-PBLKL	Square	Leveler	3/4" Phenolic
MCT-SC4824-PBLKL	Square	Caster	3/4" Phenolic
MCT-SC4830-PBLKL	Square	Caster	3/4" Phenolic
MCT-SC4836-PBLKL	Square	Caster	3/4" Phenolic
MCT-SC6024-PBLKL	Square	Caster	3/4" Phenolic
MCT-SC6030-PBLKL	Square	Caster	3/4" Phenolic
MCT-SC6036-PBLKL	Square	Caster	3/4" Phenolic
MCT-SC7224-PBLKL	Square	Caster	3/4" Phenolic
MCT-SC7230-PBLKL	Square	Caster	3/4" Phenolic
MCT-SC7236-PBLKL	Square	Caster	3/4" Phenolic

REV-003


*Butcher Block Tops Depths are **25"**, 30 and 36"



Black



Ash




Storm



Formal Walnut





Natural Ash



Bianco Kos

*Additional Laminates Available Upon Request.



Scan for more info